



Bed Techs Inc.

Specializes in Reconditioned Hill Rom® Beds & Stretchers

Hill Rom® Century Plus 835/ 837

Hill Rom® Beds for Safety, Quality, Comfort. Great Products—Intelligent, Well-Designed, Clinically Effective and Easy to Use.

Our Process Includes

- Professionally refurbished medical beds to manufacturer's specifications
- JCAHO required documentation
- New paint & labels
- Fully function tested
- Basic bed function (Head, Knee, Hi-low)
- Spring sleep surface (Standard) or Orthopedic/Hard Pan Surface (Upgrade)
- Reconditioned PC boards, 1 year warranty
- Fully rebuilt side rails w/new plastics and switches
- New head/foot boards, your choice of laminate colors, lifetime warranty
- New paint, Taupe (Standard) or Light Neutral (Upgrade)
- New motor couplings and capacitors
- Rebuilt casters with new components and wheels
- Delivery/Installation performed by our technicians (Domestic only)
- 1 year warranty (Domestic only)
- Quality & Price

Specifications

- Sleep Surface: 36" x 80"
- Overall Length (w/out bumpers): 36.5" x 88"
- Bed Frame Clearance w/ 5" casters: 7.5"
- Low position: 17"
- High position: 30"
- Central Brake and Steer

Options

- Radio/TV
- Bed exit
- Nurse call
- New PC boards, 5 year warranty
- Parallelogram (new style) side rails
- Orthopedic (Hard pan) surface
- Mattress (built to OEM spec)
- Roller bumpers
- Service Contracts for preventative maintenance
- 220-volt transformers

