

650 Ford Street Colorado Springs, CO 80915 Toll:(800)433-4064 * Ph: (719)637-0300 Fax:(719)597-8273 * www.harloff.com

Six Drawer Mini Line Emergency Cart with Breakaway Lock, Specialty Accessory Package



Drawer Configuration:

- (3) 3-inch (7.6 cm)
- (2) 6-inch (15.2 cm)
- (1) 9-inch (22.9 cm)

Sixteen available color choices:

*Antimicrobial Finish Available

Burgundy, Hammer Red, Red, Yellow, Beige, Sand, White, HammerTone Green, Teal, HammerTone Blue, Navy, Light Blue, HammerTone Gray, Light Gray, Mauve, Hammer Black

Specifications

Construction:

· Painted steel cabinet and drawers

Features:

- Space-saver design fits into smaller areas
- Allows supply storage in areas with space limitations
- Cabinet with 30 vertical inches (76.2 cm) drawer space
- Break-away locking system with package of 100 plastic seals
- 3" (7.6 cm) single-wheel stem caster set, full swivel two with brakes
- Pre-threaded accessory mounting holes
- · Ball-bearing full extension drawer slides
- Removable plastic top
- Latex free
- Twelve year warranty
- Cart requires caster installation and accessory assembly after delivery
- Durable hard-baked powder coat finish

Accessory Package Includes:

- 2-prong extendable IV pole (680401)
- Oxygen tank holder (AL680403)
- Cardiac board and brackets (680406)
- Aluminum push handle / side rail (68161)

Internal Drawer Dimensions:

• 12.6" W x 16.75" D (32 cm x 42.5cm)

Assembled Cart Dimensions

- 66.25" H x 24.5" W x 20.5" D (168.3 cm H x 62.2 cm W x 52 cm D) height includes IV Pole - extends to 84.75" / 215.3 cm H
- Cart net weight: 96 pounds (43.5 kg)
- Approximate shipping weight: 116 pounds (52.6 kg)