

Pediatrics > Head Supports

Tumble Forms²

Carrie[®] Collar with Tumble Forms 2[®] Antimicrobial Protection. This hinged collar supports the jaw line and occipital region to help improve head control, neck flexion and midline positioning. Can be used by itself or used with other seating systems. Designed to fit the Elementary size Carrie Seat System. Attaches in the front with a quick-release hook & loop fastener. The Carrie Collar is odor resistant with increased durability. The seamless covering is impervious to fluids and easy to clean. Fits 14" neck circumference.

Note: The Carrie Collar should not be used as a child restraint in transportation.

2774C



Danmar Products
Equipping people. Enhancing lives.

Hensinger Head Support. Soft foam wrap-around collar supports the head by using base support areas at the occiput and jaw line. After the rest of the body is positioned, the support can help maintain mid-line head positioning, improve breathing and swallowing, and simplify feeding and social interaction. Built with strong yet pliable reinforcement material around the outside, and a quick-release buckle for easy on/off. Soft sock covers are easy and inexpensive to replace and clean.

Black	Blue	Size	Circumference
081627660	924550	XS	10" / 1" (Neck Circumference / Thickness*)
081627678	924551	S	12" / 1½" (Neck Circumference / Thickness*)
081627686	924552	M	14" / 2" (Neck Circumference / Thickness*)
081627694	924553	L	16" / 2½" (Neck Circumference / Thickness*)

*Measure bottom of ear lobe to top of shoulder (minus ½" so that it does not rub) for thickness.



Pediatrics

Pediatric Cervical Collar. Soft foam treated with Microban[®] antimicrobial protection for a fresher, more hygienic product. Low contour material maintains head in neutral position. Simple application with hook & loop closures. Latex free.



Item No.	Size	Approx. Age	Neck Circ.	H x L
56474901	Infant	0–24 mos	9"–10½"	2" x 13"
56474902	Ped	2–7 yrs	10"–12"	2½" x 15"
56474903	Youth	8–12 yrs	12"–14"	2½" x 17"

Anterior Head Support.

Closed front and reinforced sides allow firm but gentle lifting of the chin and lateral support of the head (not to be attached or mounted to a seat or headrest). Padded roll around the rear strap ensures a comfortable fit. Vinyl-covered foam with soft cloth cover included. Size by neck circumference: **A.** Neck circumference **B.** Shoulder to bottom of ear, allowing 1" between ear-lobe and top of head support.



Item No.	Size	A	B
924557	Small	9"–12"	2–3"
924558	Medium	11"–14"	2½"–4"
924559	Large	13"–17"	3½"–5"