

Hill-Rom® Electric Stretcher



Go Beyond Mobility™

Hill-Rom's Electric Stretcher provides the comfort of a hospital bed for those unexpected long stays. With positioning controls integrated into the siderails, patients are in control of their own comfort.

A rush to the operating room is not just about running down hallways and turning corners. It's about the minute details that contribute to a safe and comfortable experience for both you and your patient.

Our stretchers **provide more than mobility** for your patients:

- A safe and comfortable patient experience
- Simple and safe maneuvering
- Durability to withstand your demands


Enhancing outcomes for
patients and their caregivers:

Hill-Rom


A Safe and Comfortable Patient Experience

Patients are spending more time on stretchers, and as a result, the patient experience, including pressure injury prevention, has become increasingly important.

- **Help prevent pressure injuries** through your choice of pressure redistribution surfaces.
- **Customize comfort** with electric patient positioning controls.
- **Prevent patient migration** with electric knee articulation.
- **Protect your patient's hands** with bumper protected siderails.
- **Provide a gentle patient ride** with spring-loaded Steering Plus™ system that absorbs shock of driving over bumps and thresholds.
- **Allow for safe patient ingress and egress** with low height.
- **Provide a comfortable option** for patients to recover with Recovery Position.



Recovery position



Simple and Safe Maneuvering

Whether you're making crisp turns down narrow halls, getting the stretcher rolling from a stop, or nestling it into a crowded elevator, our stretchers are ideal in any setting. Hill-Rom Stretchers are easy to maneuver to help reduce your risk of injuries.

- **Simple activation** requires less force to activate than Stryker's Big Wheel.
- **360° maneuvering** with Steering Plus™ system grips the floor and reduces your turning radius.
- **Ease the burden off your back** with electric positioning controls.
- **Zero Transfer Gap** — OneStep® Tuck-Away Siderails for seamless patient transfers.
- **Obtain an immediate weight without moving your patient** with scale that memorizes last weight.

Durability to Withstand your Demands

At Hill-Rom, we understand that because caregivers have to do more with less, having a reliable and durable product is more important than ever.

- **700 lb weight capacity** for all stretchers
- **Compatibility** with lifts and overbed tables
- **Easy to clean** smooth deck and siderail design

Go Beyond Mobility™


A safe and comfortable experience for both you and your patient starts on a Hill-Rom® stretcher.



Patient Controls



Caregiver Control Panel



AccuMax® VPC Surface



Dual-Display Scale



Steering Plus™ System



Active Hand Brake

A SAFE & COMFORTABLE PATIENT EXPERIENCE	STANDARD FEATURES	OPTIONAL FEATURES
30" width	X	
Electric Fowler Articulation	X	
Electric Knee Articulation	X	
Bumper Protected Siderails	X	
Trendelenburg/Reverse Trendelenburg	X	
Six Transport Strap Locations	X	
PRESSURE REDISTRIBUTION MATTRESS OPTIONS		
4" Comfortline® foam mattress	X	
5" Comfortline® foam mattress		X
5" AccuMax® VPC Surface (Elastimax or Hypolex)		X
SIMPLE & SAFE MANEUVERING		
Steering Plus™ System	X	
OneStep® Tuck-Away Siderails	X	
Caregiver Electric Fowler Control	X	
Caregiver Electric Knee Control	X	
Caregiver Fowler & Knee Lockout Controls	X	
Brake/Steer – head and foot	X	
Brake/Steer – all four sides		X
Dual-Sided hi/lo Hydraulic Control Pedals	X	
Ergonomic Push Handles	X	
Ergonomic Push Handles with Integrated IV Pole Transport Holder		X
Active Hand Brake		X
3/4" Standard IV Pole Sockets	X	
Removable IV Pole – 2 Stage		X
Permanent IV Pole – 3 Stage (40 lb weight capacity)		X
Integrated Dual-Display Scale		X
DURABILITY TO WITHSTAND YOUR DEMANDS		
700 lb weight capacity	X	
Full Perimeter Bumpering	X	
Smooth Deck Design	X	
Decals Options	X	
Bumper Color Options – choice of neutral, blue, red, teal or purple	X	

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

©2016 Hill-Rom Services, Inc. ALL RIGHTS RESERVED.
197892 rev 1 29-JUN-2016 ENG – US

For further information about this product or a service, please contact your local Hill-Rom representative or visit our webpage:

USA 800-445-3730
Canada 800-267-2337

www.hill-rom.com

Technical Specifications

Maximum weight limit	700 lbs (317.5 kg)
Low position (floor to top of sleep deck)	22.5" (57.2 cm)
High position (floor to top of sleep deck)	34.25" (87 cm)
Overall length	83" (210.8 cm)
Overall width (siderails up)	36" (91.44 cm)
Overall width (siderails stored)	30.5" (77.5 cm)
Siderail length	48" (121.9 cm)
Siderail height above sleep deck	14.5" (36.8 cm)
Comfortline® mattress size	26" (66 cm) x 75" (190.5 cm) or 30" (76.2 cm) x 75" (190.5 cm)
AccuMax® VPC mattress size	30" (76.2 cm) x 75" (190.5 cm)
Maximum head elevation	90°
Maximum Trend/Reverse Trend	18°
Floor to base clearance	3.5" (8.9 cm)
Caster size	8" (20.3 cm)

Warranty

Frames & welds for life of product

Three-year parts

Two-year mattress

One-year labor

Accessory Options — Please confirm accessory availability for model of stretcher ordered

Accessory IV pole	Liquid oxygen tank holder
Ankle stirrups	Oxygen tank holder
Armboard	Paper roll dispenser
Chart holder	Patient tray
Convertible footboard	Restraint straps
Detachable footboard	Siderail covers
Foot extender	Stretcher pillow
Infusion support system	Transport straps
IV transporter	Utility tray
Lateral X-ray cassette holder	

Enhancing outcomes for patients and their caregivers:

Hill-Rom