



**DEVICE SPECS**

# Ethicon Linear Cutter Portfolio



**ETHICON™ Linear Cutter**



**PROXIMATE™ Linear Cutter**

<b>Codes</b>	NTLC55, NTLC75	TVC55, TLC55, TCT55, TLC75, TCD75, TCT75, TLC10, TCT10
<b>Associated Reloads</b>	SR55, SR75	TVR55, TCR55, TRT55, TCR75, TRD75, TRT75, TCR10, TRT10
<b>Description</b>	Selectable Linear Cutter with 3D Staples	Linear Cutter with Reload
<b>Device Technology</b>	Selectable Staple Height	Cam Closure Mechanism
<b>Size</b>	55mm, 75mm	55mm, 75mm, 100mm
<b>Staple Rows</b>	6	4
<b>Closed Staple Height</b>	1.5mm, 1.8mm, 2.0mm	1.5mm, 1.8mm, 2.0mm
<b>Max Firings</b>	12	8
<b>Staple Technology</b>	3D	B-Form
<b>Qty./Box</b>	3	3