

## SOOTHE-A-CISER™ CERVICAL TRACTION CUSHION

- Helps alleviate headaches and neck pain
- Gently stretches upper chest and neck muscles
- Lumbar extension comfortably supports upper thoracic spine
- Helps restore the natural curve of the neck
- Promotes recovery from whiplash injury
- Ideal for patients with tech neck or military neck



Available From:

## HELP RELIEVE HEADACHES AND NECK PAIN

So many daily activities affect our posture, placing the head and neck in a detrimental forward head position. This can be attributed to looking down at handheld devices, sleeping with the head elevated excessively, poor posture, or an acute injury such as whiplash. Using the Soothe-A-Ciser™ helps alleviate tension headaches and neck pain by gently stretching the upper chest and anterior neck muscles.

Helps provide relief of headaches and neck pain associated with whiplash injury, tech neck, and other neck injuries by counteracting forward head posture. The lumbar extension provides firm support, placing the neck in a proper position to help reestablish proper cervical lordosis, reducing forward head posture.

## COMFORTABLE DESIGN

The contoured precision cut foam and unique lumbar extension provide support to the upper thoracic spine for comfort and also allow the shoulders to drop, stretching the chest and neck muscles. Opening the chest in this manner promotes breathing and facilitates relaxation. The Soothe-A-Ciser is an effective tool to assist healthcare professionals in restoring a patient's proper cervical lordosis.

Massage therapists appreciate the contoured design that elevates the head, shoulders, and neck for better access during massage. This cervical traction cushion provides comfortable support for clients while placing the neck in mild traction.

## **SPECIFICATIONS**

Product Selection: Soothe-A-Ciser™ Cervical Traction Cushion BAK-440, Black

Product Weight: 1.1 lbs

**Dimensions:** 11" x 20.5" x 6.75"

Materials: Polyurethane foam, cotton,

polyester

Case Quantity: 4

Case Dimensions: 15" x 14" x 22.5"

Case Weight: 6 lbs

**Latex:** This product is not made with natural rubber latex

Country of Origin: Made in USA

California Residents Only:

▲WARNING: Cancer and Reproductive Harm — www.P65Warnings.ca.gov