

SWEDE-O® STRAP LOK® ANKLE BRACE

- Easily control support with external figure eight straps
- Stays tight and secure with exclusive lacing system and top strap
- Cool, comfortable, breathable tongue
- Durable, lightweight ballistic nylon is ideal for sports
- Full elastic back reduces irritation and potential for blistering
- Bilateral design fits both right and left ankle





Available From:



MAXIMUM PROTECTION & SUPPORT

The Swede-O® Strap Lok® ankle brace has external figure eight straps which allow you to easily control and customize the amount of support as well as simulate professional taping.

The exclusive dual eyelet lacing system is offset to create a "pulley effect", keeping the brace tighter longer. The top strap covers the laces and locks the figure eight straps in position. The full elastic back is comfortable against the Achilles tendon and helps reduce irritation and the potential for blistering.

The breathable tongue helps keep the foot cool and comfortable. The Strap Lok is backed by an industry leading I year warranty and constructed from a durable ballistic nylon which is comfortable thin, and lightweight.

INDICATIONS FOR USE

Use the Strap Lok after injury or to protect your ankle from future injury. This ankle brace helps prevent and relieve aches, pains, sprains, and other ankle conditions.

SPECIFICATIONS

Product Selection: Swede-O® Strap

Lok® AKL-6334, Black

Product Weight: .35 lbs

Dimensions:

XSmall: 7.375" x 8"

Small: 8.75" x 8.75" Medium: 9" x 9"

Large: 9" x 9.5" XLarge: 9" x 9.5"

2XLarge: 9.5" x 10"

Materials: Polyester, nylon, polyurethane foam, neoprene, brass

Case Quantity: 12

Case Dimensions: 11" x 14" x 13"

Case Weight: 7 lbs

 $\textbf{Latex:} \ This \ product \ is \ not \ made \ with$

natural rubber latex

Country of Origin: Made in USA

California Residents Only:

▲WARNING: Cancer and Reproductive Harm — www.P65Warnings.ca.gov