

# Spec Sheet

## LINKTEMP<sup>®</sup>

### Non-Contact Infrared Thermometer

<b>Device Type</b>	Non-contact infrared thermometer
<b>Measurement Distance</b>	1.2" to 2" (3 to 5 cm) from forehead
<b>Measuring Range</b>	Body temperature: 89.6°F to 109°F (32°C to 42.9°C) Surface temperature: 32°F to 140°F (0°C to 60°C) Ambient temperature: 32°F to 104°F (0°C to 40°C)
<b>Displayed Temperature Range</b>	89.6°F to 109.2°F (32°C to 42.9°C)
<b>Display Type</b>	Large, 2-color LCD with highly accurate 1/100 degree reading
<b>Display Resolution</b>	0.1°F (0.1°C)
<b>Accuracy</b>	89.6°F to 94.8°F=±0.6°F (32°C to 34.9°C=±0.3°C) 95°F to 107.6°F=±0.4°F (35°C to 42°C=±0.2°C) 107.8°F to 109.2°F=±0.6°F (42.1°C to 42.9°C=±0.3°C)
<b>Operating Temperature</b>	50°F to 104°F (10°C to 40°C)
<b>Operating Humidity</b>	Up to 85% non condensing
<b>Storage Temperature</b>	-4°F to 131°F (-20°C to 55°C)
<b>Storage Humidity</b>	Up to 95% non condensing
<b>Pressure Altitude</b>	700 hPa - 1060 hPa
<b>Consumption</b>	≤450 mW
<b>Calibration</b>	Preset at factory. Device re-calibrates each time device is turned ON
<b>Memory</b>	Stores up to 32 readings
<b>Automatic Power Off</b>	30 seconds after last measurement
<b>Batteries</b>	2 AA batteries
<b>Unit Weight</b>	3.48 oz without batteries
<b>Unit Size</b>	6.81" x 2.05" x 1.77" (173 x 52 x 45 mm) in L x W x H
<b>Warranty</b>	18 month limited warranty



For more information on the **LinkTemp™** Non-Contact Infrared Thermometer, please contact your local distributor or Links Medical Products customer service at (888) 425-1149.



Links Medical Products Inc.  
9247 Research Drive  
Irvine, CA 92618 USA  
[www.linksmed.com](http://www.linksmed.com)