



LT-980

High strength Aluminum ultra-lightweight wheelchair

- Ultra-light weight Aluminum Frame
- Weighs only 24 lbs. (w/o riggings)
- Folding backrest for easy travel or storage
- Fixed armrests with molded arm pads
- Composite side guards
- Swing-away/detachable footrest
- Height adjustable footplate
- Adjustable push to lock wheel brakes
- 6" x 1 " Flat free front polyurethane casters
- 22" x 1 3/8" Spoke wheels with high profile/polyurethane tires
- Easy-Grip push-rims
- Carry pocket on backrest
- Easy to clean Nylon upholstery
- Seat height adjustability
- Weight capacity: 250 lbs
- Frame colors available: Black or Merlot Mica
- Lifetime Limited Warranty

The Karman Ultra Lightweight LT-980 is by far the lightest wheelchair at its price point. Weighing in at only 24 lbs, this chair is ideal for those who are traveling, caretakers who are tired of lifting a standard steel wheelchair, or any wheelchair user who is looking for something easier to propel and live life independently.

SPECIFICATIONS (subject to change without notice): LT-980-BK (Black Frame) LT-980-BD (Merlot Mica)

Model No.	seat width (in.)	seat depth (in.)	armrest height (in.)	seat height (in.)	back height (in.)	overall height (in.)	overall open width (in.)	folded width (in.)	overall length (in.)	weight w/o riggings (lb.)	weight capacity (lb.)
LT-980-BK	18	16	9	18-19	19	38	26	13	28	24	250
LT-980-BD	18	16	9	18-19	19	38	26	13	28	24	250
LT-980W-BK	20	16	9	18-19	19	38	28	13	28	24	250



KARMAN HEALTHCARE Inc.

19255 San Jose Avenue . City of Industry . CA 91748-1418

Order. 1-800-805-2762 Tel. 626-581-2235 Fax. 626-581-2335