



FOLDING WALKERS

colorSelect[®] ADULT WALKER WITH WHEELS

Features

- Sturdy 1" aluminum tubing provides maximum strength while remaining lightweight
- Dual-release folding mechanism enables user to fold walker sides independently
- 5" front fixed wheels
- Complete with glide tips that replace rear leg rubber tips to enhance walker maneuverability on all surfaces
- Designed to meet the requirements of HCPCS code: E0143
- Maximum Weight Capacity 300 lb - EVENLY DISTRIBUTED
- One (1) year limited warranty



Color Options

		Width inside handgrips		Width inside rear legs		overall Width		Depth front to rear		Depth folded		minimum Height		maximum Height		Weight Capacity EVENLY DISTRIBUTED		UOM
		inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	lb	Kg	
716270B-2	Blue	17	43.18	20	50.80	25.5	64.77	20.25	51.44	4	10.16	32	81.28	39	99.06	300	136.08	2/cs
716270BK-2	Black	17	43.18	20	50.80	25.5	64.77	20.25	51.44	4	10.16	32	81.28	39	99.06	300	136.08	2/cs
716270G-2	Gold	17	43.18	20	50.80	25.5	64.77	20.25	51.44	4	10.16	32	81.28	39	99.06	300	136.08	2/cs
716270P-2	Plum	17	43.18	20	50.80	25.5	64.77	20.25	51.44	4	10.16	32	81.28	39	99.06	300	136.08	2/cs
716270PK-2	Pink	17	43.18	20	50.80	25.5	64.77	20.25	51.44	4	10.16	32	81.28	39	99.06	300	136.08	2/cs
716270R-2	Red	17	43.18	20	50.80	25.5	64.77	20.25	51.44	4	10.16	32	81.28	39	99.06	300	136.08	2/cs
716270B-2S	Silver	17	43.18	20	50.80	25.5	64.77	20.25	51.44	4	10.16	32	81.28	39	99.06	300	136.08	2/cs



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2011, GF Health Products, Inc. All Rights Reserved.
Lumex is a registered trademark of GF Health Products, Inc.

graham
field

GF.1100059Rev.D22