



Mayo Stand Large

4366

1/EA

Product Features



WARNING: This product can expose you to chemicals including chromium (hexavalent compounds), which are known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov.

The 4366 features a stainless steel tray that offers 25% more surface space

Ideal for any market where a larger work space is needed

High grade stainless steel tray (18/8,304)

Base minimizes interference with other equipment

Secure height adjustments 37" - 53"

Tray size: 16-3/4" x 21-1/2"

Base width: 20", depth: 25"

Weight Capacity Standard: 50 lbs

Replacement tray available (4266)

Depth at base: 24-1/3"

Inside Base Area: 10-1/2" at widest end and 6-1/4" at narrow end

Attributes

Brand Name	Tech-Med
Sterile	No
Weight	18.29837 lbs

Packaging

Packaging	1/EA
Case Packaging	1/EA
Total Case Quantity	1
Case Packaging	(01)20665973022295