



Mayo Stand

4352

1/EA

Product Features



WARNING: This product can expose you to chemicals including chromium (hexavalent compounds), which are known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov.

Features a high grade stainless steel tray (18/8,304)

U-shaped base minimizes interference with other equipment

Secure height adjustments 28" - 48"

Tray Size: 12-1/2" x 19"

Base width: 20"

Weight Capacity: 11 lbs

Replacement tray available (4264)

Attributes

Brand Name	Tech-Med
Sterile	No
Weight	11.99315 lbs

Packaging

Packaging	1/EA
Case Packaging	1/EA
Total Case Quantity	1
Case Packaging Barcode	(01)20665973022240