

# RL-MIS-20

Mechanical Baby Scale

HEALTH SCALES



## Standard Features

- Painted beam
- Sturdy, enameled steel body
- Easy-to-clean ABS plastic tray
- Unit of measure: kg
- Built-in measuring tape 0 to 58 cm

## Specifications

**Capacity:**

70 kg x 20 g

**Tray Dimensions:**

(L x W x H)

61 x 41 x 7.5 cm

## Part Number/Price

Part #	Model #	Capacity	Price
131471	RL-MIS-20	70 kg x 20 g	Consult

## Shipping Weight

Part #	Model #	L x W x H	Package Wgt. (lb)	Dim Wgt. (lb)	Shipping Wgt. (lb)
131471	RL-MIS-20	28 x 20 x 14 in	30	57	57

## Dimensions

Centimeters	
A	61.0
B	68.6
C	54.5

