

MICRO-TOUCH[®] Plus

Essential Series

Latex, powder-free, medical exam glove with excellent tactile sensitivity and dexterity

- **Powder-free:** These non-sterile gloves are powder-free to eliminate powder-related glove complications
- **Secure grip:** The fully textured outer surface ensures a secure grip in wet or dry conditions and prevents slippage
- **Heightened sensitivity:** Provides the wearer with superior tactile sensitivity
- **Easy donning:** The inner polymer coating allows for ease of donning
- **Comfortable fit:** Due to the natural elastic properties of latex, this glove has excellent elasticity and memory, offering a more secure fit that provides greater comfort without bagginess
- **Smart packaging:** The 150-count packaging is smartly packed for reduced waste

CAUTION: This product contains natural rubber latex, which may cause allergic reactions

Recommended For

- Brief, low-risk examination procedures
- Environmentally friendly facilities



MICRO-TOUCH[®] Everyday, Everywhere Protection

MICRO-TOUCH[®] premium medical examination gloves are scientifically designed to complement every need of medical professionals. Sophisticated design features include Ansell ergonomic technologies that ensure comfort during extended procedures, wet grip innovations that prevent harmful mishandling, and micro-thin films for maximum sensitivity. Caregivers who are fully prepared and comfortably protected by versatile MICRO-TOUCH styles have the peace of mind to focus on what matters most: the health and protection of themselves and their patients.

MICRO-TOUCH® Plus

Essential Series

PRODUCT INFORMATION

| | | | |
|---|----------------------|------------------------------------|--|
| Material | Natural Rubber Latex | Not Made From Natural Rubber Latex | No |
| Color | Cream | Cuff Length | Standard |
| Powder Content | Powder-Free | External Glove Surface | Fully Textured |
| Freedom from Holes (Inspection level I) | 1.5 AQL | Palm Thickness (mm/mil) | 0.12 / 4.7 |
| Finger Thickness (mm/mil) | 0.16 / 6.3 | Available Sizes | XS (5.5 - 6), S (6.5 - 7), M (7.5 - 8), L (8.5 - 9), XL (9.5 - 10) |
| Tested For Use With Chemotherapy Drugs | No | Sterile | No |
| Antistatic | Not Tested | Glove Length (mm/inches) | 230 / 9.1 |
| Product Segmentation | Essential Series | | |

ORDERING INFORMATION

| Size | XS (5.5 - 6) | S (6.5 - 7) | M (7.5 - 8) | L (8.5 - 9) | XL (9.5 - 10) |
|--------------|--------------|-------------|-------------|-------------|---------------|
| Product Code | 6015300 | 6015301 | 6015302 | 6015303 | 6015304 |

PACKAGING AND STORAGE

| | |
|----------------------|--|
| Packaging | 150 gloves per dispenser 10 dispensers per case 1500 gloves per case |
| Shelf Life | 3 Years |
| Storage Instructions | Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition. |

STANDARDS AND CERTIFICATIONS

ASTM D3578, EN 420:2003 + A1:2009, FDA21 CFR 177-2600-US Food Contact Approved, ISO 9001, ISO 10993-10:2002, ISO 13485

Contact your Ansell representative for ordering or more information.

North America
US Tel: 800 952 9916
CA Tel: 1-844-494-7854
Email: insidesalesus@ansell.com

Latin America & Caribbean
Tel: +52(442) 296 20 50
Email: cslac@ansell.com
Brazil Tel: +55-11-3356-3100
Brazil Email: luvas.medicas@ansell.com

Europe, Middle East & Africa
Tel: +32 (0) 2 528 74 00
Email: info@ansell.eu

Asia Pacific
Tel: +603 8310 6688
Email: apac.medical@ansell.com

Australia & New Zealand
Tel: +61 3 9270 7270
Email: protection@ap.ansell.com