

DIGITAL HOLTER BY MIDMARK



IQ
holter



The Midmark IQholter® System combines a compact digital recorder with three powerful software applications. Our three tiers of Holter solutions are specifically designed to meet different practice needs, from family physicians to cardiologists.

IQholter® Technology by Midmark

All the features you need in a solution tailored just for you.

A POWERFUL HOLTER SYSTEM DESIGNED FOR YOUR PRACTICE

IQholter® Recorder:

- Small and lightweight, the recorder weighs only 2.8 oz.
- Well-lit display and intuitive keypad make clinician navigation and menu control easy
- On-screen leads and channel identification for speed and simplicity during patient hookup
- Automatic startup self-test with warning alarm ensures reliable recording
- 5-lead hookup for 3-channel recording; 7-lead option available
- 24-hour standard recording (48 and 72-hour recording capabilities available)
- Removable Secure Digital (SD) card allows for immediate patient turnaround

IQholter®:

- Color-coded arrhythmia events to identify different ECG morphologies
- Comprehensive reports with auto narrative summary to effectively reveal important test data and statistics
- Ability to edit and review test results and full disclosure locally or on the network without printing

IQholter® EX:

- Provides the ability to preview and adjust criteria settings to match patient's recording before the analysis
- Identifies Atrial Fibrillation (AF) events based on parameters set by the clinician

- Smart page scan to highlight specified events and reveal entire episodes at a glance
- Arrhythmia templates with intelligent sorting to facilitate quick reviews and edits
- QT and QTc analysis for identifying prolonged intervals

IQholter® EP:

- Offers all the standard features of IQholter® and IQholter® EX
- Includes pacemaker detection to identify single and dual-chamber paced beats for non-captures and sensing failures
- Automatic pacemaker detection eliminates setting changes between patients



The IQholter® Recorder can be worn as comfortably as an MP3 player or cell phone.



Comprehensive report with auto narrative summary

Designed with a well-lit display and intuitive keypad for ease of use



Easy to use SD card allows for immediate patient turnaround



Fully automated analysis and reporting process

Make your practice more efficient with an IQ Digital Diagnostic System

Workflow efficiency and improvements in patient care are just a few of the benefits that can result from implementing an IQ Digital Diagnostic System. Call us today and we'll put together a customized suite of solutions to help grow as your practice grows.

Diagnostic Products				
	EHR/EMR Connectivity	Single Platform (IQmanager®)	Standalone	Thin-client (IQpath®)
IQecg®	■	■		■
IQholter®	■	■		■
IQholter® EX	■	■		■
IQholter® EP	■	■		■
IQstress®	■	■		
IQspiro®	■	■		■
IQvitals®	■	■	■	■
IQvitals® PC	■	■		■

Minimum System Requirements

- Windows® 8, Windows® 7; Professional and Enterprise; 32-bit and 64-bit
- Intel® Core™ 2 Duo Processor E4300 or 64-bit processor or faster
- 2 GB of RAM
- 2 GB available hard disk space
- 2 USB Ports
- Mouse, keyboard and high resolution display monitor

Contact us or visit midmark.com for more detailed system requirements.

Windows is a registered trademark of Microsoft Corporation.
Intel Core is a trademark of Intel Corporation in the U.S. and/or other countries.

Thin-client (IQpath®) System Requirements

Please contact your local Midmark Sales Representative for additional system requirements.

Ordering Information

4-000-0110	IQholter® Digital Holter with Recorder	1-370-0020	IQholter® Upgrade Kit (24, 48, 72 hours)
4-000-0113	IQholter® EX Digital Holter with Recorder	2-100-0080	IQholter® Prep Kit, 5-lead
4-000-0116	IQholter® EP Digital Holter with Recorder	2-100-0090	IQholter® Prep Kit, 7-lead
4-000-0119	IQholter® Digital Holter Recorder (5-lead standard)	2-100-0060	IQholter® Electrodes (box of 50)
4-100-1200	IQmanager® Software - Patient data management of IQecg®, IQspiro®, IQholter®, IQstress®, IQvitals® and IQvitals® PC	2-100-0050	IQholter® Recorder Pouches, Disposable (pack of 25)
1-370-0010	IQholter® Secure Digital Memory Card	2-100-0051	IQholter® Recorder Pouch with Belt, Reusable
		90-30-0205	IQholter® Lead Wire, 5-lead
		90-30-0206	IQholter® Lead Wire, 7-lead

For information, pricing or to place an order, contact us at:

1-800-624-8950 or
1-310-516-6050 or
midmark.com

Midmark is an ISO 13485 Certified Company

Manufactured and/or distributed by Midmark Corporation
690 Knox St, Suite 100, Torrance, CA 90502
midmark.com
1-800-624-8950

© 2011 Midmark Corporation
60 Vista Drive Versailles, Ohio 45380
Products subject to improvement changes without notice
Litho in U.S.A. 35-80-0001 Rev B. (7/15)



RETHINK THE CLINICAL SPACE

At the core of Midmark is a new way to think about caregiver/patient interaction and the delivery of efficient patient care. We are dedicated to working with you to redefine the future of the clinical space.

midmarkclinicalsolutions.com

