

# Integra® Miltex®

Female Patient Care Products


Limit uncertainty with our comprehensive selection of female patient care products



INTEGRA®  
LIMIT UNCERTAINTY

## Pessaries

- Soft and pliable
- Made of medical grade silicone\*
- Foldable for insertion and removal
- Variety of shapes and sizes available



### Ring

- Used to treat first or second-degree uterine prolapse
- Available with or without supporting membrane


	#0	#1	#2	#3	#4	#5	#6	#7	#8	#9
Size	1.75"	2.00"	2.25"	2.50"	2.75"	3.00"	3.25"	3.50"	3.75"	4.00"
With Support	30-RSo	<b>30-RS1</b>	<b>30-RS2</b>	<b>30-RS3</b>	<b>30-RS4</b>	<b>30-RS5</b>	<b>30-RS6</b>	<b>30-RS7</b>	<b>30-RS8</b>	30-RS9
Without Support	30-Ro	30-R1	<b>30-R2</b>	<b>30-R3</b>	<b>30-R4</b>	<b>30-R5</b>	<b>30-R6</b>	30-R7	30-R8	30-R9



### Ring with Knob

- Used to treat first or second-degree uterine prolapse
- May be used to treat stress urinary incontinence
- Available with or without supporting membrane


	#0	#1	#2	#3	#4	#5	#6	#7
Size	1.75"	2.00"	2.25"	2.50"	2.75"	3.00"	3.25"	3.50"
With Support	30-RKSo	30-RKS1	<b>30-RKS2</b>	<b>30-RKS3</b>	<b>30-RKS4</b>	<b>30-RKS5</b>	<b>30-RKS6</b>	30-RKS7
Without Support	30-RKo	30-RK1	30-RK2	<b>30-RK3</b>	<b>30-RK4</b>	30-RK5	30-RK6	30-RK7



### Gellhorn

- Used to treat second or third-degree uterine prolapse, cystocele or rectocele
- Drainage feature

	#0	#1	#2	#3	#4	#5	#6	#7	#8
Size	1.50"	1.75"	2.00"	2.25"	2.50"	2.75"	3.00"	3.25"	3.50"
Standard Stem	30-GDo	<b>30-GD1</b>	<b>30-GD2</b>	<b>30-GD3</b>	<b>30-GD4</b>	<b>30-GD5</b>	<b>30-GD6</b>	<b>30-GD7</b>	30-GD8



### Gellhorn Short-Stem

- Used to treat second or third-degree uterine prolapse, cystocele or rectocele
- Shorter stem may provide a better fit for some patients
- Drainage feature


	#0	#1	#2	#3	#4	#5	#6	#7	#8
Size	1.50"	1.75"	2.00"	2.25"	2.50"	2.75"	3.00"	3.25"	3.50"
Short Stem	30-GSo	30-GS1	<b>30-GS2</b>	<b>30-GS3</b>	<b>30-GS4</b>	<b>30-GS5</b>	<b>30-GS6</b>	<b>30-GS7</b>	30-GS8



### Donut

- Used to treat second or third-degree uterine prolapse, cystocele or rectocele
- Easily compressed for insertion and removal

	#0	#1	#2	#3	#4	#5	#6	#7
Size	2.00"	2.25"	2.50"	2.75"	3.00"	3.25"	3.50"	3.75"
Donut	30-Do	<b>30-D1</b>	<b>30-D2</b>	<b>30-D3</b>	<b>30-D4</b>	<b>30-D5</b>	<b>30-D6</b>	30-D7



### Cube

- Used to treat second or third-degree uterine prolapse, cystocele or rectocele
- Available with or without drains

	#0	#1	#2	#3	#4	#5	#6	#7	#8	#9	#10
Size	25 mm	29 mm	33 mm	37 mm	41 mm	45 mm	50 mm	56 mm	63 mm	70 mm	75 mm
With Drain	30-CUDo	<b>30-CUD1</b>	<b>30-CUD2</b>	<b>30-CUD3</b>	<b>30-CUD4</b>	<b>30-CUD5</b>	30-CUD6	30-CUD7	30-CUD8	30-CUD9	30-CUD91
Without Drain	30-CUo	30-CU1	30-CU2	30-CU3	30-CU4	30-CU5	30-CU6	30-CU7	30-CU8	30-CU9	30-CU91

Items in **Bold** represent most popular patterns



## Dish

- Used to treat second or third-degree uterine prolapse, cystocele or rectocele
- May be used to treat stress urinary incontinence
- Available with or without supporting membrane

	#0	#1	#2	#3	#4	#5	#6	#7
Size	50 mm	55 mm	60 mm	65 mm	70 mm	75 mm	80 mm	85 mm
With Support	30-DSHS0	30-DSHS1	<b>30-DSHS2</b>	<b>30-DSHS3</b>	<b>30-DSHS4</b>	<b>30-DSHS5</b>	30-DSHS6	30-DSHS7
Without Support	30-DSH0	30-DSH1	30-DSH2	30-DSH3	30-DSH4	30-DSH5	30-DSH6	30-DSH7



## Marland

- Used to treat second or third-degree uterine prolapse or cystocele
- Elevated platform may provide relief of stress urinary incontinence symptoms
- Available with or without supporting membrane

	#2	#3	#4	#5	#6	#7	#8
Size	2.25"	2.50"	2.75"	3.00"	3.25"	3.50"	3.75"
With Support	30-MS2	30-MS3	30-MS4	30-MS5	30-MS6	30-MS7	30-MS8
Without Support	30-M2	30-M3	30-M4	30-M5	30-M6	30-M7	30-M8



## Oval

- Used to treat first or second-degree uterine prolapse or cystocele
- Available with or without supporting membrane

	#1	#2	#3	#4	#5	#6	#7	#8	#9
Size	2.00"	2.25"	2.50"	2.75"	3.00"	3.25"	3.50"	3.75"	4.00"
With Support	30-OVS1	30-OVS2	30-OVS3	30-OVS4	30-OVS5	30-OVS6	30-OVS7	30-OVS8	30-OVS9
Without Support	30-OV1	30-OV2	30-OV3	30-OV4	30-OV5	30-OV6	30-OV7	30-OV8	30-OV9



## Shaatz

- Used to treat first or second degree uterine prolapse
- Convex shape
- Drainage feature

	#0	#1	#2	#3	#4	#5	#6	#7	#8
Size	1.50"	1.75"	2.00"	2.25"	2.50"	2.75"	3.00"	3.25"	3.50"
Shaatz	30-SH0	30-SH1	30-SH2	30-SH3	30-SH4	30-SH5	30-SH6	30-SH7	30-SH8



## Gehrung

- Used to treat second or third-degree uterine prolapse, cystocele or rectocele
- Malleable to allow for individualized fitting . \*Pessary contains metal and must be removed prior to MRI, ultrasound and x-ray procedures.
- With support membrane

	#0	#1	#2	#3	#4	#5	#6	#7	#8	#9
Size	50 mm	55 mm	60 mm	65 mm	70 mm	75 mm	80 mm	85 mm	90 mm	95 mm
Gehrung	30-GH0	30-GH1	<b>30-GH2</b>	<b>30-GH3</b>	30-GH4	30-GH5	30-GH6	30-GH7	30-GH8	30-GH9



## Cup

- Used to treat second or third-degree uterine prolapse, cystocele or rectocele
- May be used to treat stress urinary incontinence
- Available with or without supporting membrane

	#0	#1	#2	#3	#4	#5	#6	#7
Size	50 mm	55 mm	60 mm	65 mm	70 mm	75 mm	80 mm	85 mm
With Support	30-CPS0	30-CPS1	30-CPS2	30-CPS3	30-CPS4	30-CPS5	30-CPS6	30-CPS7
Without Support	30-CP0	30-CP1	30-CP2	30-CP3	30-CP4	30-CP5	30-CP6	30-CP7



## Hodge

- Used to treat first or second-degree uterine prolapse
- Used to relieve the symptoms of retro displacement or an incompetent cervix
- Malleable to allow for individualized fitting . \*Pessary contains metal and must be removed prior to MRI, ultrasound and x-ray procedures.
- Available with or without supporting membrane

	#0	#1	#2	#3	#4	#5	#6	#7	#8	#9
Size	65 mm	70 mm	75 mm	80 mm	85 mm	90 mm	95 mm	100 mm	105 mm	110 mm
With Support	30-HDS0	30-HDS1	30-HDS2	30-HDS3	30-HDS4	30-HDS5	30-HDS6	30-HDS7	30-HDS8	30-HDS9
Without Support	30-HD0	30-HD1	30-HD2	30-HD3	30-HD4	30-HD5	30-HD6	30-HD7	30-HD8	30-HD9

Items in **Bold** represent most popular patterns. \***IMPORTANT Note:** Gehrung and Hodge Pessaries contain metal and must be removed prior to MRI, ultrasound and x-ray procedures. See DFU for important information on pessaries.

## Pessary Fitting Set

- Helps to determine style and size of pessary needed
- Includes six popular sizes of Ring pessary
- Embossed cross-reference chart on top of case
- Autoclavable

Reference	Description
30-FS1000	Pessary Fitting Set



## Pessary Remover

- Stainless steel
- Assists with the removal of any style pessary

Reference	Description
30-900	Pessary Remover



## Pessary Selection Chart

The following chart shows numerous types and degrees of POP (pelvic organ prolapse) in addition to other conditions which may be treatable with the use of a pessary. Some pessaries are designed to provide better symptom relief for certain conditions. This chart is intended to provide guidance for trained medical professionals only.

	Pessary	Prolapse		Cystocele	Rectocele	Stress Incontinence	Retro Displacement	Incompetent Cervix
		(1-2 Degree)	(2-3 Degree)					
	Cube		✓	✓	✓			
	Cup		✓	✓	✓	✓		
	Dish		✓	✓	✓	✓		
	Donut		✓	✓	✓			
	Gellhorn		✓	✓	✓			
	Gellhorn (Short Stem)		✓	✓	✓			
	Gehrung *		✓	✓	✓			
	Hodge *	✓					✓	✓
	Marland		✓	✓		✓		
	Oval	✓		✓				
	Ring	✓						
	Ring w/ Knob	✓				✓		
	Shatz	✓						

**\*IMPORTANT Note:** Gehrung and Hodge Pessaries contain metal and must be removed prior to MRI, ultrasound and x-ray procedures.

HCPCS Code for Pessary ..... A4562  
Code for the Office Visit ..... 99214  
CPT Procedure Code for Fitting & Insertion ..... 57160

## Silicone Vaginal Dilator Sets

- Set includes four sizes of dilators
- Made of medical grade silicone
- Dilators are pliable to facilitate ease of use
- For vaginal use only
- Not available individually, only sold in sets

Reference	Size
30-3001	Large Set*
30-3002	Medium Set*
30-3003	Small Set*

\* Each set includes four sizes, see chart below

Large Set	Length	Small OD	Large OD
30-3001/XL*	150 mm	30 mm	38 mm
30-3001/L*	138 mm	25.5 mm	32 mm
30-3001/M*	128 mm	22.5 mm	28 mm
30-3001/S*	120 mm	19.5 mm	24 mm

Medium Set	Length	Small OD	Large OD
30-3002/XL*	95 mm	26 mm	33 mm
30-3002/L*	86 mm	22.5 mm	28 mm
30-3002/M*	78 mm	19.5 mm	23 mm
30-3002/S*	75 mm	12.5 mm	17 mm

Small Set	Length	Small OD	Large OD
30-3003/XL*	75 mm	13.5 mm	18 mm
30-3003/L*	65 mm	12 mm	16 mm
30-3003/M*	60 mm	10.5 mm	14 mm
30-3003/S*	50 mm	9.5 mm	13 mm

\* Note: For measurement reference only. Product not sold individually.



## Reusable Os Cervical Dilator Set


**Os Locator** (single-ring around handle) - Locates and dilates the external os of the cervix.

**Canal Dilator** (double-ring around handle) - Dilates the length of the cervical canal.

**Fundus Sound** (triple-ring around handle) - Dilates the internal os.

- Flexible and smooth
- Made of PTFE
- Autoclavable

Reference	Description
30-3000	Os Cervical Dilator Set contains one of each:
	Overall Length    Width    Tip Width    Length
Os Locator	215 mm    4 mm    1.8 mm    23 mm
Canal Dilator	215 mm    4 mm    2 mm    45 mm
Fundus Sound	215 mm    4 mm    2 mm    85 mm




## Os Locator, Canal Dilator and Fundus Sound

Cervical dilators are useful in locating and dilating the external os, cervical canal and internal os of the uterus. Dilation is often required prior to endometrial sampling and other diagnostic testing procedures.

- Sterile, single-use
- Flexible and smooth
- Individually packaged


Reference	Description	Length	Quantity
30-5001	Os Locator	8-1/2" (21.6 cm)	25/box
30-5002	Canal Dilator	8-1/2" (21.6 cm)	25/box
30-5003	Fundus Sound	8-1/2" (21.6 cm)	25/box



## Uterine Sound

- Sterile, single-use
- Flexible
- Individually packaged

Reference	Description	Length	Quantity
30-6000	Uterine Sound	9-3/4" (24.7 cm)	25/box




## Endometrial Sampling Set and Endometrial Biopsy Cannula

The Endometrial Sampling Set and Endometrial Biopsy Cannula are used to remove samples from the uterus and mucosal lining of the uterus by scraping or vacuum suction. These devices obtain tissue samples for the purpose of biopsy, pre-cancer screening or menstrual extraction.

### Endometrial Sampling Set

- Sterile, single-use
- 10 ml standard twist-and-lock syringe
- Curette working length = 8-1/8" (20.64 cm)
- Stopette provides protection against perforation

Reference	Description	Quantity
30-3015	Endometrial Sampling Set (includes Curette with a Twist-and-Lock Syringe)	10/box



### Endometrial Biopsy Cannula


#### SoftFlex™ Endometrial Biopsy Cannula

- Sterile, single-use
- Used for routine sampling
- Flexible cannula

#### FirmFlex™ Endometrial Biopsy Cannula


- Sterile, single-use
- Used for sampling of stenosis or atrophic endometrium
- Rigid cannula allows for controlled insertion

Reference	Description	Quantity
30-3010	SoftFlex™ Endometrial Biopsy Cannula	25/box
30-3012	FirmFlex™ Endometrial Biopsy Cannula	25/box



## Standard Colposcopy Set

Reference	Description	Bite Size/Shape	Jaws/Tips	Length
30-15	Graves Vaginal Speculum, Medium	-	-	1-3/8" (3.5 cm) x 4" (10.2 cm)
30-55	Pederson Vaginal Speculum, Narrow Blades, Medium	-	-	7/8" (2.2 cm) x 4" (10.2 cm)
30-951	Emmett Uterine Tenaculum Hook, Style 1	-	sharp, half-curved tip	9" (22.9 cm)
30-1350	Kogan Endocervical Speculum with Gauge & Set Screw	-	fenestrated, 7mm x 25 mm	9-1/2" (24.1 cm)
30-1382	Kevorkian-Younger Endocervical Curette with Basket	-	loop, 3 mm x 12 mm	12" (30.5 cm)
301442WL	Tischler Biopsy Forceps with Lock	7mm x 3mm, oblong	single tooth on both jaws	9-1/4" (23.5 cm)



\* Sets may be customized to meet procedural needs

\* Tolerances in product specifications are subject to minor variations in dimensions and finishes due to the hand craftsmanship of the instruments.



## Delicate Colposcopy Set

Reference	Description	Bite Size/Shape	Jaws/Tips	Length
30-10	Graves Vaginal Speculum, Small	-	-	3/4" (1.9 cm) x 3" (7.6 cm)
30-50	Pederson Vaginal Speculum, Narrow Blades, Small	-	-	1/2" (1.3 cm) x 3" (7.6 cm)
30-951	Emmett Uterine Tenaculum Hook, Style 1	-	sharp, half-curved tip	9" (22.9 cm)
30-1351	Miltex-Kogan Endocervical Speculum with Ratchet & Gauge	-	fenestrated, 3 mm x 25 mm	11" (27.9 cm)
30-1383	Kevorkian-Younger Endocervical Curette without Basket	-	loop, 3 mm x 12 mm	12" (30.5 cm)
301443WL	Burke (Modified Baby Tischler) Biopsy Forceps with Lock	5mm x 3mm, oval	single tooth on both jaws	7-3/4" (19.7 cm)



\* Sets may be customized to meet procedural needs

\* Tolerances in product specifications are subject to minor variations in dimensions and finishes due to the hand craftsmanship of the instruments.

## Instrument Care Kit System

A great way to try various Integra® Miltex® instrument care products!

Reference	Description
3-800	Spray Lube (3-700), Surgical Instrument Cleaner (3-720), Surgical Instrument Stain Remover (3-740), Instrument Care DVD (3-810), Nylon & Stainless Steel Brushes (3-1000, 3-1001)



3-800

## Instrument Care Training

This 25-minute Instructional DVD assists practitioners and staff to properly care for instruments and to keep them in good working condition. The DVD offers step-by-step guidance on handling, cleaning, rinsing, inspecting, lubricating and sterilizing surgical instruments.

The step-by-step poster (SURGHIPPOST) provides easy instructions to clean and prepare instruments for sterilization.

Try a variety of our popular cleaning products with a H.I.P. sample box.

Reference	Description
3-810	Instrument Care DVD
SURGHIPPOST	H.I.P. Poster
SURGHIPBOX	H.I.P. Sample Box



SURGHIPBOX  
(contents)



3-810



SURGHIPPOST

**For a complete listing of our female patient care products, visit [integralife.com/integra-miltex](https://integralife.com/integra-miltex)**

- Biopsy Punch Forceps
- Vaginal and Endocervical Specula
- Obstetrical Forceps
- Gynecological Forceps
- Curettes and Hooks
- Dilators and Sounds
- Scissors and Needle Holders



For more information or to place an order, please contact:  
Integra ■ 589 Davies Drive, York, PA 17402  
866-854-8300 USA ■ 717-840-2763 outside USA ■ 717-840-9347 fax  
[integralife.com/integra-miltex](http://integralife.com/integra-miltex)

