



# mayo stands

## foot-pedal operated

- Base designed to move under low clearance equipment with maximum stability
- Convenient foot-pedal operated lock adjusts height from 36" to 62"
- **co•re casters™** corrosion-resistant, maintenance-free
- All welded stainless steel construction



**P-1065-SS**

- Removable extra large 20" x 25" stainless steel tray for easy cleaning
- Two 2" corrosion-resistant Gemini **co•re casters™**
- 50 lbs. weight capacity evenly distributed



**P-1066-SS**

- Removable 16¼" x 21¼" stainless steel tray for easy cleaning
- Two 2" corrosion-resistant Gemini **co•re casters™**
- 50 lbs. weight capacity evenly distributed



**P-1068-SS**

- Removable 12<sup>5</sup>/<sub>8</sub>" x 19<sup>1</sup>/<sub>8</sub>" stainless steel tray for easy cleaning
- Two 2" corrosion-resistant Gemini **co•re casters™**
- 50 lbs. weight capacity evenly distributed



**P-1069-SS**

- Removable large 16¼" x 21¼" stainless steel tray for easy cleaning
- Four 2" corrosion-resistant Gemini **co•re casters™**
- 50 lbs. weight capacity evenly distributed



**P-1066-A-SS**

- Removable 16 $\frac{1}{4}$ " x 21 $\frac{1}{4}$ " stainless steel tray for easy cleaning
- Adjust height 36" to 60" - friction lock
- Two 2" corrosion-resistant Gemini **co•re casters**<sup>™</sup>
- 50 lbs. weight capacity evenly distributed



**P-1068-A-SS**

- Removable 12 $\frac{5}{8}$ " x 19 $\frac{1}{8}$ " stainless steel tray for easy cleaning
- Adjust height 36" to 60" - friction lock
- Two 2" corrosion-resistant Gemini **co•re casters**<sup>™</sup>
- 50 lbs. weight capacity evenly distributed



# mayo stands

**hand operated**

- Base designed to move under low clearance equipment with maximum stability
- **co•re casters**<sup>™</sup> corrosion-resistant, maintenance-free
- All welded stainless steel construction



**P-1069-A-SS**

- Removable large 16 $\frac{1}{4}$ " x 21 $\frac{1}{4}$ " stainless steel tray for easy cleaning
- Adjust height 36" to 60" - friction lock
- Four 2" corrosion-resistant Gemini **co•re casters**<sup>™</sup>
- 50 lbs. weight capacity evenly distributed



**P-66**

- Removable 19 $\frac{1}{8}$ "L x 12 $\frac{5}{8}$ "W stainless steel tray on powder-coated frame
- Adjust height from 32" to 53"
- Two 3" corrosion-resistant **co•re casters**<sup>™</sup>
- 20 lbs. weight capacity evenly distributed



**P-65 Medical Tray Unit**

- Removable 19 $\frac{1}{8}$ "L x 12 $\frac{5}{8}$ "W stainless steel tray on chrome plated base
- Adjust height from 30" to 40"
- Four 3" corrosion-resistant **co•re casters**<sup>™</sup>
- 40 lbs. weight capacity evenly distributed