



# BD Medical

## Medical Surgical Systems Catalogue Canada



Helping all people  
live healthy lives



Needles  
and  
Syringes



Infusion  
Therapy  
Devices



Pharmacy  
Solutions



Surgical  
Products



Anesthesia  
Products



Sharps  
Disposal  
Products

# Sharps Collectors

## Chemotherapy

Chemotherapy Sharps Collectors feature a wide entry port design for disposal of bulky chemotherapy waste. In addition, these collectors are designed to offer increased leak resistance.

### Features:

- Point-first disposal
- Wide entry port for disposal of bulky chemotherapy waste items
- Exclusive safety plug for temporary closure
- Single-piece moulding—no assembly required
- Sturdy puncture-resistant construction
- Attached locking lid for secure closure and handles for safe transport
- Not autoclaveable

## Chemotherapy Sharps Collectors

PRODUCT #	SIZE (L)	DESCRIPTION	PKG. QTY.
300459	6.6	Plug cap	12
300162	10.3	Large open cap	12
300161	18.7	Large open cap	8



Plug cap



Large open cap



300459



300162



300161



# Sharps Collectors

## Extra-Large

These stackable and sturdy sharps collectors feature large openings for bulky medical waste.

### Features:

- Large access openings for up to 10-inch diameter medical waste
- Temporary and permanent closures
- Sturdy puncture-resistant construction
- Nesting for efficient storage
- Built-in handles for easy handling and disposal

## Extra-Large Yellow Sharps Collectors

PRODUCT #	SIZE (L)	DESCRIPTION	PKG. QTY.
305579	34	Slide top	8
305647	72	Slide top	5
305604	34	Slide top w/ gasket	8
305603	34	Hinge top w/ gasket	8
305614	64	Hinge top w/ gasket	5
305613	72	Slide top w/ gasket	5



## Extra-Large Red Sharps Collectors

PRODUCT #	SIZE (L)	DESCRIPTION	PKG. QTY.
305602	34	Slide top w/ gasket	8
305666	72	Slide top w/ gasket	5
305615	34	Hinge top	8
305610	34	Hinge top w/ gasket	8



# Sharps Collectors

## Phlebotomy

Clear top provides visualization of collector fill levels. Increased fill-level identification helps prevent collector overflow, a leading cause of disposal injuries.

### Features:

- Strong and durable puncture-resistant construction
- Two-way safety slot securely grips needle hub and permits smooth, one-handed needle removal
- Temporary and permanent closures allow for safer transport and reduced risk of tampering and accidental needlesticks
- Larger all-purpose opening permits disposal of bulkier medical wastes, such as swabs and needle holders



367216



**3-port design for use with:**  
**Standard needles**  
**Safety needles**  
**Pen needles**

### Phlebotomy Yellow Sharps Collector

PRODUCT #	SIZE (L)	DESCRIPTION	PKG. QTY.
367216	1.0	Universal unwinder	60

### Phlebotomy Red Sharps Collectors

PRODUCT #	SIZE (L)	DESCRIPTION	PKG. QTY.
305487	1.4	Small open top	36
305635	1.0	Universal unwinder	60



305487



305635



# Sharps Collectors

## Point-of-Use

One-hand, one-step technique, with a counterbalanced door. Container provides security against tampering.<sup>1</sup>

### Features:

- Counterbalanced door returns automatically for one-hand, one-step disposal
- Counterbalanced door stays shut to deter overfilling
- Safety transport handle
- Sturdy puncture-resistant construction

## Point-of-Use Sharps Collectors

PRODUCT #	SIZE (L)	DESCRIPTION	PKG. QTY.
300475	5.1	Side entry	12
300974	5.1	Horizontal entry	20
303050	7.6	Horizontal entry	10
303051	11.3	Horizontal entry	10
305517	5.1	Red Horizontal entry	20



<sup>1</sup> Selecting, Evaluating, and Using Sharps Disposal Containers. CDC/NIOSH Publication 1997-111.



# Sharps Collectors

## One-Piece

With convenient, vertical “point-first” drop, these robust sharps collectors are preassembled and offer a patented one-way funnel valve to minimize needlesticks and overfilling.

### Features:

- Point-first disposal
- One-way tamper-resistant entry
- Secure cap closure
- Single-piece moulding—no assembly required
- Sturdy puncture-resistant construction
- Attached locking lid for secure closure and safe transport

## One-Piece Sharps Collectors

PRODUCT #	SIZE (L)	DESCRIPTION	PKG. QTY.
300460	1.4	Funnel cap	36
300462	1.4	Plug cap	36
300466	3.1	Funnel cap	24
300452	10.3	Large open cap	12



**Funnel cap**



**Plug cap**



**Large open cap**



300460



300462



300466



300452



# Sharps Collectors

## Nestable

With an exclusive funnel-shaped opening and safety neck, these sharps collectors are nestable for easy storage. They allow point-first drop and reduce the risk of tampering and spilling.

### Features:

- Clear top helps deter overfilling
- Point-first disposal
- One-way tamper resistant entry and safety neck
- Nestable for efficient storage
- Sturdy carrying handles
- Sturdy puncture-resistant construction
- Secure cap closure for safe transport

## Nestable Yellow Sharps Collectors

PRODUCT #	SIZE (L)	DESCRIPTION	PKG. QTY.
300450	3.1	Funnel cap	36
300440	7.6	Funnel cap	24
300443	13.2	Funnel cap	20
300434	22.7	Open cap	12
300483	22.7	X-large funnel cap	12



Funnel cap



Open cap



X-large funnel cap



300450

300440

300443

300483

300434



305457



305343

## Nestable Red Sharps Collectors

PRODUCT #	SIZE (L)	DESCRIPTION	PKG. QTY.
305343	7.6	Open cap	24
305457	22.7	Open cap	12



# Sharps Collectors

## Pharmaceutical

BD Pharmaceutical Collectors feature a variety of lid designs and sizes to collect and dispose of unused, expired or other pharmaceutical waste.

### Features:

- Gaskets or liquid-absorbing materials to prevent leakage
- Visual, audible, tactile (VAT) perimeter locks to ensure that the lid is securely attached to the base and to add security against tampering
- 5.1L and 11.3L Pharmaceutical Collectors with counterbalanced door:
  - Facilitate hands-free disposal and are designed to reduce the potential for sharps injuries
  - Designed to improve tamper resistance by preventing hand access to the contents
  - Compatible with locking cabinets to provide security and prevent unwanted access to collector contents

## Pharmaceutical Sharps Collectors

PRODUCT #	SIZE (L)	DESCRIPTION	PKG. QTY.
305058	5.1	Counterbalanced door	20
305622	11.3	Counterbalanced door	10
305633	11.3	Hinge cap	12



305633

305058

305622





## Sharps Collectors Accessories

FOR USE WITH PRODUCT #	ACCESSORY #	DESCRIPTION	PKG. QTY.
300161	305495	Security Bracket	12
	305497	Quick Release Strap	12
300162	305495	Security Bracket	12
	305497	Quick Release Strap	12
300434	305485	Nestable Bracket	12
300440	305485	Nestable Bracket	12
	305570	Stability Bracket	5
	305018	Wire Bracket	6
300443	305485	Nestable Bracket	12
300450	305475	3.2 qt. Cabinet	12
300452	305495	Security Bracket	12
	305497	Quick Release Strap	12
300459	305496	Security Bracket	12
300460	305451	Cup Bracket	12
300462	305451	Cup Bracket	12
300466	305495	Security Bracket	12
	305497	Quick Release Strap	12
300475	305447	Locking Bracket	10
	305021	Stability Bracket	5
300483	305485	Nestable Bracket	12
300974	305447	Locking Bracket	10
	305021	Stability Bracket	5
	305017	Cabinet Beige	1
303050	305552	Stability Bracket	5
	305409	Locking Bracket	10
	305588	Wire Bracket	6
	305420	Cabinet Beige	1
303051	305552	Stability Bracket	5
	305409	Locking Bracket	10
	305588	Wire Bracket	6
	305420	Cabinet Beige	1
305343	305485	Nestable Bracket	12
	305570	Stability Bracket	5
	305018	Wire Bracket	6
305457	305485	Nestable Bracket	12
305487	305970	Cup Bracket	12
305579	305090	Trolley Basic	1
	305091	Trolley, Foot Operated	1
305602	305090	Trolley Basic	1
	305091	Trolley, Foot Operated	1
305635	305970	Cup Bracket	12
305647	305092	Trolley Basic	1
	305093	Trolley, Foot Operated	1
305648	305970	Cup Bracket	12
305666	305092	Trolley Basic	1
	305093	Trolley, Foot Operated	1
367216	305970	Cup Bracket	12



305447



305093



Interlink® is a trademark of Baxter. © 2015

LifeShield® is a trademark of Abbott. © 2015

Safeline® is a trademark of Braun. © 2015

Teflon® is a trademark of E.I. duPont de Nemours and Company. © 2015

BD, BD Logo and all other trademarks are property of Becton, Dickinson and Company. © 2015 BD  
BDM0150-CE Rev1 04/15



2100 Derry Road W., Suite 100  
Mississauga, Ontario L5N 0B3  
[www.bd.com/ca](http://www.bd.com/ca)