

# Patient Interfaces

## Roscoe Nasal Mask

The Roscoe Medical Nasal CPAP Mask with Comfort Cushion™ is extremely lightweight and comes complete with headgear. The headgear buckle has a quick-snap design making it easy to fasten and unfasten. The Mask is constructed of a one piece, flexible silicone cushion that allows for a comfortable secure fit.



Replacement Seal

### Product Specifications

- One piece shaping: flexible, lightweight silicone allows for a comfortable, secure fit
- Headgear with Velcro® fasteners provide maximum versatility and movement
- Unique 360° swivel connector design with standard 22 mm tubing attachment port
- Headgear buckle with quick-snap design making it quick and easy to fasten and unfasten
- Comfortable cushion provides stable, superior and absolute comfort
- Comes complete with sampling port
- Forehead support is flexible allowing for self-adjustment
- Latex-free
- Mask is dishwasher safe for easy cleaning
- 90 day limited warranty
- HCPCS Code - A7034\* (mask)
- HCPCS Code - A7032\* (replacement seal)
- HCPCS Code - A7035\* (replacement headgear)



**Roscoe Medical Offers a 30-Day Money-Back Guarantee on all Roscoe® CPAP Masks, effective from the date of purchase.**

If your patient is not 100% completely satisfied and it is within 30 days from your Roscoe purchase date, we will refund the purchase price to you, no questions asked.

Item #	Description	Qty
CPM-NMS	Nasal Mask with Headgear, small	each
CPM-NMM	Nasal Mask with Headgear, medium	each
CPM-NML	Nasal Mask with Headgear, large	each
90085	Nasal Mask Replacement Seal, small	each
90087	Nasal Mask Replacement Seal, medium	each
90089	Nasal Mask Replacement Seal, large	each
CPH-FF	Replacement Headgear, black	each

\* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.



**To Order:**  
 (800) 3-ROSCOE  
 fax: (440) 572-4261  
 roscoemedical.com