

Cardinal Health™ Nasopharyngeal Airways

A trusted solution for your
airway management needs.

Safe • Reliable • High-quality • Clinically and cost effective




Nasopharyngeal Airways

Features and advantages:

- **Anatomical design with thin wall construction** – Facilitates easy insertion and acceptance of maximum size catheter for safe suctioning and patient comfort.
- **Soft PVC airway with smooth satin finish** – Shields and protects the patient's nasopharyngeal mucosa from adverse effects of repeated passage, as well as minimizing coughing and gagging during suctioning.
- **Labeled with product size/I.D. in Millimeters (mm) and French gauge (Fr)** – For quick and easy size identification.
- **Available in a full range of sizes** – A full range of sizes meets the needs of patients and clinicians.
- **Not made from natural rubber latex**
- **Disposable, single-use**
- **Individually packaged, sterile**

Product specifications and ordering information:

Cat. no.	Size/I.D. (Fr)	Size/I.D. (mm)	Qty.
NA12FR	12	3.0mm	10 ea/bx, 10 bx/cs
NA14FR	14	3.5mm	10 ea/bx, 10 bx/cs
NA16FR	16	4.0mm	10 ea/bx, 10 bx/cs
NA18FR	18	4.5mm	10 ea/bx, 10 bx/cs
NA20FR	20	5.0mm	10 ea/bx, 10 bx/cs
NA22FR	22	5.5mm	10 ea/bx, 10 bx/cs
NA24FR	24	6.0mm	10 ea/bx, 10 bx/cs
NA26FR	26	6.5mm	10 ea/bx, 10 bx/cs
NA28FR	28	7.0mm	10 ea/bx, 10 bx/cs
NA30FR	30	7.5mm	10 ea/bx, 10 bx/cs
NA32FR	32	8.0mm	10 ea/bx, 10 bx/cs
NA34FR	34	8.5mm	10 ea/bx, 10 bx/cs
NA36FR	36	9.0mm	10 ea/bx, 10 bx/cs



Adjustable Nasopharyngeal Airways

Features and advantages:

- Anatomical design with thin wall construction
- Soft PVC airway with smooth satin finish
- Adjustable flange- allows clinician to adjust the flange to the appropriate length
- Labeled with product size/I.D, in millimeters (mm) and French gauge (Fr)
- Available in a full range of sizes
- Not made with natural rubber latex
- Disposable, single-use
- Individually packaged, sterile

Product specifications and ordering information:

Cat. no.	Size/I.D. (Fr)	Size/I.D. (mm)	Qty.
ADJ12FR	12	2.5mm	10 ea/bx, 10 bx/cs
ADJ14FR	14	3.0mm	10 ea/bx, 10 bx/cs
ADJ16FR	16	3.5mm	10 ea/bx, 10 bx/cs
ADJ18F	18	4.0mm	10 ea/bx, 10 bx/cs
ADJ20FR	20	4.5mm	10 ea/bx, 10 bx/cs
ADJ22FR	22	5.0mm	10 ea/bx, 10 bx/cs
ADJ24FR	24	5.5mm	10 ea/bx, 10 bx/cs
ADJ26FR	26	6.0mm	10 ea/bx, 10 bx/cs
ADJ28FR	28	6.5mm	10 ea/bx, 10 bx/cs
ADJ30FR	30	7.0mm	10 ea/bx, 10 bx/cs
ADJ32FR	32	7.5mm	10 ea/bx, 10 bx/cs
ADJ34FR	34	8.0mm	10 ea/bx, 10 bx/cs
ADJ36FR	36	8.5mm	10 ea/bx, 10 bx/cs

To place an order, please contact your Cardinal Health sales representative or customer service at 800.964.5227.



Cardinal Health
General Anesthesia
Waukegan, Illinois

cardinalhealth.com/anesthesia