

Sterile

Sutures

Suture Material:

- Variety of Nylon, Polyglycolic Acid, Chromic Gut, Silk and Polypropylene Sutures
- Easy to handle for smooth tissue passage and minimal tissue reaction
- High tensile strength reduces breakage

Needles:

- Needles crimped to thread using swaging allowing for a smooth transition from needle to thread
- Sharp, fine point needle ensures smooth penetration of tissue
- Most popular needle sizes
- 3/8 circle reverse cut



Black Nylon Sutures



A non-absorbable suture offering high tensile strength. Monofilament for a smooth passage through tissue. Sterile. 12/box.



Black Nylon Monofilament – Non-Absorbable				
Item #	Needle	Needle Description	Suture Size	Suture Length
P420661	FS-2	3/8 Circle Reverse Cut 19mm	5-0	18"
P420662	FS-2	3/8 Circle Reverse Cut 19mm	4-0	18"
P420663	FS-1	3/8 Circle Reverse Cut 24mm	3-0	18"
P420697	P-1	3/8 Circle Reverse Cut 11mm	6-0	18"
P420698	P-3	3/8 Circle Reverse Cut 13mm	5-0	18"
P420699	P-3	3/8 Circle Reverse Cut 13mm	4-0	18"
P421666	PS-2	3/8 Circle Reverse Cut 19mm	5-0	18"
P421667	PS-2	3/8 Circle Reverse Cut 19mm	4-0	18"
P421698	P-3	3/8 Circle Reverse Cut 13mm	6-0	18"

Braided Polyglycolic Acid (PGA) Sutures

An absorbable, undyed suture made from polyglycolic acid. Braided for strength and easy handling and coated to allow for easy passage through tissue. Sterile. 12/box.





Braided, Undyed PGA – Absorbable				
Item #	Needle	Needle Description	Suture Size	Suture Length
P420422	FS-2	3/8 Circle Reverse Cut 19mm	4-0	27"
P420493	P-3	3/8 Circle Reverse Cut 13mm	5-0	18"
P420494	P-3	3/8 Circle Reverse Cut 13mm	4-0	18"
P420496	PS-2	3/8 Circle Reverse Cut 19mm	4-0	18"
P420497	PS-2	3/8 Circle Reverse Cut 19mm	3-0	18"

Chromic Gut Sutures


An absorbable suture made from purified connective tissue derived from cow or sheep intestines. Treated with chromic salt to prolong the absorption process and improve tensile strength. Packaged in alcohol to prevent drying. Sterile. 12/box.





Chromic Gut – Absorbable				
Item #	Needle	Needle Description	Suture Size	Suture Length
P420635		3/8 Circle Reverse Cut 19mm	4-0	27"
P420636		3/8 Circle Reverse Cut 19mm	3-0	27"

Braided Black Silk Sutures


A non-absorbable suture made from natural silk that is degummed to remove extraneous materials. Braided for strength and easy handling. Sterile. 12/box










Braided Black Silk – Non-Absorbable				
Item #	Needle	Needle Description	Suture Size	Suture Length
P420683		3/8 Circle Reverse Cut 19mm	4-0	18"
P420684		3/8 Circle Reverse Cut 24mm	3-0	18"

Polypropylene Sutures

A non-absorbable suture made of polypropylene and dyed blue for enhanced visibility. Monofilament for a smooth passage through tissue. Sterile. 12/box



Blue Polypropylene, Monofilament – Non-Absorbable				
Item #	Needle	Needle Description	Suture Size	Suture Length
P428661		3/8 Circle Reverse Cut 19mm	5-0	18"
P428683		3/8 Circle Reverse Cut 19mm	4-0	18"
P428684		3/8 Circle Reverse Cut 24mm	3-0	18"
P428695		3/8 Circle Reverse Cut 13mm	6-0	18"
P428697		3/8 Circle Reverse Cut 11mm	6-0	18"
P428698		3/8 Circle Reverse Cut 13mm	5-0	18"
P428699		3/8 Circle Reverse Cut 13mm	4-0	18"

Label Description

P420662

Color Code

Item Code

Pro ADVANTAGE[®]
by **NDC**

**MONOFILAMENT
NON-ABSORBABLE**

Suture Material — **BLACK NYLON**

Suture Size — **4 – 0** Metric 1.5

Suture Length — **18" (25 cm)**

19mm 3/8 Circle
Reverse Cutting

Needle

FS-2

STERILE | **EO**

Distributed by: