



## B. Braun Needles and Syringes

Consistent, dependable performance in a wide range of options

### Needles

#### Optimized Patient Comfort and Convenience

- Extremely smooth surface with light silicone coating to minimize discomfort
- Color-coded hubs for easy identification
- Wide range of diameters and lengths



### Syringes

#### Consistent Performance

- Smooth plunger action and double-seal stopper
- High visibility – transparent barrel with bold graduations
- Secure backstop to prevent accidental plunger withdrawal
- Double seal stopper to prevent medication leakage
- Large, stable finger grips
- DEHP-free and latex-free

#### Product Configurations

- 1mL to 100mL
- With and without graduations
- Standard, fine dose, insulin and catheter-fitting syringes
- Luer lock, Luer slip and Eccentric tips

#### Product Materials

- Barrel and plunger rod: Polypropylene
- Plunger tip: Polyisoprene
- DEHP-free and latex-free

# Ordering Information

## Needles

REF Number	Description	Units Per Case
4665118-02	18G x 1 in.	1000
4665120-02	18G x 1½ in.	1000
4667123-02	18G x 2 in.	1000
4665317-02	19G x 1¼ in.	1000
4657799-02	19G x 1½ in. RB	1000
4665112-02	19G x 1½ in. SB	1000
4665511-02	19G x 2 in.	1000
4657500-02	20G x 1 in.	1000
4657519-02	20G x 1½ in.	1000
4667093-02	20G x 2 in.	1000
4665791-02	20G x 2¾ in.	1000
4657543-02	21G x 1 in.	1000
4657527-02	21G x 1½ in.	1000
4665503-02	21G x 2 in.	1000
4665465-02	21G x 3⅛ in.	1000
4665643-02	21G x 4¾ in.	1000
4657624-02	22G x 1¼ in.	1000
4660021-02	22G x 1½ in.	1000
4650018-02	22G x 2 in.	1000
4657659-02	22G x 1 in.	1000
4657667-02	23G x 1 in.	1000
4657640-02	23G x 1¼ in.	1000
4665600-02	23G x 2⅜ in.	1000
4665635-02	23G x 3⅛ in.	1000
4657675-02	24G x 1 in.	1000
4657853-02	25G x 5⁄8 in.	1000
9186158-02	25G x 1 in.	1000
9186166-02	25G x 1½ in.	1000
4657683-02	26G x 1 in.	1000
4665457-02	26G x ½ in.	1000
4657705-02	27G x ¾ in.	1000
4665406-02	27G x ½ in.	1000
9186174-02	27G x 1 in.	1000
9186182-02	27G x 1½ in.	1000
9180109-02	21G x 7⁄8 in. Blunt Cannula	1000
9180117-02	27G x 1 in. Blunt Cannula	1000

## Tuberculin Syringes

9161406V-02	1mL TB LS	1800
9161422-02	1mL TB 25G x 5⁄8 in. DET	1600
9161414-02	1mL TB 26G x 1½ in. DET	1600
9161502SUS	1mL TB 30G x ½ in. PN	2400

## Insulin Syringes

9151117SUS	0.5mL Insulin 30G x 5⁄16 in. PN	2400
9151125SUS	0.5mL Insulin 30G x ½ in. PN	2400
9151133SUS	1mL Insulin 30G x 5⁄16 in. PN	2400
9151141SUS	1mL Insulin 30G x ½ in. PN	2400
9161732VUS	1mL Insulin 27G x ½ in. DET	1600
9161708VUS	1mL Insulin	1800

For more information please contact us at Customer Support in the US at **1-800-848-2066**.

Visit our website at [www.bbraunusa.com](http://www.bbraunusa.com)

## Syringes

REF Number	Description	Units Per Case
4610303-02	3mL LL	2400
4610302-02	3mL LS	2400
4610304-02	3mL LL 18G x 1 in.	1600
4610305-02	3mL LL 18G x 1½ in.	1600
4610306-02	3mL LL 20G x 1 in.	1600
4610307-02	3mL LL 20G x 1½ in.	1600
4610308-02	3mL LL 21G x1 in.	1600
4610309-02	3mL LL 21G x 1½ in.	1600
4610310-02	3mL LL 22G x 1 in.	1600
4610311-02	3mL LL 22G x 1½ in.	1600
4610312-02	3mL LL 23G x 1 in.	1600
4610313-02	3mL LL 23G x 1¼ in.	1600
4610314-02	3mL LL 25G x 5⁄8 in.	1600
4610315-02	3mL LL 25G x 1 in.	1600
4610316-02	3mL LL 27G x 1½ in.	1600
4610518-02	5mL LS	2000
4617053V-02	5mL LL	2000
4616057V-02	5mL Eccentric	2000
4610519-02	5mL LL 18G x 1½ in.	900
4610520-02	5mL LL 20G x 1 in.	900
4610521-02	5mL LL 20G x 1½ in.	900
4610522-02	5mL LL 21G x 1 in.	900
4610523-02	5mL LL 21G x 1½ in.	900
4610524-02	5mL LL 22G x 1 in.	900
4610525-02	5mL LL 22G x 1¼ in.	900
4610526-02	5mL LL 22G x 1½ in.	900
4610527-02	5mL LL 23G x 1 in.	900
4610528-02	5mL LL 25G x 1½ in.	900
4610529-02	5mL LL 27G x 1½ in.	900
4611030-02	10mL LS	1200
4617100V-02	10mL LL	1200
4616103V-02	10mL Eccentric	1200
4611031-02	10mL LL 18G x 1 in.	900
4611033-02	10mL LL 18G x 1½ in.	900
4611034-02	10mL LL 20G x 1½ in.	900
4611035-02	10mL LL 21G x1 in.	900
4611036-02	10mL LL 21G x 1½ in.	900
4611037-02	10mL LL 22G x 1 in.	900
4611038-02	10mL LL 22G x 1½ in.	900
4611039-02	10mL LL 23G x 1 in.	900
4611040-02	10mL LL 25G x 1½ in.	900
4612041-02	20mL LS	800
4617207V-02	20mL LL	800
4616200V-02	20mL Eccentric	800
4612042-02	20mL LL 21G x 1½ in.	800
4613043-02	30mL LS	600
4617304F-02	30mL LL	600
4616308F-02	30mL Eccentric	600
4615044F-02	50mL LS	50
4617509F-02	50mL LL	50
4616502F-02	50mL Eccentric	50
4613503F-02	50mL CATHETER	50
4613554F-02	50mL CATHETER with Luer Adaptor	50
4614003F-02	100mL Syringe with Luer Adaptor	50