



Neocussor™

Technical Specifications



Output Speed -	40 cycles per second
Power Source -	1.5 volt DC
Speed Control -	Single Switch On/Off
Housing Material -	ABS Cicolac
Configurations -	Hand-held
Unit Dimensions -	6" H x 1 ¹ / ₄ " W x 1 ¹ / ₄ " D
Unit Weight -	6 oz.
Shipping Dimensions -	10" x 7" x 5"
Shipping Weight -	1 Lbs.



GENERAL PHYSIOTHERAPY, INC.
13222 LAKEFRONT DRIVE - EARTH CITY (ST. LOUIS), MISSOURI 63045-1504 - USA
PHONE: (314) 291-1442 TOLL FREE: (800) 237-1832 FAX: (314) 291-1485
www.g5.com

