Vinyl Powder Free Examination Glove

The OmniTrust #312 Series Vinyl Powder Free Examination Glove is made from the highest quality polyvinyl chloride material and formulation available. Our vinyl gloves provide you with the best possible value because they're made on the industry's most efficient automated manufacturing lines and tested to the most stringent international quality standards. All of our gloves are produced with DOTP, a safe alternative to the harmful plasticizers found in other vinyl gloves. All of this adds up to a superior vinyl glove in quality, comfort and performance:



- California Proposition 65 Compliant
- Ambidextrous, Beaded Cuff, Single Use Only
- Smooth Finish for Improved Feel / Tactile Sensitivity
- 100% Synthetic Vinyl
- Not Made With Natural Rubber Latex

Tensile Strength - Aged/Unaged	Elongation - Aged/Unaged
min. 11 MPa / min. 11 MPa	min. 350% / min. 350%

Dimensions	Finger	Palm	Cuff	Length	Width
Small	min. 0.08 mm	min. 0.08 mm	min. 0.05 mm	min. 240 mm	85 ± 5 mm
Medium	min. 0.08 mm	min. 0.08 mm	min. 0.05 mm	min. 240 mm	95 ± 5 mm
Large	min. 0.08 mm	min. 0.08 mm	min. 0.05 mm	min. 240 mm	105 ± 5 mm
Extra Large	min. 0.08 mm	min. 0.08 mm	min. 0.05 mm	min. 240 mm	115 ± 5 mm

Item Numbers	Small	Medium	Large	X-Large
312 Series	312-11	312-12	312-13	312-14

Packaging	Boxes	Cases
packaged by weight	100 Gloves	10 Boxes



