



## Overbed Tables

### Model 125



### Features:

- Melamine laminate top
- Flush mounted vinyl t-mold edge
- Base: Spring assisted lift mechanism
- Infinite height adjustment between 28" - 45"
- Slight upward pressure allows the table to be adjusted without requiring the handle to be activated
- Available in both U and H style base design
- 24" opening
- Standard Chrome plated steel column & base
- Standard Gunstock Walnut finish
- Optional dual mirror vanity

### Specifications:

- Top: 15" x 30"
- Height: 28" - 45"
- 1.5" Casters
- Weight capacity: 100 lbs evenly distributed
- Warranty: 1 Year
- **129 same table as shown with vanity**

## Finish Options



Bannister  
Oak

Clear Maple

English Oak

Figured  
Mahogany

Fusion  
Maple



Gunstock  
Walnut

Nepal Teak

Quartered  
Oak

Solar Oak

Wild Cherry

### Standard Edge Choices



Brown

Opal

Black

### Standard Base Finish Choices



Brown

Gray

Opal

Black

Chrome

\*Additional costs may apply to finish and trim upgrades