


Splinting > Soft Strapping Materials



Rolyan® SoftStrap® Strapping Material. Similar to Beta Pile II™. Soft, nonstretchable 1/4" thick (6.4mm) strapping material with sealed side edges for increased durability. Composed of a 5-layer laminate of hook-compatible nylon fabric, open-cell polyurethane foam and a central stabilizing layer with a self-sealing edge. Mates with standard and HTH hook. Conforms anatomically without over-tightening. Recommended where edema or hypersensitive skin is a problem. Latex free.

Item No.	Dimensions	Price
551798	1/2" x 25 yd. (1.3cm x 23m) Beige	\$21.00
A34960	1" x 10 yd. (2.5cm x 9m) Beige	\$17.85
6993	1 1/2" x 10 yd. (3.8cm x 9m) Beige	\$27.85
A34961	2" x 10 yd. (5cm x 9m) Beige	\$29.10
6991	6" x 3.3 yd. (15.2cm x 3m) Beige	\$27.15
7100260	1" X 10 yd. (2.5 cm x 9 m) Black	\$16.42
7100261	1" X 10 yd. (2.5 cm x 9 m) Blue	\$16.42
7100262	2" X 10 yd. (5 cm x 9 m) Black	\$26.77
7100263	2" X 10 yd. (5 cm x 9 m) Blue	\$26.77



Splinting



Rolyan® Economy SoftStrap Strapping Material. Soft, nonstretchable 1/4" thick (6.4mm) strapping material. Composed of a five-layer laminate of hook-compatible nylon fabric, open-cell polyurethane foam and a central stabilizing layer. Mates with standard and HTH hook. Conforms anatomically without over-tightening. Latex free.



Beige	Strapping Material	Color	Price
081515238	1" x 10 yd. (2.5cm x 9m)	Beige	\$15.10
081515246	1 1/2" x 10 yd. (3.8cm x 9m)	Beige	\$25.20
081515253	2" x 10 yd. (5cm x 9m)	Beige	\$26.25
081571116	1" x 10 yd. (2.5cm x 9m)	Black	\$15.10
081571124	2" x 10 yd. (5cm x 9m)	Black	\$26.25



Rolyan® Economy Super Strap II Strapping Material. Recommended for strapping requiring strong fastening strength and longer use. 1/4" (6.4mm) thick breathable foam and soft tricot backing. Non-stretchable and hypoallergenic. Washable. Both sides mate with hook materials. Latex free.



Beige	Dimensions	Price
081515261	1" x 10 yd. (2.5cm x 9m)	\$23.20
081515279	2" x 10 yd. (5cm x 9m)	\$38.35



Rolyan® Super Strap® II Strapping Material. Recommended for strapping requiring strong fastening strength and longer use. Consists of non-fraying knitted loop fabric, 1/4" (6.4mm) thick breathable foam and soft tricot backing. Non-stretchable and hypoallergenic. Washable. Both sides mate with hook materials. Latex free. Strapping material sold in a bag.

755701 1" x 10 yd. (2.5cm x 9.14m)
\$24.15

755702 2" x 10 yd. (5cm x 9.14m)
\$41.50

755704 4" x 10 yd. (10.16cm x 9.14m)
\$81.10

546

Call toll-free 1.800.323.5547

fax 1.800.547.4333

Hook and loop in one product



Rolyan® Back-to-Back Hook & Loop. White hook and royal blue loop permanently bonded back-to-back. Excellent for quick circumferential applications to meet positioning, restraining and strapping needs. Sold in a bag. 10 (9 m) or 25 (22.8 m) yd. length. Latex free.



TIP: Reserve one pair of scissors exclusively for cutting self-adhesive strapping and padding materials. The non-stick scissors on page 560 are especially easy to clean.

Item No.	Size	Price
566114	3/8" x 10 yd. (1.58 cm x 9 m)	\$14.35
A28690	1" x 10 yd. (2.5 cm x 9 m)	\$37.54
566116	1" x 25 yd. (5 cm x 22.8 m)	\$40.21
A28691	2" x 10 yd. (5 cm x 9 m)	\$47.31
566118	2" x 25 yd. (5 cm x 22.8 m)	\$72.19



VELCRO® Brand Translucent Self-Adhesive Low Profile HTH Hook. Flexible translucent plastic hook doesn't attract lint and seldom snags clothing. Mates with all loop and soft strapping materials. Translucent material blends with any color of splint or strapping material. Sold in a box. Latex free.

Item No.	Width	Length	Price
A34024	1" (2.54 cm)	25 yds. (22.86 m)	\$58.30
A34025	2" (5 cm)	25 yds. (22.86 m)	\$115.95
A34022	1" (2.54 cm)	10 yds. (9.14 m)	\$35.15
A34023	2" (5 cm)	10 yds. (9.14 m)	\$60.55

VELCRO® Brand Back-to-Back Hook & Loop. Continuous roll is permanently bonded back-to-back, and accommodates a wide variety of applications. Velcro can be cut to size, reused, replaced and cut again to a smaller size. All of these features help you reduce strap inventory.



The hook & loop attaches to itself without adhesive, so it's not messy, not painful to patients when removed and does not leave a sticky residue. Soft edges prevent irritation and discomfort for the patient. This versatile tape is ideal for quick circumferential applications to meet positioning, restraining and strapping needs. Packaged in a clear, sealed polybag.

Item No.	Description	Dimensions	Price
081610708	Narrow	3/4" x 10 yd.	\$28.65
081610716	Regular	1 1/2" x 10 yd.	\$45.60
081610724	Wide	3" x 10 yd.	\$85.70
081610732	X-Wide	5" x 10 yd.	\$126.75

Dual Lock® Fastening System. Locks securely to itself for high shear strength application. Mates with standard loop for high-retention attachment. Self-adhesive. Similar to Poly-Lock®. Latex free.
081592799
1" x 10 yd.
(2.5 cm x 9 m)
\$99.95



Strong holding power!

SUCCESSFUL SPLINTING

HANDS-ON FABRICATION AND TROUBLESHOOTING

6 CONTACT HOURS/6 CEU

WHO SHOULD ATTEND
PT, PTA, OT, OTA

EDUCATIONAL LEVEL
Intermediate

PREREQUISITES
None

To view current PHA courses, visit PHA.credit

Splinting

Order online. It's click & easy: www.performancehealth.com