

## Pivot & Tilt Overbed Table

HCPCS: E0274



FIGURE A

- The mast of the table pivots and can be locked in one of the three positions from flat on the floor to 90° (Figure A)
- New and improved larger tubing
- The pivot feature allows the table top to be positioned closer to individuals in a bed, sitting in a wheelchair, or folded for storage
- The tabletop tilts with the angle of the mast to provide a flat surface regardless of the angle of the mast
- Four 2" casters (two locking)
- Only wheels have to be installed

PRODUCTS	Description
13000	Overbed Table,Pivot/Tilt,1/cs
SPECIFICATIONS	Description
Carton Shipping Weight	22 lbs.
Overall Dimensions at Base	25" (W) x 28" (D)
Tabletop Adjustment	28" - 37"
Tabletop Surface Dimensions	15.75" (D) x 23.5" (W)
Warranty	Limited Lifetime
Weight	21 lbs.
Weight Capacity	40 lbs