




Regu
PROCARE®

XTEND 174 cervical collar
Proven Design. Improved Comfort.




The ProCare® XTEND 174 cervical collar helps provide substantial motion restriction while helping to minimize painful pressure points that can lead to skin maceration and/or breakdown. Removable, open cell foam padding wicks moisture away from the skin to help promote proper patient hygiene and enhanced comfort while minimizing skin irritation. X-ray, CT, and MRI lucent.

Innovative design features that help encourage patient compliance

- Helps restrict range-of-motion
- Helps reduce Occipital Pressure Ulcers
- Aids in proper fit and optimum immobilization*

Enhanced patient comfort using soft, breathable foam padding

- Anti-Microbial
- Designed to help support appropriate skin care and aid in proper pressure distribution*
- Fabric performs at highest rating for moisture management properties*



Pediatric cervical collar

- Same innovative design as adult sizing
- Kid-friendly design elements help make wearing a cervical collar less intimidating





XTEND 174 cervical collar

Proven Design. Improved Comfort.

PART NUMBER	SIZE	HEIGHT*	NECK CIRCUMFERENCE	AGE
ADULT COLLARS ONLY				
79-83201	SUPER SHORT	1"	10" - 20"	
79-83922	XSMALL	1.25"	8" - 14"	
79-83203	SHORT	2"	10" - 20"	
79-83205	REGULAR	2.5"	10" - 20"	
79-83207	TALL	3"	10" - 20"	
79-83208	STOUT	1.75"	14" - 24"	
ADULT COLLARS WITH REPLACEMENT PADS				
79-83221	SUPER SHORT W/REPL. PAD	1"	10" - 20"	
79-83232	XSMALL W/REPL. PAD	1.25"	8" - 14"	
79-83233	SHORT W/REPL. PAD	2"	10" - 20"	
79-83235	REGULAR W/REPL. PAD	2.5"	10" - 20"	
79-83237	TALL W/REPL. PAD	3"	10" - 20"	
79-83238	STOUT W/REPL. PAD	1.75"	14" - 24"	
ADULT REPLACEMENT PAD SETS				
79-83942	XSMALL REPL. PAD SET			
79-83224	STANDARD REPL. PAD SET			
79-83225	STOUT REPL. PAD SET			
PEDIATRIC COLLARS ONLY				
79-83281	PED 0			0 - 6mo.
79-83251	PED 1	1.125"	≤9"	6mo. - 2yrs.
79-83252	PED 2	1.25"	8" - 11"	2 - 6yrs.
79-83202	PED 3	1.75"	9" - 13"	6 - 12yrs.
PEDIATRIC COLLARS WITH REPLACEMENT PADS				
79-83206	PED 0 W/REPL. PAD			0 - 6mo.
79-83227	PED 1 W/REPL. PAD	1.125"	≤9"	6mo. - 2yrs.
79-83228	PED 2 W/REPL. PAD	1.25"	8" - 11"	2 - 6yrs.
79-83222	PED 3 W/REPL. PAD	1.75"	9" - 13"	6 - 12yrs.
PEDIATRIC REPLACEMENT PAD SETS				
79-83972	PED 0 REPL. PAD SET			0 - 6mo.
79-83253	PED 1 REPL. PAD SET	1.125"	≤9"	6mo. - 2yrs.
79-83254	PED 2 REPL. PAD SET	1.25"	8" - 11"	2 - 6yrs.
79-83204	PED 3 REPL. PAD SET	1.75"	9" - 13"	6 - 12yrs.

*Height - measurement from sternal notch to chin



T 800.526.8785 D 760.727.1280 F 800.457.4221
2900 Lake Vista Drive | Dallas, TX 75067 | U.S.A.
DJOGlobal.com

* Data on file.

Individual results may vary. Neither DJO Global, Inc. nor any of its subsidiaries dispense medical advice. The contents of this sheet do not constitute medical, legal, or any other type of professional advice. Rather, please consult your healthcare professional for information on the courses of treatment, if any, which may be appropriate for you.