

SHARPS CONTAINERS

by **P** LASTI – **P** RODUCTS, INC.

Made in the U.S.A.

Because You Can Never Be Too Safe

WALL MOUNTED SHARPS DISPOSAL SYSTEM

The 5 qt. Sharps Container is designed to be used in the Plasti-Products Wall Cabinet or the Kendall®/Covidien®/Cardinal Health™ 5 qt. Wall Cabinet.



Touchless Deposit of Sharps

The automatic rotating chamber is designed to safely deposit sharps into the puncture resistant container. If necessary, sharps can be dropped manually by using the actuating lever.

- Optional glove box may be attached to the wall cabinet
- Additional glove boxes may be attached to the wall



141020



145002

The 5.4 qt. Sharps Container is designed to be used in the Plasti-Products Wall Bracket 149154 or BD® 5.4 qt. Wall Cabinet or with the BD® 5.4 qt. Locking Wall Bracket. Container has an automatic rotating chamber for touchless deposit of most sharps.



149154 Locking Wall Bracket for 5.4 qt.



143154 5.4 qt. Sharps Container

The Locking Wall Bracket is designed to be used with the 5.4 qt. Sharps Container

HORIZONTAL ENTRY SHARPS CONTAINERS



Horizontal entry sharps containers operate the same as the 5 qt. wall mounted system for touchless deposit of sharps. The automatic rotating chamber is designed to safely deposit sharps into the puncture resistant container. If necessary, sharps can be dropped manually by using the actuating lever.

Other Features

- Horizontal drop maximizes the use of the container
- Available in sizes small (4 qt.), medium (8 qt.) and large (14 qt.)
- All sizes may be mounted to the wall using optional wall mounting bracket
- 4 qt. and 8 qt. containers are also ideal for countertop use

BIG MOUTH SHARPS CONTAINERS

These containers feature a large opening to accommodate bulky sharps and contaminated materials.








The Big Mouth Containers should never be used in areas accessible to patients or to the public

Other Features

- Large screw cap closure with locking device
- Available in sizes medium (8 qt.) and large (14 qt.)
- Both sizes may be mounted to the wall using optional wall mounting bracket



Ordering Information

PRODUCT NUMBER	AMOUNT PER CASE	DESCRIPTION	PRODUCT NUMBER	AMOUNT PER CASE	DESCRIPTION
	2	Locking Wall Cabinet (for 5 qt. Container)		2	Glove Box (with Wall Bracket)
143002	2	Glove Box (for 5 qt. Wall Cabinet)	148002	2	Glove Box (with Wall Bracket)
144002	1	Combo (Wall Cabinet, Glove Box, Red Container)	145004	25	4 qt. Horizontal Entry Sharps Container, Red
145002	2	Glove Box (with Wall Bracket)	145008	20	8 qt. Horizontal Entry Sharps Container, Red
148002	2	Glove Box (with Wall Bracket)	145014	10	14 qt. Horizontal Entry Sharps Container, Red
	20	5 qt. Sharps Container, Red		20	8 qt. Big Mouth Sharps Container, Red
141020	20	5 qt. Sharps Container, Clear	146008	10	14 qt. Big Mouth Sharps Container, Red
142020	20	5 qt. Sharps Container, Clear packed 10 per dispenser box	146014	10	14 qt. Big Mouth Sharps Container, Red
	20	5.4 qt. Sharps Container, Red		5	Wall Mounting Bracket
143154	20	5.4 qt. Sharps Container, Clear	147005	5	Wall Mounting Bracket
143254	20	5.4 qt. Sharps Container, Clear packed 10 per dispenser box	100020	10	Security Lock
	2	Locking Wall Bracket (for 5.4 qt. Container)			for 4, 8, 14 qt. Containers
148154	2	Locking Wall Bracket (for 5.4 qt. Container)			

PLASTI-**P**RODUCTS, INC.

14315 "C" Circle – Omaha, NE 68144 – 800-527-5396

www.plastiproducts.com