

Jackson Original Model – Short Length

- Full 90 degree curve
- ¼" shorter than standard Jackson length
- Complete with inner and outer cannula
- Non-fenestrated
- Tucked in distal end
- Manufactured using Stainless Steel

1395-1#4 1 ea
 1395-1#5 1 ea
 1395-1#6 1 ea
 1395-1#8 1 ea

TUBE DIMENSIONS		
O.D. mm	Actual Size	Size
3		00
4		0
5		1
6		2
7		3
8		4
9		5
10		6
11		7
12		8
13		9
14		10
16		12



TRACHEA AND LARYNGECTOMY TUBE LENGTHS
 LENGTH IN INCHES MEASURED ON OUTSIDE CURVE

JOHN BUNN® ITEM NO.	00	0	1	2	3	4	5	6	7	8	9	10	12
	1 ³ / ₄	2	2 ¹ / ₄	2 ¹ / ₂	2 ³ / ₄	3	3 ³ / ₈	3 ⁵ / ₈	4	4 ¹ / ₄	4 ³ / ₈	4 ³ / ₈	4 ³ / ₈
	1 ¹ / ₂	1 ³ / ₄	2 ¹ / ₈	2 ¹ / ₄	2 ³ / ₈	2 ¹ / ₂	2 ⁵ / ₈	2 ⁷ / ₈	3	3 ¹ / ₈	3 ¹ / ₄	3 ³ / ₈	

