

Everyday Raised Toilet Seat

- Adds 4 ½" to toilet seat height reducing bending effort
- Lightweight and portable
- Maximum Weight Capacity 250 lb - EVENLY DISTRIBUTED
- Designed to meet the requirements of HCPCS code: E0244
- Lifetime limited warranty

6909A 6/cs
6909A-1 1 ea



Specifications

Width Between Arms 18 ¼"
Overall Width 18 ½"
Seat Height 16"
Seat Depth 16"
Weight Capacity 250 lb EVENLY DISTRIBUTED

