

The Everest & Jennings Transport Chair is designed as an affordable patient transfer device. Its compact size and lighter weight allow a companion to lift and store the chair with ease.



7 stylish color options

**Standard Features:**

- Available in 17" and 19" seat widths
- Aluminum frame in stylish color options of Red, Blue, Silver, Green, Pink, Aqua and Black
- Fold-down back makes transport and storage easy
- Black padded nylon upholstery is flame retardant
- Carry straps on seat allow for ease of lifting and folding
- Full length padded arms provide user comfort
- Tool-free adjustable detachable swingaway footrests
- Comes standard with 3" webbing black nylon leg strap with velcro adjustable loops to slide over installed footrests
- Black padded nylon leg strap fits 17" and 19" seat width transport chair frames
- Attendant wheel locks on both rear casters for better control and access by companion
- Weighs 20 lb (w/o front rigging)
- Maximum Weight Capacity 250 lb - EVENLY DISTRIBUTED
- Designed to meet the requirements of Medicare code E1038

**Warranty:**

One (1) year limited warranty



<b>Aluminum Transport Chair</b>		<i>Dimensions and Weights</i>									
		Seat			Back	Overall			Weight		
		Width	Depth	Height	Height	Depth	Height	Width	Folded Width	Without front rigging	Max Capacity
<b>EJ781-1</b>	Lightweight Aluminum - Red	17"	16"	19 ¼"	18 ½"	23"	39 ¼"	22 ¼"	10 ¾"	20 lb	250 lb
<b>EJ782-1</b>	Lightweight Aluminum - Red	19"	16"	19 ¼"	18 ½"	23"	39 ¼"	24 ¼"	10 ¾"	20 lb	250 lb
<b>EJ783-1</b>	Lightweight Aluminum - Silver	17"	16"	19 ¼"	18 ½"	23"	39 ¼"	22 ¼"	10 ¾"	20 lb	250 lb
<b>EJ784-1</b>	Lightweight Aluminum - Silver	19"	16"	19 ¼"	18 ½"	23"	39 ¼"	24 ¼"	10 ¾"	20 lb	250 lb
<b>EJ787-1</b>	Lightweight Aluminum - Blue	17"	16"	19 ¼"	18 ½"	23"	39 ¼"	22 ¼"	10 ¾"	20 lb	250 lb
<b>EJ788-1</b>	Lightweight Aluminum - Blue	19"	16"	19 ¼"	18 ½"	23"	39 ¼"	24 ¼"	10 ¾"	20 lb	250 lb
<b>EJ789-1</b>	Lightweight Aluminum - Green	17"	16"	19 ¼"	18 ½"	23"	39 ¼"	22 ¼"	10 ¾"	20 lb	250 lb
<b>EJ790-1</b>	Lightweight Aluminum - Pink	17"	16"	19 ¼"	18 ½"	23"	39 ¼"	22 ¼"	10 ¾"	20 lb	250 lb
<b>EJ791-1</b>	Lightweight Aluminum - Pink	19"	16"	19 ¼"	18 ½"	23"	39 ¼"	24 ¼"	10 ¾"	20 lb	250 lb
<b>EJ793-1</b>	Lightweight Aluminum - Aqua	19"	16"	19 ¼"	18 ½"	23"	39 ¼"	24 ¼"	10 ¾"	20 lb	250 lb
<b>EJ785-1</b>	Lightweight Aluminum - Black	17"	16"	19 ¼"	18 ½"	23"	39 ¼"	22 ¼"	10 ¾"	20 lb	250 lb
<b>EJ786-1</b>	Lightweight Aluminum - Black	19"	16"	19 ¼"	18 ½"	23"	39 ¼"	24 ¼"	10 ¾"	20 lb	250 lb

All dimensions are ±1/4"  
\*Without Front Rigging.

<b>Common Accessories &amp; Parts</b>		<b>HCPCS Code</b>
<b>EJ76-BELT</b>	Positioning belt (plastic buckle) 42" for Transport chairs	E0978