



BLADES

OPHTHALMOLOGY STERILE MICRO SCALPELS

Feather sterile ophthalmology micro scalpels offer a lightweight aluminum handle with non-slip finger grips combined with an angled, ultra-sharp stainless steel blade.

Feather, a blade manufacturer of global renown, utilizes high-precision technology and carefully selected materials. Each Feather scalpel is made in Japan and meets strict ISO standards, resulting in exceptionally high quality and consistent strength and sharpness.

Features

- Sterilized with gamma radiation and hermetically sealed in blister package
- Non-slip finger grip ensures scalpel remains in user's grasp
- Knife can be autoclaved within included heat-resistant storage case
- Color-coded handles for easy identification of blade type
- Ideal for hospital operating rooms and surgical facilities, especially for ophthalmic surgeries involving semilaminar incisions and centesis incisions of the sclera and cornea
- Available in 3 angle options (15°, 30° and 45° angles)
- Individually packaged 5 per dispensing box



Item No.	Angled, Ultra Sharp Blade Paired With Stainless Steel Handle	Blade Angle	Color Code	UOM
2979#15		15°	Blue	5/box dispenser
2979#30		30°	Black	5/box dispenser
2979#45		45°	Gold	5/box dispenser



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2021, GF Health Products, Inc. All Rights Reserved.
GF Health Products, Inc. is the US distributor for FEATHER.
FEATHER is a registered trademark of FEATHER Safety Razor Co., Ltd

NOTE: Dispose of scalpel in an approved sharps container in accordance with your facility's sharps disposal procedures

graham field

GF2100005Rev/C22