



ROLLATORS

WALKABOUT CONTOUR DELUXE FOUR-WHEEL ROLLATOR

The Lumex Walkabout Four-Wheel Contour Deluxe Rollator offers a comfortable, high end and affordable solution to mobility. Designed with an aluminum frame, this lightweight rollator weighs 16.5 lbs and supports up to a 300 lb weight capacity. The Walkabout folds quickly and easily into a compact unit for storage and transport.

Features

- Ergonomic hand grips
- Easy to operate locking and "quick stop" loop brakes (pull up to stop; push down to lock)
- Padded seat for added comfort
- Rear reflectors to increase visibility at night
- Removable padded contoured backbar for extra support
- Adjustable handle height to accommodate broad range of user heights
- Durable 8" wheel for indoor and outdoor use
- Basket that can be placed underneath seat or in front for carrying of personal items
- Includes easy to read, retail-friendly user manual
- Designed to meet the requirements of HCPCS codes: E0143 & E1056
- Maximum Weight Capacity 300 lb - EVENLY DISTRIBUTED
- Five (5) year limited warranty



		Width inside handgrips		overall Width		overall Depth		Seat Width		Seat Depth		Seat Height	
		inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm
RJ4805B	Royal Blue	16	40.64	25	63.50	27.5	69.85	14	35.56	13.5	34.29	22.5	57.15
RJ4805CH	Champagne	16	40.64	25	63.50	27.5	69.85	14	35.56	13.5	34.29	22.5	57.15
RJ4805K	Black	16	40.64	25	63.50	27.5	69.85	14	35.56	13.5	34.29	22.5	57.15
RJ4805L	Lavender	16	40.64	25	63.50	27.5	69.85	14	35.56	13.5	34.29	22.5	57.15
RJ4805R	Burgundy	16	40.64	25	63.50	27.5	69.85	14	35.56	13.5	34.29	22.5	57.15
RJ4805P	Pink	16	40.64	25	63.50	27.5	69.85	14	35.56	13.5	34.29	22.5	57.15

		Wheel diameter		minimum Handle Height		maximum Handle Height		Product Weight		Weight Capacity EVENLY DISTRIBUTED		UOM
		inch	cm	inch	cm	inch	cm	lb	Kg	lb	Kg	
RJ4805B	Royal Blue	8	20.32	34	86.36	39	99.06	16.5	7.48	300	136.08	1 ea
RJ4805CH	Champagne	8	20.32	34	86.36	39	99.06	16.5	7.48	300	136.08	1 ea
RJ4805K	Black	8	20.32	34	86.36	39	99.06	16.5	7.48	300	136.08	1 ea
RJ4805L	Lavender	8	20.32	34	86.36	39	99.06	16.5	7.48	300	136.08	1 ea
RJ4805R	Burgundy	8	20.32	34	86.36	39	99.06	16.5	7.48	300	136.08	1 ea
RJ4805P	Pink	8	20.32	34	86.36	39	99.06	16.5	7.48	300	136.08	1 ea



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2009, GF Health Products, Inc. All Rights Reserved.
Lumex is a registered trademark of GF Health Products, Inc.

Rollators are not intended to be propelled while seated.
Do NOT use this rollator as a wheelchair or transport device.

graham
field

CF0900094Rev.D02z