



## HOMECARE BEDS

# ★★★ PATRIOT

Homecare Bed Series

## FULL-ELECTRIC/LOW HOMECARE BED

The single-motor design homecare bed is now available in a full-electric/low version. The Patriot Full-Electric Homecare Bed offers the same features as the Patriot Semi-Electric Homecare Bed with the additional feature of electric height adjustment. This makes it easier for caregivers to raise and lower the bed for patient care and linen changes. The electric height adjustment also enables the patient to easily position the bed for maximum comfort. The full-electric model can be adjusted to be a conventional height homecare bed OR a low bed for better inventory management!



Slots in frame ensure proper rail placement



Locking pins secure frame to sleep deck

Designed and engineered by the same team who develop the GF Long Term Care beds.

Sleep surface and motors are **MADE IN THE USA!**

- Total assembled bed weighs only 178 pounds (bed and motor)
- Powerful DC motor for effortless height adjustment
- Hi-Lo motor comes pre-installed on foot section for easy set-up and disassembly
- All functions fully-integrated into convenient handheld pendant
- Welded motor-guard protects motor and universal joint from damage during shipment and deliveries
- Available with Fiberboard Decorative Walnut Bed-ends or High-impact Plastic Bed-ends
- Designed to meet the requirements of applicable HCPCS codes
- Maximum patient weight 350 lbs
- 450 lb evenly distributed total weight capacity including patient, mattress, rails, and accessories
- Height Adjustment:
  - Casters in Low Position:*  
Min. 9.5", Max. 20" (at deck level)
  - Casters in Standard Position:*  
Min. 13", Max. 23.5" (at deck level)
- Overall Bed Size: 87" L x 36" W
- Lift Range: Head - 60°, Knee - 30°
- No external pinch points!
- Slots in frame to ensure proper rail placement
- Locking pins secure bed ends to sleep deck. Attached lanyards prevent loss while moving or storing
- Only homecare bed to offer a 4-inch extension kit that can extend either the head, foot, or both head and foot decking. Bed length will be 84" with one kit, 88" with 2 kits
- The only homecare bed to offer optional corner mattress retainers

### Grid Deck offers the following advantages:

- Much stronger frame
- Better mattress support
- No spring sag or bottoming out with heavier patients
- Easier to clean
- No connecting links or springs to lose.



GF Health Products, Inc. is an  
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.  
The most current and complete product information can be found on our website.

[www.grahamfield.com](http://www.grahamfield.com)





## HOMECARE BEDS



### FULL-ELECTRIC/LOW HOMECARE BED

Fiberboard Decorative Walnut Bed-Ends	High Impact Plastic Bed-Ends	Description	UOM	Meets HCPCS
US0458	US0458PL	Full-Electric Bed Only	1/ea	E0297
US0458-RFPKG	US0458PL-RFPKG	Full-Electric Bed with 1633-Rolled Foam Mattress and Full Chrome Rails	1/ea	E0265
US0458-RFPKGHRA	US0458PL-RFPKGHRA	Full-Electric Bed with 1633-Rolled Foam Mattress and Clamp-On Half Chrome Rails	1/ea	E0265
US0458-RPKG	US0458PL-RPKG	Full-Electric Bed with No Mattress and Full Chrome Rails	1/ea	E0265
US0458-RPKGHRA	US0458PL-RPKGHRA	Full-Electric Bed with No Mattress and Clamp-On Half Chrome Rails	1/ea	E0265
US0458	US0458PL	Full-Electric Bed Only	1/ea	E0297
US0468-RFPKG	US0468PL-RFPKG	Full-Electric Bed with 1633- Rolled Foam Mattress and FDA Full Rails	1/ea	E0265
US0468-RFPKGHR	US0468PL-RFPKGHR	Full-Electric Bed with 1633- Rolled Foam Mattress and FDA Half Rails	1/ea	E0265
US0468-RFPKGQR	US0468PL-RFPKGQR	Full-Electric Bed with 1633- Rolled Foam Mattress and FDA Quarter Rails	1/ea	E0265
US0468-RPKG	US0468PL-RPKG	Full-Electric Bed with No Mattress and FDA Full Rails	1/ea	E0265
US0468-RPKGHR	US0468PL-RPKGHR	Full-Electric Bed with No Mattress and FDA Half Rails	1/ea	E0265
US0468-RPKGQR	US0468PL-RPKGQR	Full-Electric Bed with No Mattress and FDA Quarter Rails	1/ea	E0265

#### Warranty:

- Two year warranty for defects in workmanship and materials of mechanical components, frame, and electronics.
- During the warranty period, defective items will be repaired or replaced at manufacturer's option. The warranty does not include any labor charges incurred in replacement part(s) installation or any associated freight or shipping charges to the manufacturer.

#### Bed shipped in three cartons:

1. Foot section and motor.
2. Head section and Hi/Lo rod.
3. Bed end set.



GF Health Products, Inc. is an  
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.  
The most current and complete product information can be found on our website.

[www.grahamfield.com](http://www.grahamfield.com)

