

# Parallel jaw closure. **Even tissue compression.**

## PROXIMATE® Reloadable Linear Stapler (TX)



### Consistent staple formation

- Parallel jaw closure helps to compress tissue evenly within the jaws of the instrument for consistent staple formation

### Easy to use

- Intermediate locking position facilitates proper placement of stapler on tissue
- One handed operation to close, fire and release

### Applicable specialties

- General
- Colorectal
- Thoracic

## Product Overview

- Two staggered rows of titanium alloy staples in 30mm and 60mm sizes.
- Three different reloads accommodates a variety of tissue thicknesses

## Device Specifications

Description	Device Code	Tissue	Staple rows	Staple line length	Number of staples	Max firings	Qty/box
Proximate Linear Stapler 30mm Reload-Vascular	TX30V	<input type="checkbox"/> Vascular	3	30mm	23	8	3
Proximate Linear Stapler 30mm	TX30B TX30G	<input checked="" type="checkbox"/> Regular <input checked="" type="checkbox"/> Thick	2	30mm	11	8	3
Proximate Linear Stapler 60mm	TX60B TX60G	<input checked="" type="checkbox"/> Regular <input checked="" type="checkbox"/> Thick	2	60mm	21	8	3



## Reload Specifications

Description	Reload code	Tissue	Closed staple height (mm)	Wire diameter	Open staple height	Qty/box
Proximate Linear Stapler 30mm Reload-Vascular	XR30V	<input type="checkbox"/> Vascular	1.0mm	0.20mm	2.4mm	12
Proximate Linear Stapler 30mm Reload	XR30B XR30G	<input checked="" type="checkbox"/> Regular <input checked="" type="checkbox"/> Thick	1.5mm 2.0mm	0.23mm 0.27mm	3.6mm 4.7mm	12
Proximate Linear Stapler 60mm Reload	XR60B XR60G	<input checked="" type="checkbox"/> Regular <input checked="" type="checkbox"/> Thick	1.5mm 2.0mm	0.23mm 0.27mm	3.6mm 4.7mm	12



## How to Order

All purchase orders are made to Johnson & Johnson Health Care Systems, Inc. (JJHCS) or your authorized distributor.

### Electronic Ordering Options

(Note: Placing orders electronically avoids minimum order fees for your hospital.)

#### Johnson & Johnson Gateway: [www.jnjgateway.com/commerce](http://www.jnjgateway.com/commerce)

For questions about your order, visit the website or call 1-866-JNJ-GATE (565-4283).

#### Global Healthcare Exchange (GHX): [www.ghx.com](http://www.ghx.com)

For questions about your order, visit the website or call 1-800-YOUR-GHX (968-7449).

#### Electronic Data Interchange (EDI): JJHCS EDI Help Line: 1-800-262-2888

### Nonelectronic/Manual Ordering Options

Call JJHCS at 1-800-255-2500 between 8:30 a.m. and 9:00 p.m. Eastern Standard Time or fax your order to 1-732-562-2212.

## Customer Support

For product use assistance, clinical guidelines, service and repair, emergency assistance, copy of a 510(k) clearance letter or complaints, please contact our Customer Support Center at [customersupport@eesus.jnj.com](mailto:customersupport@eesus.jnj.com) or by calling **1-877-ETHICON** (384-4266). Our Support Center is staffed 24 hours a day, 7 days a week by qualified nurses to answer your product-related questions.