



Electronic UBE Upper Body Ergometer

Electronic Resistance WATTS & METS Commercial Grade

PhysioTrainer PRO	
Technical Details	
Resistance System	Electro Magnetic
Resistance Levels	1 to 16
WATT Range	8 to 220
Programs	21
Handles	360° Swivel Handles
Pedals	OrthoPedic Heel Cups
Cranks	Adjustable Cranks
Motion	Bi-Directional Resistance
Max User Weight	Unlimited
Power	AC 110V-220V 50/60Hz
Product Type	Commercial Grade
Model #	PT-PRO
Computer Details	
Display	High Definition 9" LCD Display
Backlight	Easy to Read High Contrast Blue
USB Summary	Export Workout Data Including Peaks and Averages
Feedback	Time, Distance, Pulse, Steps Per Minute (SPM), Calories, WATTS, METs, & Resistance Level "1-16"
Programs	21 Programs: Manual, Warm Up, Interval, Ramps, Mountain, Rolling, Climb, Random, Hill, Fat Burn, Cardio, HRC 65%, HRC 85%, Constant SPM, Constant Watts, Goal Calories, Goal Distance, User 1, User 2
Heart Rate	Built in Receiver Compatible with Polar and other Manufactures *Heart Rate Receiver Not Included
Warranty	
Frame	Lifetime
Belt Drive	3 Years
Mechanical Parts	3 Years
Electronic Parts	3 Years
Unit Dimensions	
Assembled Dimensions	26" L X 19" W X 31" H
Unit Weight	55 lbs.




9" Color-Coded Display



Transport Handle



Ortho Foot Pedals Included



Optional HCI Tables



HEALTHCARE
INTERNATIONAL

www.hcfitness.com

