

# Rescue Sheet

THE RESCUE SHEET IS WIDELY ACCEPTED AS A VITAL AID IN THE RAPID EVACUATION OF NON-AMBULATORY PATIENTS DURING EMERGENCY SITUATIONS.

The Rescue Sheet is located permanently and unobtrusively under the mattress and does not interfere with bed making. The patient is rescued with their bed mattress to offer protection under the patient, especially when going down stairs.



## Features

### 1 Long pulling straps

Located at both head and feet ends to allow for multiple rescuers

### 2 Strong easy slide material

For easy and efficient patient movement

### 3 Rubber side panels

Prevents the patient slipping when getting in to or out of bed, which lift clear of the floor when secured for an evacuation

### 4 Reinforcing cross straps

Sewn all the way across the sheet adding reinforcement when being pulled across a rough surface

### 5 Installation label

'This Side Up' label insuring correct installation on the bed, with space for ward name or number

### 6 Plastic strap buckles

Safer than metal buckles, protecting against injury and fully washable

### 7 Cross strap press studs

Press studs contain the cross straps which help prevent straps falling and catching in bed frame or rails

### 8 Pulling strap press studs

Press studs for the pulling straps help keeping them tucked under the mattress, and easily accessible



## Includes

- 660lbs Weight capacity
- Fits all bed models and mattresses
- Heavy duty construction for maximum efficiency
- Designed to enable one or more persons to evacuate a person in their mattress during an emergency

## Dimensions

Height	75"
Width	36"
Weight Capacity	660lbs

**VIDEO LINK:** [evac-chair.com/evacuation-sleds/rescue-sheet/](http://evac-chair.com/evacuation-sleds/rescue-sheet/)

For over 35 years EVAC+CHAIR® is America's Original and World's Number One Selling Evacuation Chair. Don't Risk Lives, Beware of Poor Imitations and Inferior Substitutions. There is only one EVAC+CHAIR®.

Unsurpassed Quality and Workmanship, Successful in Saving Lives in 9/11, 2 YEAR MANUFACTURERS WARRANTY.

For further information please contact:

**Evac+Chair North America LLC**

3000 Marcus Ave, Suite #3E6, Lake Success, NY, 11042-1012

T 516-502-4240 W [www.evac-chair.com](http://www.evac-chair.com) E [sales@evac-chair.com](mailto:sales@evac-chair.com)

**EVAC+CHAIR®**

THE WORLD'S NO.1 STAIRWAY EVACUATION CHAIR SINCE 1982