

# Mechanical Chair Scale

## Standard Features

- Capacity: 440 lb x 0.25 lb (200 kg x 100 g)
- Swing-up arms for easy access and exit
- Fold up footrest
- Sturdy enameled, steel frame
- Dual read die cast beam
- ABS plastic seat 17 in (D) x 18.3 in W
- Swing-up armrests
- Heavy duty wheels - two back wheels lock, two front wheels swivel

## Specifications

### Capacity:

450 lb x 0.25 lb (200 kg x 0.1 kg)

### Seat Size:

17 in (D) x 18.3 in W

### Warranty:

Two-year limited warranty

## Part Number/Price

PART #	MODEL	CAPACITY	EST. SHIP WEIGHT	PRICE
111326	RL-MCS	440 lb (200 kg) x 4 oz (100 g)	85 lb	Consult



Model RL-MCS

