


Modalities

Catalog



Richmar[®]



OUR PROMISE



INNOVATIVE



HIGH QUALITY



CUSTOMER-CENTRIC



PRODUCT SUPPORT

To the physical therapists, chiropractors, and athletic trainers, and to our distributors: Richmar is the only dealer's manufacturer that offers high-quality, innovative medical products supported by superior customer service, marketing support and clinical product training. This is our promise to our partners.

TABLE OF CONTENTS

1 Clinical Combination Therapy

4 Clinical Electrotherapy

7 Clinical Ultrasound

10 Deep Oscillation

12 Diathermy

13 Laser & Light Therapy

16 Hot/Cold Therapy

22 Portable Electrotherapy

32 Portable Ultrasound

35 Electrodes

41 Accessories

47 Topical Analgesics

48 Kinesio Tape

52 Bracing & Support



WINNER EVO COMBO

The Winner EVO combination unit provides the clinician with the most comprehensive therapies in one unit with a reliable platform that is simple to use. Winner EVO CM4 can run 4 channels of electrical stimulation, ultrasound and laser and / or light all at the same time ensuring maximum utilization and return on investment in a busy clinic.

PRODUCT FEATURES

- Dual Frequency (1 & 3 MHz) Clinical Ultrasound with Industry Leading BNR
- Continuous or Pulsed Modes
- 5 cm², 2 / 5 cm² Therapy Hammer, 10 cm², or AutoSound Transducers Available
- Combination, Ultrasound, Electrotherapy, and/or Laser in One
- Tru-Stim Electrotherapy with Five Waveforms On-Board
- Therapy Profiles for Commonly Used Treatments
- Lead Cord Continuity Test on Board All Winner EVO Series
- Accessory Holders for Easy Applicator Organization
- Available AutoSound and Light Cluster Probe Allow for Hands-Free Therapy
- Optional Storage Cart Available
- Three-Year Warranty
- Proudly Made in USA

ORDERING DETAILS

410-012	Winner EVO CM4 4 Channel Combo 120V
410-012-220	Winner EVO CM4 4 Channel Combo 220V
410-010	Winner EVO CM2 2 Channel Combo 120V
410-010-220	Winner EVO CM2 2 Channel Combo 220V

STANDARD ACCESSORIES

410-160	5 cm ² Applicator
A1717A	Lead Cord Set 1.8m (6 ft.) CM2 (2 x); CM4 (4 x)
203-016	US Hospital Grade Power Cord 2.4m (8 ft.)
LS5266	TheraSonic Ultrasound Gel, 16 ounces
410-250	Instruction Manual
	MultiStim Electrodes, 2pk

OPTIONAL ACCESSORIES

410-040	StoreMore Therapy Cart*
DQCART-2	Therapy Cart
410-200	AutoSound Hands-Free Ultrasound Applicator
410-115	Therapy Hammer 2 cm ² / 5 cm ² (2 in 1)
410-176	10 cm ² Applicator (1 MHz)
600-123	Light Cluster Probe, Hands-Free (500 mW - 870 nm LED / 633 nm SLD)
LCT450	LaserPrism Laser Probe, Tethered (450 mW / 808 nm)
LCS450	LaserPrism Laser Probe, Cordless (450 mW / 808 nm)
410-060A	Laser Module - Winner EVO (Required for Tethered Laser / Light Probes)
K36031H5183BB-D	UK Power Cord
K5A031H5183BB	Euro Power Cord

TECHNICAL SPECIFICATIONS

US Frequency	1 & 3 MHz
US Output	Continuous & Pulsed

*Product discontinued 2019; available while supplies last.



THERATOUGH® CX4 COMBO

The CX4 is an advanced four-channel combination therapy system providing electrotherapy and ultrasound treatment options in a simple-to-use, ergonomic design.

The CX4's bilingual interface intuitively groups and displays clinical protocols to guide clinicians step-by-step to ensure an ideal therapeutic treatment, including an anatomical library with electrode placement recommendations to ensure treatment efficacy and clinic efficiency. Therapeutic ultrasound can be combined with stimulation using the dedicated stim channel with any supplied waveform.

PRODUCT FEATURES

- Dual Frequency (1 & 3 MHz) Clinical Ultrasound
- Pulsed and Continuous Operating Modes
- Ergonomically Designed 5 cm² Applicator Included with a Maximum BNR Ratio of 5:1
- Rapid Agonist Antagonist Sequencing (RAAS) Is a Unique Stimulation Waveform Which Replicates the Typical Firing Pattern of Normal Muscle Contractions, and Used to Effectively Re-Educate Weak, Inhibited Muscles Following Surgery, Injury, or Trauma
- Eight Preset Clinical Protocols for Stimulation, with Ability to Customize Parameters and Save Up to 99 User-Defined Protocols for Future Use
- Bilingual Interface Selectable Between English and Spanish
- Three-Year Warranty
- ETL Certified

ORDERING DETAILS

DQ8200	CX4 Clinical Electrotherapy and Ultrasound System
DQ8001	CX4 Clinical Electrotherapy and Ultrasound System with Therapy Cart

STANDARD ACCESSORIES

DQ8000X	Medical Grade Power Cord
DQ8002W5	5 cm Applicator
LS5288	TheraSonic Conductive Gel, 8 oz.
DQSWITCH	Patient Interrupt Switch
400-8771	MultiStim Electrodes, 2" x 2", 4/pk (2x)
400-8791	MultiStim Electrodes, 3" x 5", 2/pk
400-8521	MultiStim Electrodes, 2" x 3.5", 4/pk (2x)
WQ8112	Lead Wires, 110", 4/pk
WQ8112-SNGL	Combination Lead Wire

OPTIONAL ACCESSORIES

DQ8002W1	1 cm ² Ultrasound Applicator
DQCART-2	Therapy Cart (only standard with DQ8001)

TECHNICAL SPECIFICATIONS

Power	100 / 240 VAC
Frequency	50 / 60 Hz
Waveforms	Biphasic, NMS, High Volt, PreMod, IFC, RAAS, Russian, Microcurrent
US Frequency	1 & 3 MHz
US Output	Continuous & Pulsed



COMBOCARE®

The ComboCare is a professional two-channel device that combines clinical ultrasound therapy with five electrotherapy treatment waveforms: IFC, Premod, Russian, EMS and TENS waveform. Its advanced operating system ensures a complete range of therapy.

The large easy-to-read LCD panel allows the clinician to quickly set treatment parameters. A dual frequency 5 cm² applicator comes standard with the device. Its portable design lets you take it to your patient or keep it at the clinic. It's easy to use, and is up and running within minutes.

PRODUCT FEATURES

- Portable Therapy - Provide Ultrasound and E-Stim Simultaneously, Great for At the Clinic or On the Go
- Two Channels of Electrical Stimulation Output
- Large Easy-to-Read Graphical LCD Display with Exact Dial Performance and Programmable Timer
- Ultrasound Has 10% - 100% Duty Cycle Selections
- 10 Presets Per Waveform
- 20 Ultrasound Presets
- Several Different Amplitude Modulation Options
- Interferential and Premodulated Modes Offer Frequency Modulation As Well As a Static Frequency Option
- User Defined Memory Positions for User Protocols
- Designed for Tabletop, Therapy Cart or Mobile Use
- Includes 5cm² Applicator, Two Sets of Electrodes, Two Sets of Rubber Electrodes with Sponges, One Cable Connector, Two Elastic Straps, Two Sets of Lead Wires, One Combination Therapy Lead Wire, Power Cord and Instruction Manual
- Two-Year Warranty

TECHNICAL SPECIFICATIONS

- Power: 100 / 240 VAC
- Frequency: 50 / 60 Hz
- Waveforms: TENS, EMS, Premodulated, IFC, Russian
- Ultrasound Frequency: 1 & 3 MHz
- Ultrasound Output: Continuous & Pulsed

ORDERING DETAILS

DQ7844 ComboCare

STANDARD ACCESSORIES

EIP2020WC2	Polybag White Cloth Electrodes, 2" x 2"
EIP2035WC2	Polybag White Cloth Electrodes, 2" x 3.5"
ER2535B2	Rubber Electrodes, 2.5" x 3.5", 4/pk
ER2743B2	Rubber Electrodes, 2.75" x 4.3", 4/pk
DQ9275W5	5 cm Applicator
ES2740Y2	Electrode Sponges, 2.75"x4", 4/pk
ES3047Y2	Electrode Sponges, 3"x4.75", 4/pk
DQ9276MGC	SoundCare/Combo Care Medical Grade Cord
EW3023BW2	Elastic Strap, 3"x23.5", 2/pk
EW3047B2	Elastic Strap, 3"x47", 2/pk
WW1001PLUG	Combination Lead Wire
WQ6005PLUG-2	Lead Wires, 78", 4/pk
WQ8450CPLUG	Cable Connector for Lead Wires, 2/pk
LS5288	Therasonic 8oz ultrasound gel

OPTIONAL ACCESSORIES

DQ9275W1 1 cm² Applicator



WINNER EVO STIM

For those clinicians only interested in utilizing electrical stimulation within their practice, the Winner EVO Stim device maintains the rugged durability clinicians have come to expect. The Winner EVO Stim has proven that it will provide years of worry-free treatments with up to four channels of electrotherapy.

PRODUCT FEATURES

- EVO Tru-Stim Electrotherapy with Five Waveforms on Board
- User-Defined Profiles for Commonly Used Treatments
- Easily Incorporate Laser/Light Therapy Now or Upgrade Later
- Lead Cord Continuity Test on Board All Winner EVO Series
- Accessory Holders for Easy Applicator Organization Included
- Optional Storage Cart Available
- Three-Year Warranty
- Proudly Made in USA

ORDERING DETAILS

410-016	Winner EVO ST4 4 Channel Stim 120V
410-016-220	Winner EVO ST4 4 Channel Stim 220V
410-014	Winner EVO ST2 2 Channel Stim 120V
410-014-220	Winner EVO ST2 2 Channel Stim 220V

STANDARD ACCESSORIES

A1717A	Lead Cord Set 6ft. (1.8 m) ST2 (2x); ST4 (4x)
203-016	US Hospital Grade Power Cord 8ft. (2.4 m)
410-250	Instruction Manual
	MultiStim Electrodes, 2pk

OPTIONAL ACCESSORIES

410-040	StoreMore Therapy Cart*
600-123	Light Cluster Probe, Hands-Free (500mW - 870nm LED/633nm SLD)
LCT450	LaserPrism Laser Probe, Tethered (450 mW / 808 nm)
LCS450	LaserPrism Laser Probe, Cordless (450 mW / 808 nm)
410-060A	Winner Laser Module (Required for Laser / Light Emitter)
K36031H5183BB-D	UK Power Cord
K5A031H5183BB	Euro Power Cord

TECHNICAL SPECIFICATIONS

Power	115 / 230 VAC
Frequency	50 / 60 Hz
Waveforms	Interferential, Premodulated, High Volt, Russian, Microcurrent

*Product discontinued 2019; available while supplies last.



THERATOUGH® EX4 STIM

The EX4 is an advanced four-channel electrotherapy system offering a wide range of treatment options in one user-friendly, ergonomic design.

Equipped with a bilingual menu-driven touchscreen interface, the EX4 guides clinicians step-by-step through treatment protocols to ensure the ideal therapeutic treatment – including electrode placement for better treatment outcomes.

PRODUCT FEATURES

- Rapid Agonist Antagonist Sequencing (RAAS) Is a Unique Stimulation Waveform Which Replicates the Typical Firing Pattern of Normal Muscle Contractions, and Used to Effectively Re-Educate Weak, Inhibited Muscles Following Surgery, Injury, or Trauma
- Eight Preset Clinical Protocols for Stimulation, with Ability to Customize Parameters and Save Up to 99 User-Defined Protocols for Future Use
- Bilingual Interface Selectable Between English and Spanish
- Three-Year Warranty
- ETL Certified

ORDERING DETAILS

DQ7200	EX4 Clinical Electrotherapy System
DQ7001	EX4 Clinical Electrotherapy System with Therapy Cart

STANDARD ACCESSORIES

DQ8000X	Medical Grade Power Cord
DQSWITCH	Patient Interrupt Switch
400-8771	MultiStim Electrodes, 2" x 2", 4/pk, (2x)
400-8791	MultiStim Electrodes, 3" x 5", 2/pk
400-8521	MultiStim Electrodes, 2" x 3.5", 4/pk
WQ8112-SNGL	Clinical Grade Lead Wires, Pin / Pin Connectors, 110"
WQ8112	Clinical Grade Lead Wires, Pin / Pin Connectors, 110", 4/pk

OPTIONAL ACCESSORIES

DQCART-2	Therapy Cart (only standard with DQ7001)
----------	--

TECHNICAL SPECIFICATIONS

Power	110 / 240 VAC
Frequency	50 / 60 Hz
Waveforms	RAAS, Interferential, Premodulated, Biphasic, NMS, High Volt, Russian, Microcurrent

QUATTRO 2.5™



The microprocessor controlled Quattro 2.5 is a professional four-channel electrotherapy device that provides IFC and Premodulated, Russian, EMS and TENS waveforms with two separate timers. Choose between several different amplitude modulation options. The interferential and pre-modulated modes offer frequency modulation as well as a static frequency option. The portable design lets you take it to your patient or keep it at the clinic.

PRODUCT FEATURES

- 4-Channel Electrotherapy with Separate Timers for Simultaneous Treatment for Two Treatment Areas or Two Patients At Once
- Large Easy-to-Read Graphical LCD Display with Exact Dial Performance and Two Programmable Timers
- Preset Protocols As Well As Customized User Defined Protocols
- Waveforms Include TENS, EMS, Premod (IF-2P), IFC (IF-4P) and Russian
- Constant Current and Constant Voltage Modes
- 10-Presets Per Waveform Totaling 50-Presets in All (Plus Professional Mode)
- Designed for Tabletop, Therapy Cart or Mobile Use
- Includes Four Sets of Electrodes in Two Different Sizes, Four Sets of Carbon Rubber Electrodes with Sponges in Two Different Sizes, Four Straps, Four Sets of Lead Wires, Four Cable Connectors and a User Manual with Protocol Guide
- Two-Year Warranty

ORDERING DETAILS

DQ8450 Quattro 2.5 Professional Electrotherapy Device

STANDARD ACCESSORIES

EIP2020WC2	Polybag White Cloth Electrodes, 2" x 2"
EIP2035WC2	Polybag White Cloth Electrodes, 2" x 3.5"
ER2535B2	Rubber Electrodes, 2.5" x 3.5", 4/pk
ER2743B2	Rubber Electrodes, 2.75" x 4.3", 4/pk
ES2740Y2	Electrode Sponges, 2.75" x 4", 4/pk
ES3047Y2	Electrode Sponges, 3" x 4.75", 4/pk
DQ2001MGC	Medical Grade Power Cord
EW3023BW2	Elastic Strap, 3" x 23.5", 2/pk
EW3047B2	Elastic Strap, 3" x 47", 2/pk
WQ6005PLUG-2	Lead Wires, 78", 4/pk
WQ8450CPLUG	Cable Connector for Lead Wire, 2/pk

TECHNICAL SPECIFICATIONS

Power	110 / 240 VAC
Frequency	50 / 60 Hz
Waveforms	TENS, EMS, Premodulated, IFC, Russian



THERASOUND EVO

Richmar's heritage is ultrasound. The legacy products from Richmar were among the first ultrasound devices engineered more than forty years ago. That same knowledge, craftsmanship and quality reside in today's TheraSound. A dual frequency clinical ultrasound device that has an industry leading BNR, TheraSound can also become a hybrid ultrasound/laser modality with two independent channels. This field upgrade can be done now or later.

PRODUCT FEATURES

- Dual Frequency Clinical Ultrasound (1 & 3 MHz) with Industry Leading BNR
- Continuous or Pulsed Duty Cycles
- Patented Therapy Hammer 2 & 5 cm² Combined in One Transducer
- TheraSound EVO Upgradeable to Accept Laser
- Both AutoSound and Light Cluster Probe Allow for Hands-Free Therapy
- With the Inclusion of a TheraSound Cradle, Integration with a StoreMore Cart Is Available
- Three-Year Warranty
- Proudly Made in USA

ORDERING DETAILS

410-020	TheraSound EVO (120V)
410-028	TheraSound EVO (220V)
410-026	TheraSound 3 Standard (not shown)

STANDARD ACCESSORIES

TheraSound EVO	
410-115	Therapy Hammer 2 cm ² / 5 cm ² Ultrasound Applicator, (2 in 1)
410-125	Instruction Manual
LS5266	TheraSonic Ultrasound Gel, 16 oz
TheraSound 3 Standard	
410-160	5 cm ² Applicator
410-125	Instruction Manual
LS5266	TheraSonic Ultrasound Gel, 16 oz

OPTIONAL ACCESSORIES

410-040	StoreMore Therapy Cart*
410-042	TheraSound Cradle (Mount unit to cart)
410-200	AutoSound Hands-Free Ultrasound Applicator
410-160	5 cm ² Ultrasound Applicator
410-176	10 cm ² Ultrasound Applicator (1 MHz)
LCS450	LaserPrism Laser Probe, Cordless (450 mW / 808 nm)

TECHNICAL SPECIFICATIONS

Power	115 / 230 VAC
Frequency	50 / 60 Hz
US Frequency	1 & 3 MHz
US Output	Continuous & Pulsed

*Product discontinued 2019; available while supplies last.



SOUNDCARE PLUS®

The SoundCare Plus is a professional ultrasound device, designed by clinicians. With this unique design, the SoundCare Plus comes equipped with two applicators, 1cm² and 5cm², that are plugged into the device allowing the user to switch between applicators with a single button push.

The SoundCare Plus has 20 presets and operates in 1 MHz and 3 MHz frequencies. The portable design lets you take it to your patient or keep it at the clinic. It's easy to use, and up and running within minutes.

PRODUCT FEATURES

- For the Treatment of Chronic and Acute Muscular Pain
- Includes Two Applicators, 1cm² and 5cm²
- Fully Functional 1 MHz and 3 MHz Frequencies
- Simultaneously Switch Between Applicators with the Push of a Single Button
- 20 User-Defined Presets with Pulsed and Continuous Therapy
- Variable Duty Cycle Outputs
- Ergonomic Handle Design
- Two-Year Warranty

ORDERING DETAILS

DQ9275 SoundCare Plus Professional Ultrasound Device

STANDARD ACCESSORIES

DQ9275W5 5 cm Applicator
 DQ9276MGC Medical Grade Power Cord
 LS5288 TheraSonic Ultrasound Gel, 8 oz
 DQ9275W1 1 cm² Applicator

TECHNICAL SPECIFICATIONS

Power 100 / 240 VAC
 Frequency 50 / 60 Hz
 US Frequency 1 & 3 MHz
 US Output Continuous & Pulsed



AUTOSOUND

Maximizing clinical efficiency while ensuring positive clinical outcomes can be achieved through the utilization of the time-tested AutoSound.

PRODUCT FEATURES

- Patented, Hands-Free Ultrasound Applicator
- 1 MHz, 3 MHz As Well As 1 & 3 MHz Sweep
- Applicator Includes Swivel Attachment with Three FabriFoam Velcro Straps for Easy Fixation to the Body
- Sequences Between Each of the Four 3.5 cm² Transducers for 1 Second Each
- Easily Accessorize the Winner EVO CM2 / CM4 or TheraSound EVO Device for Hands-Free Therapy
- One-Year Warranty
- Proudly Made in USA

ORDERING DETAILS

410-200 AutoSound Hands-Free Ultrasound Applicator

STANDARD ACCESSORIES

A1220 Swivel Attachment
with Three FabriFoam Velcro Straps

400-046 AutoGel Pads (50x / case)

TECHNICAL SPECIFICATIONS

US Frequency 1 & 3 MHz



HIVAMAT® 200 EVIDENT

The Hivamat 200 Evident provides deep oscillation therapy, an FDA-approved treatment for pain relief, muscle relaxation, and increased circulation to help the body's natural healing processes. This clinical device gives therapists the ability to edit and store parameters within the device, or to a blank patient card. A simple, full-color interface with animated guides make for an easy-to-use therapy system.

PRODUCT FEATURES

- Full-Color Touchscreen Display with Simple-to-Use Interface
- Pre-Populated Treatment Programs
- Area-Specific Instructional Application Animations
- User-Controlled Parameters
- Two Independent Treatment Channels
- Save Treatment Programs to Device or Patient Card
- Three-Year Warranty
- Optional Storage Trolley

ORDERING DETAILS

RM02037 Hivamat 200 Evident

STANDARD ACCESSORIES

RM00379	Handheld Applicator (4x)
RM00261	Connection Cable (4x)
RM00778	US Hospital Grade Power Cord
RM00381	Oscillator Head, 5 cm (2x)
RM00386	Oscillator Head, 9.5 cm (2x)
RM00348	Powder
RM00382	Titanium Neutral Element (2x)
RM00229	Patient Therapy Card (3x)
RM00260	Patient Lead Wire (2x)
RM00262	Connection Cable Snap (4x)
201-131HV	Hivamat Electrodes 2/pk (20x)
RM00200	Treatment Gloves (Box of 100)

OPTIONAL ACCESSORIES

RM01004	Evident Trolley
201-131	Hivamat Electrodes, 2/pk (50x)

TECHNICAL SPECIFICATIONS

Power	100 / 240 VAC
Frequency	50 / 60 Hz



HIVAMAT® 200 PORTABLE

Whether for traveling athletes, easy transportation between different floors of a facility, or self-treating at home between office visits, the Hivamat 200 Portable provides deep oscillation therapy wherever it's needed. Using practice-specific Protocol Cards, Hivamat 200 Portable provides targeted treatment options to quickly begin therapy.

PRODUCT FEATURES

- Simple-to-Use Interface
- Practice-Specific Treatment Protocol Cards: Rehab, Athletic, Chiropractic, Continuing Care
- Easily Transportable / Battery Operated
- Three-Year Warranty

ORDERING DETAILS

RM01009AT	Hivamat 200 Portable - Athletic
RM01009DC	Hivamat 200 Portable - Chiropractic
RM01009CC	Hivamat 200 Portable - Continuing Care
RM01009PT	Hivamat 200 Portable - Rehab

STANDARD ACCESSORIES

RM00277	Battery Charger (US plug)
RM00348	Powder
RM00262	Connection Cable Snap (2x)
201-131HV	Hivamat Electrodes, 2/pk (20x)
RM55035	Transport Bag
RM00200	Treatment Gloves (Box of 100)

OPTIONAL ACCESSORIES

RM1012	Self Care Kit: - Handheld Applicator - Connection Cable (2x) - Oscillator Head, 5 cm (2x) - Oscillator Head, 9.5 cm (2x) - Titanium Neutral Element
201-131	Hivamat Electrodes, 2/pk (50x)

TECHNICAL SPECIFICATIONS

Power	1.2V AA NiMH Rechargeable Batteries, 4/pk
-------	---

THERATOUGH® DX2

Richmar has expanded the capability of traditional diathermy devices with the addition of several key innovations and features found only in the TheraTouch® DX2. Designed to provide safe and effective treatment options with modern conveniences for the efficient clinic, TheraTouch® DX2 is today's diathermy.



PRODUCT FEATURES

- Scientific Dosage Control (S.D.C.) Adds Two Dosimetry Settings Within Dose III, for Treatment Options That Are More Effective and Comfortable for a Wide Variety of Patients
- High Efficiency Auto-Tuning (H.E.A.T.) Ensures Maximum Energy Transfer By Constantly Adjusting the Device to Minute Changes in the Patient's Body
- Quick Links By Indication Feature 9 Commonly Used Protocols for Quick Treatment Setup
- 6 Thermal Dosimetry Settings
- Anatomical Treatment Guide for Monode Placement
- Robust Mechanical Arm with 5 Articulation Points for Easy and Reliable Monode Positioning
- Save Up To 99 User-Defined Favorites
- Color Touch Screen Interface
- Three-Year Warranty

ORDERING DETAILS

DQSWD2 TheraTouch® DX2 Shortwave Diathermy

STANDARD ACCESSORIES

DQSWD2-D Coil Monode Drum (14 cm)
 DQSWD2-DC Coaxial Cable for Coil Monode Drum
 DQSWD2-MA Mechanical Arm
 DQSWD2X Medical Grade Power Cord

TECHNICAL SPECIFICATIONS

Power	100 / 240 VAC
Frequency	50 / 60 Hz
Continuous Output	100W
Pulsed Output	200W Peak Power
Power Increments	.5 W
Treatment Duration	1 - 30 minutes, 1 minute increments



THERATOUGH® LX2

Richmar's TheraTouch LX2 simplifies laser therapy with a menu-driven touch screen interface that includes research-based clinical protocols to treat common acute and chronic conditions. An integrated dosage calculator uses multiple units of measure to display energy output based on the set treatment time. The LX2 comes standard with a 9-diode laser applicator, with single diode applicators available as additional accessories.

PRODUCT FEATURES

- Delivers Up to 35 J/cm² in a Single Treatment
- Dosage Reading Automatically Calculates and Displays Output in Joules, J/cm², and Joules Per Diode
- Programmable Treatment Time Up to 4 Minutes
- Continuous or Pulsed Output Modes
- 2 Output Channels
- 10 Predefined Clinical Protocols
- Save Up to 100 User-Defined Protocols
- Available with Therapy Cart

ORDERING DETAILS

DQLLLT TheraTouch® LX2

STANDARD ACCESSORIES

DQLLLT-CA 9 Diode Cluster Applicator

OPTIONAL ACCESSORIES

DQLLLT-SA Single Diode Applicator

TECHNICAL SPECIFICATIONS

Power	100-240V 50/60Hz
Programmable treatment time	Up to 4 minutes
Mode of operation	Continuous or pulsed
Pulsed mode	90%
Emission frequency	1Hz-10000Hz
Laser classification	3B

LASERPRISM LASER PROBE



LCS450



LCT450

PRODUCT FEATURES

- Cordless and Tethered Options Available
- Three 150 mW GaAIAS Laser Diodes with Visible Blue LED Guide Light
- 450 mW of Continuous Power At 808 nm (Frequency)
- Delivers 1 Joule Every 2.2 Seconds, Continuous Duty Cycle
- Three Preset Options: 6, 8, and 10 Joules
- Cordless LaserPrism Uses NiMH Rechargeable Battery, and Comes with Safety Key
- Tethered LaserPrism Connects to Any Winner EVO Via Laser Module Upgrade
- Two-Year Warranty

ORDERING DETAILS

LCS450	LaserPrism - Cordless
LCT450	LaserPrism - Tethered

STANDARD ACCESSORIES

600-207	Laser Goggles (2x)
500-040	LaserPrism (Cordless) 9V Universal Charger (90 - 240 VAC) US, Euro, Aus, China plug configs Nylon Carrying Bag



LIGHT CLUSTER PROBE

- Hands-Free SLD Light Therapy
- 500 mW Total Power, 870 nm LED / 633 nm SLD Frequency
- Delivers 1 Joule / cm² Every 45 Seconds
- Tethered Light Cluster Connects to Any Winner EVO Via Laser Module Upgrade
- Treatment Time Can Be Adjusted If Used in Conjunction with the Tethered Laser
- One-Year Warranty

ORDERING DETAILS

600-123 Light Cluster Probe - Hands Free



WINNER EVO LASER MODULE

- Optional Upgrade to Any Winner EVO Allows Dual Operation of Light/Laser While Performing Up to 4 Channels of Electrotherapy All Simultaneously
- Two Ports for LaserPrism and/or Light Cluster Plug-in
- One-Year Warranty

ORDERING DETAILS

410-060A Laser Module - Winner EVO



HYDRATHERM™

HydraTherm heating units set the standard for modern moist heat therapy with an updated modern design, featuring a rust-free composite tank that will not dent, and unique hot pack storage with a 12 pack capacity.

PRODUCT FEATURES

- Only Two Recommended Cleanings Per Year Thanks to Rust-Free Composite Tank, Compared to 24 for Traditional Stainless Tanks
- Digital Thermostat with Easy-to-Read Temperature Display Gives Therapist Ability to Precisely Target Treatment Temperature, Increasing Patient Safety
- Safe Operation with Low Water Cut-Off, and 24 Volt Control Unit, Compared to 110 Volt Competitors
- HydraTherm Deluxe Models Include an Electric Drainage Pump and Fold-Away Side Shelf
- Hot Pack Storage Options Include the Divider System That Provides Tong-Free Access to Hot Packs, or the Rack System That Allows for Additional Pack Storage on Top of Racks.
- Use in Conjunction with HydraHeat Packs for a Clean, Hygienic Moist Heat Therapy System
- Personal Touch - Custom Order with Facility or Team Logo
- HydraTherms Sold with Packs Include: 3 Universal, 3 Cervical, 3 Standard, and 3 Oversize HydraHeat Packs
- Two Year Warranty When Used with HydraHeat Packs; One Year Warranty When Used with Packs Other Than HydraHeat Packs
- Proudly Made in USA

ORDERING DETAILS

Divider System with Packs*

HT-R12-DW	HydraTherm Deluxe** 120V
HT-R12-DW220	HydraTherm Deluxe** 220V
HT-R12-SW	HydraTherm Standard 120V
HT-R12-SW220	HydraTherm Standard 220V

Divider System without Packs

HT-R12-D	HydraTherm Deluxe** 120V
HT-R12-D220	HydraTherm Deluxe** 220V
HT-R12-S	HydraTherm Standard 120V
HT-R12-S220	HydraTherm Standard 220V

Rack System with Packs*

HT-R12-DRW	HydraTherm Deluxe** 120V
HT-R12-DRW220	HydraTherm Deluxe** 220V
HT-R12-SRW	HydraTherm Standard 120V
HT-R12-SRW220	HydraTherm Standard 220V

Rack System without Packs

HT-R12-DR	HydraTherm Deluxe 120V
HT-R12-DR220	HydraTherm Deluxe 220V
HT-R12-SR	HydraTherm Standard 120V
HT-R12-SR220	HydraTherm Standard 220V

STANDARD ACCESSORIES

HT-R12-RCK	Rack System
HT-R12-DIV	Divider System (includes 3 clip-on shelves)
HTR121	Fold-Away Side Shelf
HT-R12-DIVS	Clip-On Shelves

TECHNICAL SPECIFICATIONS

Power	100 VAC, 240 VAC
Frequency	50 Hz, 60 Hz
Temperature Range	120 °F - 160 °F (71 °C - 49 °C)

* HydraTherm models sold with packs include: 3 Standard, 3 Universal, 3 Cervical, & 3 Oversize HydraHeat Packs. **HydraTherm Deluxe models include electric-powered drain pump, and a side-shelf.



HYDRAHEAT PACKS

HydraHeat Packs paired with HydraTherm heating units create a complete moist heat therapy system that provides patient comfort, cleanliness and a time-efficient clinic.



HP-1713-OS



HP-1905-UN



HP-1711-CX



HP-1711-ST

PRODUCT FEATURES

- The Polymer Skin of HydraHeat Packs Are Easily Disinfected with Antibacterial Wipes, Helping Reduce Cross Contamination Between Patients
- HydraHeat Packs Do Not Require High Temperatures to Inhibit Bacterial Growth As with Clay Packs, So Lower, Safer Temperatures Within the Therapeutic Range Can Be Used to Treat Patients, Minimizing Burn Risk
- Lower Treatment Temperatures Allow HydraHeat Packs to Be Covered with Only a Towel or Pillow Case, Not Bulky, Expensive, and Unsanitary Terry Covers
- The Sureseam Process Used to Seal Each Pack Lessens the Chances of Bursting or Leaking its Contents Into the Tank, or Onto the Patient
- Polymer Skin in Combination with the Non-Toxic Gel Interior Make for a Pliable Pack That Conforms to Various Treatment Areas and Allows a Greater Area of Contact for Heat Transfer
- HydraHeat Packs Maintain Heat Within the Therapeutic Range for Up to 40 Minutes, and Reheat to Treatment Temperature Within 15 Minutes
- Not Made with Natural Rubber Latex

ORDERING DETAILS

HP-1713-OS	Oversize 18" x 13" (45 x 33 cm)
HP-1905-UN	Universal 18" x 14" (46 x 10 cm)
HP-1711-CX	Neck Contour/Cervical 17" x 10.5" (44 x 27 cm)
HP-1711-ST	Standard 17" x 10.5" (49 x 27 cm)

THERAMED® MOIST HEAT PACK



Therapeutic moist heat therapy provides deeper penetration to muscle tissue, allowing temporary relief of pain due to sprains, muscle strains, back pain, bursitis, arthritis and more. Moist heat has also been established to increase tissue temperature at a faster rate than dry heat, causing increased circulation and promoting quicker recovery time. With a full range of sizes and covers available, TheraMED Professional hot packs work with most heating units on the market today and will address all your clinic needs.



PRODUCT FEATURES

- Clay-Filled Packs Provide Up to 30 Minutes of Therapeutic Relief
- Red and White Tabs Visually Assist for Easy Rotation
- Seven Sizes Available to Provide Therapy to Nearly Any Body Part
- Reusable for Cost-Effective Treatments
- Great for Clinical or At-Home Use
- Compatible with the Chattanooga™ Hydrocollator® and Other Comparable Stainless Heating Tanks
- Terry Covers Available in a Variety of Sizes (Sold Separately)
- One-Year Warranty

ORDERING DETAILS

HP0512	Half-Size Moist Heat Pack, 5" x 12"
HP1012	Standard Moist Heat Pack, 10" x 10"
HP1018	Spine Small Moist Heat Pack, 10" x 18"
HP0024	Cervical Moist Heat Pack, 24"
HP1020	Knee/Shoulder Moist Heat Pack, 10" x 20"
HP1024	Spine Large Moist Heat Pack, 10" x 24"
HP1524	Oversize Moist Heat Pack, 15" x 24"



THERAMED® MOIST HEAT TERRY COVERS

TheraMED Terry Covers are designed to cover the TheraMED Moist Heat Packs to insulate heat and protect the patient from excessive heat exposure. Great for high-volume clinic use or therapy at home.



HPC1028-TC



HPC2432-TC



HPC0926-TC



HPC1028-FF



HPC2432-FF



HPC0926-FF

PRODUCT FEATURES

- Hook and Loop Closures to Keep Hot Pack Secure and Prevent Unwrapping or Slipping
- Foam-Filled Covers Are Constructed with a Polyurethane Foam Layer Between Two Layers of Terry
- Machine Washable
- Available in a Variety of Sizes
- One Year Warranty

ORDERING DETAILS

HPC0926-FF	Cervical Terry Cover, Foam-Filled, 9" x 26"
HPC0926-TC	Cervical All-Terry Cover, 9" x 26"
HPC1028-FF	Standard Terry Cover, Foam-Filled, 10" x 28"
HPC1028-TC	Standard All-Terry Cover, 10" x 28"
HPC2432-FF	Oversize Terry Cover, Foam-Filled, 24" x 32"
HPC2432-TC	Oversize All-Terry Cover, 24" x 32"

THERAMED® DUAL MOIST-DRY HEATING PAD



The TheraMED Professional Dual Moist-Dry Heating Pad offers two applications in one for optimal therapy. It can be used dry or moist with the dampened sponge insert. The moist heat option provides deep penetrating relief from sore, stiff muscles, while the dry option provides a comforting warm treatment. The unique pad design delivers a consistent therapeutic heat in the recommended range and eliminates any spikes to maintain an even comfortable level at all times.

PRODUCT FEATURES

- Provides Soothing Relief From Chronic Back Pain, Muscle Pain, Sports Injuries, Post-Surgery Pain, Menstrual Cramps and Arthritis
- Dual Application - Can Be Used Dry or Moist with Dampened Sponge Insert
- Flexible Pad with Extra-Soft Machine Washable Cover (100% Polyester) for a Comfortable Feel
- Four Heat Settings - Warm, Low, Medium and High
- Ergonomically Designed Hand Controller with Simple One-Button Operation and Color-Coded LED Indicators
- 40 Minute Auto Shut-Off Feature
- 8-Foot Extra Long Cord for Easy Positioning and Freedom of Movement During Use
- Available in Two Sizes - Medium and King

ORDERING DETAILS

HP1215	TheraMED Pro Dual Moist-Dry Heat Pad Medium, 12" x 15"
HP1224	TheraMED Pro Dual Moist-Dry Heat Pad King, 12" x 24"



REUSABLE HOT/COLD PACKS

Soothe strains and pain with our reusable, microwaveable hot/cold packs. Stays flexible, even when frozen. Cost effective and long lasting. Can be used in conjunction with other therapies. Three sizes to choose from.



BG5010



BG7511



BG1114

PRODUCT FEATURES

- Made of Soft Touch Material
- Can Be Cleaned Hundreds of Times
- For Treating Injuries Requiring Hot or Cold Therapy
- Microwave for Heat Therapy; Freeze for Cold Therapy

ORDERING DETAILS

BG5010	Hot and Cold Gel Pack, 5" x 10"
BG7511	Hot and Cold Gel Pack, 7.5" x 11"
BG1114	Hot and Cold Gel Pack, 11" x 14"



INTENSITY™ SELECT COMBO II

The InTENSity Select Combo II offers patients four modalities in one user-friendly, ergonomic design - Transcutaneous Electrical Nerve Stimulation (TENS), Neuromuscular Electrical Stimulation (NMES), Interferential Current (IFC), and Russian (RUSS).

The wide-ranging capability of the InTENSity Select Combo II provides the most up-to-date treatment technology and waveforms. Equipped with preset body--part-specific settings as well as a customizable option, it intuitively guides the patient step-by-step to ensure the ideal therapeutic treatment throughout the healing cycle.

PRODUCT FEATURES

- User-Friendly, Sleek Design
- Easy-to-Read Large LCD Display and Convenient Easy-to-Use Controls
- 2 Channels of Independent Stimulation Output
- 4 Modalities in One Portable Device: TENS | NMES | IFC | RUSS
- Extensive Range of Programs with Manual Customize Settings and Body-Part-Specific Settings (Neck, Shoulder, Hand, Lower Back, Knee, Foot) — 18 TENS, 18 NMES, 18 IFC, and 18 RUSS (Total of 72 Programs)
- Electrode Load Indicator
- 4 Stimulation Modes for TENS
- Synchronous Mode for NMES with Adjustable Pulse Amplitude, Pulse Frequency, and Contraction/Relaxation Time Parameter
- Up to 10,000 Hz of IFC Therapy with Sweep Frequency to Prevent Accommodation of Nerves and Provide Maximum Treatment Benefits
- 2500 Hz of RUSS for a Deeper and More Comfortable Stimulation
- Preset and Adjustable Treatment Time
- Compliance Meter - Store Up to 90 Sets of Treatment Records and Total Treatment Time of 100 Hours
- Rechargeable Lithium Ion Battery with USB/ Wall Charger and LED Charging Indicator Light
- One-Year Warranty
- Rx Only

ORDERING DETAILS

DI2195 InTENSity Select Combo Portable TENS, NMES, IFC & RUSS Pain Relief System

STANDARD ACCESSORIES

DI2000-CASE	Carrying Case
WW3005	Lead Wires, 45", 2/pk
TA2000	Lithium Ion Battery
EP2020WC2-INTM	InTENSity Electrodes, 2" x 2", Mini 4/pk
DI2000-AC	AC Charger

OPTIONAL ACCESSORIES

TA2000C	Battery Charger for Lithium Ion Battery
---------	---

TECHNICAL SPECIFICATIONS

Power	3.7 V Lithium Rechargeable Battery Pack
Waveforms	TENS, NMES, IFC, Russian



INTENSITY™ IF COMBO II

The InTENSity IF Combo II combines TENS and Interferential waveforms in one sleek rechargeable design. Equipped with preset body-part-specific settings, the InTENSity IF Combo II intuitively guides the patient step-by-step to ensure the ideal therapeutic treatment throughout the healing cycle.

Simply choose a preset program by body part or customize a program to meet your specific therapy needs. With the powerful rechargeable lithium ion battery, this portable device can be used at home, at work or on-the-go.

PRODUCT FEATURES

- User-Friendly, Sleek Design
- 2 Channels of Independent Stimulation Output
- Easy-to-Read Large LCD Display and Convenient Easy-to-Use Controls
- Extensive Range of Programs with Manual Customize Settings and Body-Part-Specific Settings — 18 TENS and 18 IFC (Total of 36 Programs)
- Electrode Load Indicator
- Up to 10,000 Hz of IFC Therapy with Sweep Frequency to Prevent Accommodation of Nerves and Provide Maximum Treatment Benefits
- Preset and Adjustable Treatment Time
- Compliance Meter - Store Up to 90 Sets of Treatment Records and Total Treatment Time of 100 Hours
- Low Battery Indicator
- Rechargeable Lithium Ion Battery with USB/ Wall Charger and LED Charging Indicator Light
- One-Year Warranty
- Rx Only

ORDERING DETAILS

DI2738 InTENSity IF Combo II

STANDARD ACCESSORIES

DI2000-CASE	Carrying Case
WW3005	Lead Wires, 45", 2/pk
TA2000	Lithium Ion Battery
EP2020WC2-INTM	InTENSity Electrodes, 2" x 2", Mini 4/pk
DI2000-AC	AC Charger

OPTIONAL ACCESSORIES

TA2000C	Battery Charger for Lithium Ion Battery
---------	---

TECHNICAL SPECIFICATIONS

Power	3.7 V Lithium Rechargeable Battery Pack
Waveforms	TENS, IFC



INTENSITY™ TWIN STIM IV

The InTENSity Twin Stim IV offers patients two modalities in one user-friendly, ergonomic design - Transcutaneous Electrical Nerve Stimulation (TENS) and Neuromuscular Electrical Stimulation (NMES).

Equipped with preset body-part-specific settings as well as a customizable option to ensure the ideal therapeutic treatment throughout the healing cycle.

PRODUCT FEATURES

- User-Friendly, Sleek Design
- 2 Channels of Independent Stimulation Output
- Easy-to-Read Large LCD Display and Convenient Easy-to-Use Controls
- Extensive Range of Programs with Manual Customize Settings and Body-Part-Specific Settings — 18 TENS and 18 NMES (Total of 36 Programs)
- The NMES Setting Sends Gentle Electrical Pulses in a Contraction and Relaxation Rhythm to Help Relieve Muscle Spasms, Re-Educate Muscles, and Help Prevent Muscle Atrophy.
- 4 Stimulation Modes for TENS
- Synchronous Mode for NMES with Adjustable Pulse Amplitude, Pulse Frequency, and Contraction/Relaxation Time Parameters
- Preset and Adjustable Treatment Time
- Compliance Meter - Store Up to 90 Sets of Treatment Records and Total Treatment Time of 100 Hours
- Electrode Load Indicator
- Low Battery Indicator
- Rechargeable Lithium Ion Battery with USB/ Wall Charger and LED Charging Indicator Light
- One-Year Warranty
- Rx Only

ORDERING DETAILS

DI2717 InTENSity Twin Stim IV

STANDARD ACCESSORIES

DI2000-CASE	Carrying Case
WW3005	Lead Wires, 45", 2/pk
TA2000	Lithium Ion Battery
EP2020WC2-INTM	InTENSity Electrodes, 2" x 2", Mini 4/pk
DI2000-AC	AC Charger

OPTIONAL ACCESSORIES

TA2000C	Battery Charger for Lithium Ion Battery
---------	---

TECHNICAL SPECIFICATIONS

Power	3.7 V Lithium Rechargeable Battery Pack
Waveforms	TENS, NMES



INTENSITY™ 12

The InTENSity 12 is the newest standard in portable electrotherapy offering patients the latest in treatment options and stimulation programs in a user-friendly, ergonomic design.

Equipped with preset body-part-specific settings as well as a customizable option, the InTENSity 12 intuitively guides the patient step-by-step to ensure the ideal therapeutic treatment.

This highly customizable TENS unit has modes to reach most recurring and chronic pains. Patients will love the convenience of the powerful rechargeable lithium ion battery. Take it to work, use it at home or while traveling.

PRODUCT FEATURES

- User-Friendly, Sleek Design
- 2 Channels of Independent Stimulation Output
- Easy-to-Read Large LCD Display and Convenient Easy-to-Use Controls
- 6 Body-Part-Specific Settings (Neck, Shoulder, Hand, Lower Back, Knee, Foot) Each Equipped with Two (2) Preset Programs, and One (1) Manual Customize Setting - 18 Program Settings Total
- Preset and Adjustable Treatment Time (5-90 Minutes)
- Compliance Meter - Store Up to 90 Sets of Treatment Records, and Total Treatment Time of 100 Hours
- Electrode Load Indicator
- Low Battery Indicator
- Rechargeable Lithium Ion Battery with USB/ Wall Charger and LED Charging Indicator Light
- One-Year Warranty
- Rx Only

ORDERING DETAILS

DI1212 InTENSity 12

STANDARD ACCESSORIES

DI2000-CASE	Carrying Case
WW3005	Lead Wires, 45", 2/pk
TA2000	Lithium Ion Battery
EP2020WC2-INTM	InTENSity Electrodes, 2" x 2", Mini 4/pk
DI2000-AC	AC Charger

OPTIONAL ACCESSORIES

TA2000C	Battery Charger for Lithium Ion Battery
---------	---

TECHNICAL SPECIFICATIONS

Power	3.7 V Lithium Rechargeable Battery Pack
Waveforms	TENS



INTENSITY SELECT COMBO™

The InTENSity Select Combo is a powerful, easy-to-use portable electrotherapy device. Complete with four waveforms, the unit provides the flexibility to treat a wide variety of musculoskeletal issues. With a strong current of 105 mA, this innovative device combines TENS, EMS, Interferential and Microcurrent and comes complete with an AC adapter making home use easier than ever

PRODUCT FEATURES

- Intended for Use of Pain Relief, Muscle Reeducation, Increasing Local Blood Circulation, and Maintaining Range of Motion
- Simple-to-Operate, Portable Electrotherapy Device for Pain Management
- Four Waveforms: TENS, EMS, Interferential and Microcurrent
- Comes Complete with Carrying Case, Pair of Lead Wires, 9-Volt Battery, AC Adapter, Manual, and Pack of Four Self-Adhesive Reusable Electrodes
- One-Year Warranty
- Rx Only

ORDERING DETAILS

DI8195 InTENSity Select Combo

STANDARD ACCESSORIES

CC5081	Carrying Case
WW3005	Lead Wires, 45", 2/pk
TA5013-I	Alkaline 9V Battery
EP2020WC2-INTM	InTENSity Electrodes, 2" x 2", Mini 4/pk
DI0001X	AC Adapter

OPTIONAL ACCESSORIES

TA9265-I	9V Charger (Holds Two 9V Batteries)
TA9225-2	TENS Cell Rechargeable 9V Battery

TECHNICAL SPECIFICATIONS

Power	9V Alkaline DC
Waveforms	TENS, EMS, IFC, Microcurrent



INTENSITY™ IF COMBO

The InTENSity IF Combo combines TENS and Interferential treatment, giving you the best combination for treating patients that need something stronger than TENS therapy alone. The IF Combo has four modes of operation in both TENS and Interferential, ensuring that your patient is getting the most pain relief. The unit is 105 mA strong and operates on a standard 9-volt battery or with an AC adapter so the patient can get relief virtually anywhere.

PRODUCT FEATURES:

- Simple-to-Operate, Portable Electrotherapy Device for Pain Management
- Combines TENS and Interferential Treatment for Patients Needing Stronger Therapy
- Equipped with Four Treatment Modes
- Comes Complete with Carrying Case, Pair of Lead Wires, 9-Volt Battery, AC Adapter, Manual, and Pack of Four Self-Adhesive Reusable Electrodes
- Power Source: Standard 9-Volt Battery or AC Adapter
- One-Year Warranty
- Rx Only

ORDERING DETAILS

DI4738 InTENSity IF Combo

STANDARD ACCESSORIES

CC5081	Carrying Case
WW3005	Lead Wires, 45", 2/pk
TA5013-I	Alkaline 9V Battery
EP2020WC2-INTM	InTENSity Electrodes, 2" x 2", Mini 4/pk
DI0001X	AC Adapter

OPTIONAL ACCESSORIES

TA9265-I	9V Charger (Holds Two 9V Batteries)
TA9225-2	TENS Cell Rechargeable 9V Battery

TECHNICAL SPECIFICATIONS

Power	9V Alkaline DC
Waveforms	TENS, IFC



TWIN STIM® III

The Twin Stim III is a portable digital device that combines TENS for pain treatment and management and EMS for muscle rehabilitation and reeducation. When operating in TENS mode, the unit has four modes and when operating in EMS mode, three modes. The unit is among the upper echelon of pain relief devices, with a maximum output of 105 mA, and it also comes standard with an AC Adapter and removable belt clip for easy, at-home use and convenience.

PRODUCT FEATURES

- Intended for Use of Pain Relief, Muscle Reeducation, Increasing Local Blood Circulation and Maintaining Range of Motion
- Unit Combines TENS and EMS Therapy for Optimal Treatment
- Simple-to-Operate, Portable Electrotherapy Device for Pain Management
- Four Modes for TENS Therapy: Burst, Normal, Mod 1, and Mod 2
- Three Modes for EMS Therapy: Synchronous, Alternate, and Delay
- Comes Complete with Carrying Case, Pair of Lead Wires, 9-Volt Battery, AC Adapter, Manual, and Pack of Four Self-Adhesive Reusable Electrodes
- One-Year Warranty
- Rx Only

ORDERING DETAILS

DI3717 InTENSity Twin Stim III

STANDARD ACCESSORIES

CC5081	Carrying Case
WW3005	Lead Wires, 45", 2/pk
TA5013-I	Alkaline 9V Battery
EP2020WC2-INTM	InTENSity Electrodes, 2" x 2", Mini 4/pk
DI0001X	AC Adapter

OPTIONAL ACCESSORIES

TA9265-I	9V Charger (Holds Two 9V Batteries)
TA9225-2	TENS Cell Rechargeable 9V Battery

TECHNICAL SPECIFICATIONS

Power	9V Alkaline DC
Waveforms	TENS, EMS



INTENSITY™ 10

The InTENSity 10 Digital TENS electrotherapy device allows patients to select from 10 preset body icons for targeted therapy, making it the easiest and most effective TENS device available on the market today. With the press of two buttons, the easy-to-use device is ready to operate effectively based on the area of the body that is being treated. Comes complete with AC adapter for effective therapy on the go.

PRODUCT FEATURES:

- Simple to Operate, Portable Electrotherapy Device for Pain Management
- Two-Button Operation
- Device Has Ten Preset Body Part Selections Making It Easy for Patients to Target Certain Areas in Pain
- Body Part Selections Include: Neck, Shoulder, Mid Back, Low Back, Elbow, Wrist, Hip, Knee, Ankle, and Foot
- Individual Body Part Selections Have a Customized Treatment Program Based on Researched Protocols
- Comes Complete with Carrying Case, Pair of Lead Wires, 9-Volt Battery, AC Adapter, Manual, and Pack of Four Self-Adhesive Reusable Electrodes
- One-Year Warranty
- Rx Only

ORDERING DETAILS

D11010 InTENSity 10

STANDARD ACCESSORIES


CC5081	Carrying Case
WW3005	Lead Wires, 45", 2/pk
TA5013-I	Alkaline 9V Battery
EP2020WC2-INTM	InTENSity Electrodes, 2" x 2", Mini 4/pk
D10001X	AC Adapter

OPTIONAL ACCESSORIES

TA9265-I	9V Charger (Holds Two 9V Batteries)
TA9225-2	TENS Cell Rechargeable 9V Battery

TECHNICAL SPECIFICATIONS

Power	9V Alkaline DC
Waveforms	TENS



VENTI TENS™

The Venti TENS Digital Pain Relief System is a dual-channel digital TENS device equipped with three waveforms and twenty modes of operation to choose from – twelve presets and eight manual. Used for symptomatic relief and management of chronic pain and acute post-surgical or post-trauma pain.

This incredibly advanced, pocket-size portable device is known for its body-part specific settings, which provide ease of application and set-up by patients. Both, its pulse width and pulse frequency are also adjustable.

PRODUCT FEATURES:

- Dual Channel Digital TENS Unit
- Large LCD Display with Body-Part Diagram for Targeted Pain Relief
- Three Selectable Waveforms: Symmetrical Biphasic, Asymmetrical Biphasic, and Monophasic
- 20 Treatment Modes - 12 Preset and 8 Manual
- Last Setting Recall and Settings Lock
- Treatment Timer and Compliance Meter
- Protective Cover Over Controls
- Low Battery Indicator
- Electrode Placement Suggestions
- Includes Digital TENS Device, Pair of Lead Wires, One Package of 4 Reusable Self-Adhering Electrodes, 2 AA Alkaline Batteries, Zipper Carrying Case, and Instruction Manual
- One-Year Warranty
- Rx Only

ORDERING DETAILS

DT2020 Venti TENS Digital Pain Relief System

STANDARD ACCESSORIES

DT2020-CASE	Zipper Pouch
WW3005	Lead Wires, 45", 2/pk
TA2020	TENS Cell AA Batteries, 2/pk
EP2020WC2-INTM	InTENSity Electrodes, 2" x 2", Mini 4/pk

TECHNICAL SPECIFICATIONS

Power	2 AA alkaline batteries
Waveforms	TENS



TENS 7000™ (OTC)

The TENS 7000 is one of the most popular, and best valued, digital TENS units on the market. This digital TENS unit is proven to be among the most durable and easy-to-use TENS units. The upgraded digital processor is more precise and accurate than ever before. This unit is great for users seeking a digital device that is also equipped with physical amplitude knobs. FDA approved for over-the-counter sale, no prescription required.

PRODUCT FEATURES:

- For Relief of Acute and Chronic Pain Associated Achy, Stiff or Sore Muscle, and Joint Pain
- Modes: Burst, Normal, Modulation, SD1, and SD2
- Dual Channel for Use with Two or Four Electrodes
- Pulse Width and a Pulse Frequency Are Adjustable
- System Includes Device, Lead Wires, Four Self-Adhesive Reusable Electrodes, 9-Volt Battery, Hard Plastic Carrying Case, Quick Start Guide, and Instruction Manual
- Compliance Meter, 60 Records Up to 999 Hours
- Low Battery Indicator
- One-Year Warranty

ORDERING DETAILS

DT7202 TENS 7000 (OTC)

STANDARD ACCESSORIES

CC3001	Black Carrying Case
WW3005	45" Lead Wires, 45", 2/pk
TA9050-I	TENS Cell 9V Heavy Duty Battery
EP2020WC2-INTM	InTENSity Electrodes, 2" x 2", Mini 4/pk

TECHNICAL SPECIFICATIONS

Power	9V Alkaline DC
Waveforms	TENS



THERASOUND 3.5P

The reliability of the TheraSound ultrasound device integrated into a rugged hard case makes the TheraSound 3.5P ideal for the mobile clinician, athletic trainer, physiotherapist, or home health provider.

PRODUCT FEATURES:

- Portable Ultrasound Unit with Affixed Applicator
- Patented Therapy Hammer Applicator
- 2/5 cm² in One Ultrasound Applicator
(with Dual Frequency 1 & 3 MHz)
- One of the Lowest BNRs in the Industry
- Proudly Made in USA
- Three-Year Warranty

ORDERING DETAILS

400-020 TheraSound 3.5P

STANDARD ACCESSORIES

410-115 Therapy Hammer Ultrasound Applicator
(2 & 5 cm 2 in 1)

203-016 US Hospital Grade Power Cord 8 ft. (2.4m)
AloeSound Lotion 1 oz. (30 ml)

TECHNICAL SPECIFICATIONS

Power 115 VAC at 1 Amp / 230 VAC at .5 Amp

Frequency 50/60 Hz

US Frequency 1 & 3 MHz

US Output Pulsed and Continuous



US PRO 2000™ 2ND EDITION

The US Pro 2000 2nd Edition is a clinical grade therapeutic ultrasound device that generates deep heat through ultrasound energy. The pulsed sound waves travel deep into the tissue and induce vasodilation, which helps increase blood flow to the treated area and is found to help relieve pain and reduce muscle spasms.

Equipped with three output intensities and three treatment time selections, the US Pro 2000 2nd Edition can provide a wide range of therapy treatments. The head warming feature provides added patient comfort.

PRODUCT FEATURES:

- Operates in Continuous and Pulsed Mode (100% Duty Cycle)
- Head Warming Feature Pre-Heats the Sound Head Up to 95°F for Increased Patient Comfort
- Sound Head (Metal Part Only) Can Withstand Water Immersion for Up to 30 Minutes with Patient in No More Than One Meter Depth of Water
- Three Output Intensities Selections
- Three Treatment Times (5 Min, 10 Min, or 15 Min)
- Small, Ergonomic Design That Is Easy to Hold and Portable
- Comes Complete with Device, Ultrasound Gel, AC/DC Adapter, Hard Carrying Case, Quick Start Guide, and Instruction Manual
- One-Year Warranty
- Rx Only

ORDERING DETAILS

DU3035 US Pro 2000 2nd Edition Ultrasound Unit

STANDARD ACCESSORIES

LC2828 TheraSonic Conductive Gel 2 oz.
DU2035X AC Adapter
CC2035 Carrying Case

STANDARD ACCESSORIES

LC2828 TheraSonic Conductive Gel 2 oz.
DU2035X AC Adapter
CC2035 Carrying Case

TECHNICAL SPECIFICATIONS

Power 100-240 VAC
Frequency 50/60 Hz
US Frequency 1 MHz
US Output Pulsed and Continuous

US 1000™ 3RD EDITION

The US 1000 3rd Edition offers effective ultrasound therapy by generating deep ultrasonic waves within body tissues to treat pain, muscle spasms and joint stiffness. The portable unit comes as a complete kit with carrying case, AC adapter, ultrasound gel, and manual.



PRODUCT FEATURES:

- Therapeutic Ultrasound Device That Generates Deep Heat to Effected Area
- Device Must Be Used with Ultrasound Gel to Conduct the Ultrasound Energy to the Skin
- Three Power Output Settings
- Comes Complete with Ultrasound, AC Adapter, Carrying Case, Ultrasound Gel, and Instruction Manual
- One-Year Warranty
- Rx Only

ORDERING DETAILS

DU1025 US 1000 3rd Edition Portable Ultrasound Unit

STANDARD ACCESSORIES

LC2828 TheraSonic Conductive Gel 2 oz.
DU2035X AC Adapter
CC2035 Carrying Case

TECHNICAL SPECIFICATIONS

Power 100/240 VAC
Frequency 50/60 Hz
US Frequency 1 MHz
US Output Pulsed



MICROBLOCK ANTIMICROBIAL ELECTRODES

MicroBlock antimicrobial electrodes are the first and only electrodes on the market that help kill and inhibit the growth of common bacteria on and around the zone of inhibition of the electrode surface. These electrodes can play an important part in infection prevention for hospitals and healthcare facilities, making MicroBlock the logical choice to decrease the risk of healthcare associated infections while treating your patients.

PRODUCT FEATURES:

- Contains the Benzethonium Chloride (BEC), a Bactericidal Agent That Kills Disease-Causing Microbes and Common Antibiotic-Resistant Threats Such As:
 - E. Cloacae (CRE)
 - S. Aureus (MRSA)
 - E. Coli (ESBL)
 - E. Faecalis (VRE)
 - P. Aeruginosa and A. Baumannii (Multi-Drug Resistant)
- Same Comfortable Application and Delamination Resistance As Found in Richmar SuperStim and MultiStim Electrodes

ORDERING DETAILS

Cloth

400-872-MIC	2" (5 cm) Round, 4/pk
400-877-MIC	2" x 2" (5 cm x 5 cm) Square, 4/pk
400-879-MIC	2" x 3.5" (5 cm x 9 cm) Rectangle, 4/pk
400-881-MIC	3" (7.6 cm) Round, 4/pk
400-852-MIC	3" x 5" (7.6 cm x 13 cm) Rectangle, 2/pk

Patent Pending

Rx Only



SUPERSTIM PREMIUM ELECTRODES

The SuperStim premium electrode line is manufactured to have high tensile-strength lead wires. This is thanks to Leadstay, a proprietary bonding and moisture barrier layer protecting from delamination, particularly when used with moist heat therapy. A premium hydrogel provides comfort and uniform stimulation.

PRODUCT FEATURES:

- Available in Cloth or Foam and a Variety of Shapes and Sizes
- Strong Lead Wire Pull Strength
- Premium Dual Layer 35 mil Gel
- Gel Formulation Produces Maximum Dry-Out Resistance
- Optimum Distribution of Electrical Current
- Proudly Made in the USA

ORDERING DETAILS

Cloth

400-895	2" (5 cm) Round, 4/pk
400-899	2" x 2" (5 cm x 5 cm) Square, 4/pk
400-893	3" (7.6 cm) Round, 4/pk
400-864	1.25" (3 cm) Round, 4/pk
400-858	1.5" x 3.5" (4 cm x 9 cm) Rectangle, 4/pk
400-889	2" x 3.5" (5 cm x 9 cm) Rectangle, 4/pk
400-854	3" x 5" (7.6 cm x 13 cm) Rectangle, 2/pk
400-866	1.5" x 2.5" (4 cm x 6 cm) Oval, 4/pk

Foam

400-894	2" (5 cm) Round, 4/pk
400-898	2" x 2" (5 cm x 5 cm) Square, 4/pk
400-892	3" (7.6 cm) Round, 4/pk
400-869	1.25" (3 cm) Round, 4/pk
400-857	1.5" x 3.5" (4 cm x 9 cm) Rectangle, 4/pk
400-888	2" x 3.5" (5 cm x 9 cm) Rectangle, 4/pk
400-853	3" x 5" (7.6 cm x 13 cm) Rectangle, 2/pk
400-868	1.5" x 2.5" (4 cm x 6 cm) Oval, 4/pk



MULTISTIM ELECTRODES

Made with quality construction methods and components, MultiStim electrodes are engineered for e-stim treatments. Available in cloth or foam, in a variety of common sizes.

PRODUCT FEATURES:

- Available in Cloth or Foam and a Variety of Shapes and Sizes
- Strong Lead Wire Pull Strength
- Dual Layer 35 mil Gel
- Gel Formulation Produces Maximum Dry-Out Resistance
- Optimum Distribution of Electrical Current
- Proudly Made in the USA

ORDERING DETAILS

Cloth

400-872	2" (5 cm) Round, 4/pk
400-877	2" x 2" (5 cm x 5 cm) Square, 4/pk
400-881	3" (7.6 cm) Round, 4/pk
400-863	1.25" (3 cm) Round, 4/pk
400-856	1.5" x 3.5" (4 cm x 9 cm) Rectangle, 4/pk
400-879	2" x 3.5" (5 cm x 9 cm) Rectangle, 4/pk
400-865	1.5" x 2.5" (4 cm x 6 cm) Oval, 4/pk

Foam

400-873	2" (5 cm) Round, 4/pk
400-878	2" x 2" (5 cm x 5 cm) Square, 4/pk
400-882	3" (7.6 cm) Round, 4/pk
400-862	1.25" (3 cm) Round, 4/pk
400-855	1.5" x 3.5" (4 cm x 9 cm) Rectangle, 4/pk
400-880	2" x 3.5" (5 cm x 9 cm) Rectangle, 4/pk
400-861	1.5" x 2.5" (4 cm x 6 cm) Oval, 4/pk



HIVAMAT® ELECTRODES

Treating patients with gloved hands is the preferred application method of deep oscillation therapy, in order to mimic traditional manual therapy movements with the added benefits of deep oscillation. Hivamat Electrodes are designed specifically for the requirements and recommendations of deep oscillation therapy.

PRODUCT FEATURES:

- Specifically Engineered for Deep Oscillation Therapy
- Self-Adhesive Hivamat Electrodes with Snap Connection
- For Patient & Therapist Use
- Not Made with Natural Rubber Latex
- Proudly Made in the USA

ORDERING DETAILS

201-131

Hivamat Electrodes, 1 Set (100x)





IONTO⁺

Ionto+ Iontophoresis Electrodes come in 4 sizes and capacities, providing more available drug delivery than some competing Iontophoresis electrodes. Ionto+ Iontophoresis Electrodes contain synthetic pH stabilizing buffering agents in each electrode which guard against polarity shifts in the drug solution, helping to maintain consistent and efficient drug delivery throughout the Iontophoresis treatment.



201-452



201-455



201-449



201-453

PRODUCT FEATURES:

- Larger Adhesive Area Than Competitive Iontophoresis Electrodes
- Unique Silver-Silver Chloride Conductor (As Opposed to Carbon) Yields a More Efficient Drug Delivery
- Buffering Agents for Additional PH Stabilization
- Synthetic Gel Layer Secure Adhesion to Skin; Provides Maximum Current Dispersion, and Increases Patient Comfort Level
- Synthetic Buffering Provides Consistent PH Balancing; Helps Prevent the Effects of Hydrolysis
- No Water Needed to Activate the Adhesive
- Less Output Current Needed From Phoresor (Compared to Carbon Iontophoresis Electrodes)
- Proudly Made in USA

ORDERING DETAILS

201-452	Butterfly kit - 2.0 cc, 3mA-minute includes buffered return electrode 3.75" x 3.75" (9.5 x 9.5 cm), (12/Box)
201-455	Small kit - 1.5 cc, 3mA-minute includes buffered return electrode 2.75" x 2.75" (7 x 7 cm), (12/Box)
201-449	Medium kit - 2.5 cc, 4mA-minute includes buffered return electrode 3.25" x 3.25" (8.2 x 8.2 cm), (12/Box)
201-453	Large kit - 4.0 cc, 4mA-minute includes buffered return electrode 3.5" x 3.5" (9.0 x 9.0 cm), (12/Box)

STANDARD ACCESSORIES

201-383	Instruction Manual
201-376	Alcohol Prep Pad



CONDUCTIVE GLOVE

When attached to an electrotherapy device (not included), the conductive glove stimulates the entire hand through silver meshing. An effective alternative for the hand which is difficult to treat with self-adhesive electrodes. Commonly used to treat RSD, Carpal Tunnel and wrist and hand injuries. Should be used with conductive spray. Ambidextrous and one size fits most. Six month warranty. Rx only.



CONDUCTIVE ELBOW

When attached to an electrotherapy device (not included), the conductive sleeve is used to treat the entire elbow with conductive silver meshing. Effective for an area that is difficult to treat with self-adhesive electrodes. Commonly used to treat elbow problems, such as tennis elbow. Should be used with conductive spray. Ambidextrous and one size fits most. Six month warranty. Rx only.



CONDUCTIVE SOCK

When attached to an electrotherapy device (not included), the conductive sock treats the entire foot using conductive silver meshing. Effective for an area that is difficult to treat with self-adhesive electrodes. Commonly used to treat RSD, Tarsal Tunnel, neuropathies and basic foot/ankle injuries. Should be used with conductive spray. Ambidextrous and one size fits most. Six month warranty. Rx only.

ORDERING DETAILS

GU2856	Universal Conductive Glove
GU4017	Universal Conductive Elbow
GU4025	Universal Conductive Sock

STANDARD ACCESSORIES

EF3000WC2	3" Round White Cloth Electrodes, 4/pk
-----------	---------------------------------------

OPTIONAL ACCESSORIES

LS2177-2	Signa Conductive Spray, 2 oz
LS8812-2	Signa Conductive Spray, 8 oz

POWER SUPPLIES, CHARGERS, AND BATTERIES



DQ8000X



DQ8000X



DI1000IX



TA2000



TA2000C



TA9666-I



TA9050-I



TA5013-I



TA9225-I



TA9265-I

ORDERING DETAILS

Item No.	Description	Applicable Product(s)
203-016	US Hospital Grade Power Cord, 8 ft (2.4 m)	Winner EVO, TheraSound EVO
K36031H5183BB-D	UK Power Cord	Winner EVO
K5A031H5183BB	Euro Power Cord	Winner EVO
DQ8000X	Power Supply Cord, 96" medical grade	InTENSity EX4/CX4
DQ2001MGC	Adapter and Power Supply Cord, two-piece	Quattro 2.5
DQ9276MGC	Adapter and Power Supply Cord, two-piece	SoundCare Plus and Combo Care
DI000IX	AC Adapter	DI1010, DI3717, DI4738, DI8195, DI3302, DI9698, DT4222, DT4444
DI2000AC	AC Wall Adapter	DI1212, DI2738, DI2717, DI2195
TA2000	InTENSity Lithium Ion Battery	DI1212, DI2738, DI2717, DI2195
TA2000C	InTENSity Lithium Ion Battery Charger	InTENSity Lithium Ion Battery
TA9666-I	TENS Cell 1.5V AA Heavy Duty Batteries	
TA9050-I	TENS Cell 9V Heavy Duty Battery	
TA5013-I	TENS Cell 9V Super Alkaline Premium Battery	
TA9225-I	TENS Cell 9V 250 mAh Ni-MH Rechargeable Battery	
TA9265-I	TENS Cell 9V Charger, Two 9V Battery Capacity	

STIMULATION LEAD WIRES AND CONNECTION CABLES



A1717A



A1716A



WQ8112-SNGL



WQ6005PLUG

ORDERING DETAILS

Item No.	Description	Applicable Product(s)
A1717A	Lead Cord Set, 6 ft (2 m)	Winner EVO
A1716A	Lead Cord Set, 12 ft (4 m)	Winner EVO
WQ8112-SNGL	Combination Lead Wire 110", 1/pk	TheraTouch® EX4/CX4
WQ8112	Clinical Grade Lead Wires, Pin / Pin Connectors, 110", 4/pk	TheraTouch® EX4/CX4
WQ6005PLUG	Clinical Grade Lead Wires, Pin / Plug Connectors, 78", 4/pk	Quattro 2.5, and ComboCare
DQ8450CPLUG	Cable Connector, 2/pk	Quattro 2.5 and ComboCare
WW1001PLUG	Combination Lead Wire, Compatible with DQ8432CPlug	ComboCare

SEATING

SS3277



SS7677



ORDERING DETAILS

SS3277	Pneumatic Stool with Removable Back (Height Adjustable 18" - 23.5"), Black
SS7677	Pneumatic Stool (Height Adjustable 18" - 23.5"), Black

STIMULATION ELECTRODES AND GARMENTS



MicroBlock



SuperStim



MultiStim



Elastic Strap



Electrode Sponges



Sigma Spray



Rubber Electrodes

ORDERING DETAILS

MicroBlock Cloth

400-872-MIC	2" (5 cm) Round, 4/pk
400-877-MIC	2" x 2" (5 cm x 5 cm) Square, 4/pk
400-879-MIC	2" x 3.5" (5 cm x 9 cm) Rectangle, 4/pk
400-881-MIC	3" (7.6 cm) Round, 4/pk
400-852-MIC	3" x 5" (7.6 cm x 13 cm) Rectangle, 2/pk

SuperStim Cloth

400-895	2" (5 cm) Round, 4/pk
400-899	2" x 2" (5 cm x 5 cm) Square, 4/pk
400-893	3" (7.6 cm) Round, 4/pk
400-864	1.25" (3 cm) Round, 4/pk
400-858	1.5" x 3.5" (4 cm x 9 cm) Rectangle, 4/pk
400-889	2" x 3.5" (5 cm x 9 cm) Rectangle, 4/pk
400-854	3" x 5" (7.6 cm x 13 cm) Rectangle, 2/pk
400-866	1.5" x 2.5" (4 cm x 6 cm) Oval, 4/pk

SuperStim Foam

400-894	2" (5 cm) Round, 4/pk
400-898	2" x 2" (5 cm x 5 cm) Square, 4/pk
400-892	3" (7.6 cm) Round, 4/pk
400-869	1.25" (3 cm) Round, 4/pk
400-857	1.5" x 3.5" (4 cm x 9 cm) Rectangle, 4/pk
400-888	2" x 3.5" (5 cm x 9 cm) Rectangle, 4/pk
400-853	3" x 5" (7.6 cm x 13 cm) Rectangle, 2/pk
400-868	1.5" x 2.5" (4 cm x 6 cm) Oval, 4/pk

ORDERING DETAILS

MultiStim Cloth

400-872	2" (5 cm) Round, 4/pk
400-877	2" x 2" (5 cm x 5 cm) Square, 4/pk
400-881	3" (7.6 cm) Round, 4/pk
400-863	1.25" (3 cm) Round, 4/pk
400-856	1.5" x 3.5" (4 cm x 9 cm) Rectangle, 4/pk
400-879	2" x 3.5" (5 cm x 9 cm) Rectangle, 4/pk
400-865	1.5" x 2.5" (4 cm x 6 cm) Oval, 4/pk

MultiStim Foam

400-873	2" (5 cm) Round, 4/pk
400-878	2" x 2" (5 cm x 5 cm) Square, 4/pk
400-882	3" (7.6 cm) Round, 4/pk
400-862	1.25" (3 cm) Round, 4/pk
400-855	1.5" x 3.5" (4 cm x 9 cm) Rectangle, 4/pk
400-880	2" x 3.5" (5 cm x 9 cm) Rectangle, 4/pk
400-861	1.5" x 2.5" (4 cm x 6 cm) Oval, 4/pk

ORDERING DETAILS

ER2535B2	Rubber Electrodes, 2.5" x 3.5", 4/pk
ER2743B2	Rubber Electrodes, 2.75" x 4.3", 4/pk
ES2740Y2	Electrode Sponges, 2.75" x 4", 4/pk
ES3047Y2	Electrode Sponges, 3" x 4.75", 4/pk
EW3023BW2	Elastic Strap 3" x 23.5", 2/pk
EW3047B2	Elastic Strap 3" x 47", 2/pk
WS3147	Snap Pigtail Connector for Conductive Garments
LS2177-2	Sigma Spray, 2 oz.
LS8812-2	Sigma Spray, 8 oz.


DEEP OSCILLATION



RM00260



RM00381



RM00386



RM00382



RM00261



RM00379



RM00778



RM00262



RM00200



RM00348



RM00277



RM55035

ORDERING DETAILS

Item No.	Description	Applicable Product(s)
RM00261	Connection Cable	Hivamat 200 Evident
RM00381	Oscillator Head, 5 cm	Hivamat 200 Evident
RM00386	Oscillator Head, 9.5 cm	Hivamat 200 Evident
RM00382	Titanium Neutral Element	Hivamat 200 Evident
RM00229	Patient Therapy Card	Hivamat 200 Evident
RM00260	Patient Lead Wire	Hivamat 200 Evident
RM00379	Handheld Applicator	Hivamat 200 Evident, Hivamat 200 Portable
RM00778	US Hospital Grade Power Cord	Hivamat 200 Evident, Hivamat 200 Portable
RM00262	Connection Cable Snap	Hivamat 200 Evident, Hivamat 200 Portable
RM00200	Treatment Gloves (Box of 100)	Hivamat 200 Evident, Hivamat 200 Portable
201-131HV	Hivamat Electrodes, 2/pk (20x)	Hivamat 200 Evident, Hivamat 200 Portable
RM00348	Powder	Hivamat 200 Evident, Hivamat 200 Portable
RM00277	Battery Charger (US plug)	Hivamat 200 Portable
RM55035	Transport Bag	Hivamat 200 Portable
RM1012	Self Care Kit: Handheld Applicator, Connection Cable, Oscillator Head, 5 cm (2x), Oscillator Head, 9.5 cm (2x), Titanium Neutral Element	Hivamat 200 Portable

ULTRASOUND APPLICATORS AND TREATMENT



410-160



410-176



410115



410-200



400-026



A1220



DQ8002W1



DQ8002W5



LC5288



400-039

ORDERING DETAILS

Item No.	Description
410-160	5 cm Ultrasound Applicator (1 & 3 MHz)
410-176	10 cm Ultrasound Applicator (1 MHz)
410-115	Therapy Hammer 2-in-1 Ultrasound Applicator, 2 cm and 5 cm (1 & 3 MHz)
410-200	AutoSound Hands-Free Ultrasound Applicator
400-026	AutoGel Pads
A1220	AutoSound Swivel Adapter
DQ8002W1	1 cm ² Ultrasound Applicator with Cable
DQ8002W5	5 cm ² Ultrasound Applicator with Cable
DQ9275W1	1 cm ² Ultrasound Applicator with Cable
DQ9275W5	5 cm ² Ultrasound Applicator with Cable
LS5255	TheraSonic Conductive Gel, 5 L with 8 oz. Dispensing Bottle
LC5266	TheraSonic Conductive Gel, 16 oz.
LS5288	TheraSonic Conductive Gel, 8 oz.
LC2828	TheraSonic Conductive Gel, 2 oz.
400-039	AloeSound Lotion - 1 Gallon with Pump
400-042	AloeSound Dispenser Pump
201-1045	250 mL bottles

Applicable Product(s)

Winner EVO CM4 & CM2, TheraSound EVO
 Winner EVO CM4 & CM2, TheraSound EVO
 Winner EVO CM4 & CM2, TheraSound EVO
 Winner EVO CM4 & CM2, TheraSound EVO
 AutoSound Applicator
 AutoSound Applicator
 TheraTouch® CX4
 TheraTouch® CX4
 SoundCare Plus
 SoundCare Plus and ComboCare

LASER AND LIGHT THERAPY



LCS450



LCT450



410-060A



600-123

ORDERING DETAILS

Item No.	Description	Applicable Product(s)
LCS450	LaserPrism Cordless	Laser Module
LCT450	LaserPrism Tethered	Laser Module
600-123	Light Cluster Probe	Winner EVO
410-060A	Laser Module	

CLINICAL DEVICE PERIPHERALS



410-040



410-041



DQSWITCH



DQCART-2

ORDERING DETAILS

Item No.	Description	Applicable Product(s)
410-040	StoreMore Therapy Cart*	Winner EVO, TheraSound EVO with Cradle
410-041	TheraSound Cradle	TheraSound EVO with StoreMore Therapy Cart
DQSWITCH	Patient Interrupt Switch	TheraTouch® EX4/CX4
DQCART-2	Therapy Cart	TheraTouch® EX4/CX4

*Product discontinued 2019; available while supplies last.



SUBZERO™

Sub Zero with Cat's Claw provides fast long lasting relief from minor aches and pains associated with strenuous exercise, daily activities and the stresses of everyday life. Can be used as effective pain relief for arthritis, neck and back pain, sore muscles, sprains and spasms. Contains no artificial coloring, 3.9% Menthol and a unique herbal blend of Ilex, Boswellia and Uncaria Tomentosa



LZ1050



LZ1665



LZ4045



LZ3075



LZ75115

PRODUCT FEATURES:

- Portable Pain Relief
- Deep Penetrating, Cool Pain Relieving Gel
- Crystal Clear and Greaseless - No D&C
Blue or Yellow, Coloring Added
- Only Leading Topical Pain Reliever with Cat's Claw - a Unique Herbal Blend Proven to Provide Fast, Long-Lasting Relief
- Help Relieve Pain From: Strenuous Exercise and Sore Muscles, Strains and Spasms, Joint Pain, Back Pain, Shoulder Pain, Neck Pain, and Arthritis Pain
- Available in a Variety of Sizes

ORDERING DETAILS

LZ1050	SubZero 1 Gal. with Pump
LZ1665	SubZero 16 oz. with Pump
LZ4045	SubZero 4 oz. Tube
LZ3075	SubZero 3 oz. Roll-on
LZ5115	SubZero 5 ml Sample Pack

STRENGTHTAPE® KITS



STRENGTHTAPE is perfect for everyone from hard charging athletes to weekend warriors. It provides support & stability for muscles, joints, and tendons without limiting range of motion like a traditional brace. Featuring our exclusive AllSport adhesive, STRENGTHTAPE can be worn for up to 7 days in nearly any condition, including water.



PRODUCT FEATURES:

- Made of Breathable Cotton Fabric and Now Includes 10% More Stretch for Better Performance and Comfort
- Provides Support and Stability for Muscles, Joints, and Tendons Without Limiting Range of Motion Like a Traditional Brace
- Provides Compression for Injuries
- AllSport™ Adhesive Keeps Tape in Place for Up to 7 Days in Nearly Any Condition, Including Water
- Includes 6 Precut Strips with Easy-to-Follow Instructions

ORDERING DETAILS

6300-AF	StrengthTape Kit, Ankle/Foot
6300-BN	StrengthTape Kit, Back/Neck
6300-KN	StrengthTape Kit, Knee
6300-SHLDR	StrengthTape Kit, Shoulder
6300-MUS	StrengthTape Kit, Calf/Quad



STRENGTHTAPE® 5M ROLLS

STRENGTHTAPE Pre-Cut rolls include tape for up to 10 complete applications. Uncut rolls come with 5 meters of uncut tape. Each roll includes easy-to-follow instructions for the most common injuries.



Black



Beige



Light Blue



Red



White



Pink



Green



Royal Blue

PRODUCT FEATURES:

- Made of Breathable Cotton Fabric and Now Includes 10% More Stretch for Better Performance and Comfort
- Provides Support and Stability for Muscles, Joints, and Tendons Without Limiting Range of Motion Like a Traditional Brace
- Provides Compression for Injuries
- AllSport™ Adhesive Keeps Tape in Place for Up to 7 Days in Nearly Any Condition, Including Water
- Include Easy-to-Follow Instructions

ORDERING DETAILS

Precut Rolls - 5M

6305-510PC	StrengthTape Precut Roll, Black
6310-510PC	StrengthTape Precut Roll, Beige
6315-510PC	StrengthTape Precut Roll, Light Blue
6320-510PC	StrengthTape Precut Roll, Red
6325-510PC	StrengthTape Precut Roll, White
6330-510PC	StrengthTape Precut Roll, Pink
6370-510PC	StrengthTape Precut Roll, Green
6380-510PC	StrengthTape Precut Roll, Royal Blue

Uncut Rolls - 5M

6305-5UN	StrengthTape Uncut Roll, Black
6310-5UN	StrengthTape Uncut Roll, Beige
6320-5UN	StrengthTape Uncut Roll, Red
6380-5UN	StrengthTape Uncut Roll, Royal Blue



STRENGTHTAPE® 35M ROLLS

STRENGTHTAPE 35M jumbo uncut rolls are great for athletic trainers, therapists and sport enthusiasts alike looking for a larger, longer lasting roll for multiple applications. Great for the clinical setting.




Black



Beige



Red



Royal Blue

PRODUCT FEATURES:

- Made of Breathable Cotton Fabric and Now Includes 10% More Stretch for Better Performance and Comfort
- Provides Support and Stability for Muscles, Joints, and Tendons Without Limiting Range of Motion Like a Traditional Brace
- Provides Compression for Injuries
- AllSport™ Adhesive Keeps Tape in Place for Up to 7 Days in Nearly Any Condition, Including Water
- Include Easy-to-Follow Instructions

ORDERING DETAILS

6305-35UN	StrengthTape 35 m Uncut Roll, Black
63105-35UN	StrengthTape 35 m Uncut Roll, Beige
6320-35UN	StrengthTape 35 m Uncut Roll, Red
6380-35UN	StrengthTape 35 m Uncut Roll, Royal Blue



STRENGTHTAPE® DISPLAYS

StrengthTape Point-of-Purchase Displays are attractive and informative. With a small compact footprint, these displays will increase your cash sale revenue in a tiny space. A great addition to any waiting room in the clinic setting.



PRODUCT FEATURES:

- Made of Breathable Cotton Fabric and Now Includes 10% More Stretch for Better Performance and Comfort
- Provides Support and Stability for Muscles, Joints, and Tendons Without Limiting Range of Motion Like a Traditional Brace
- Provides Compression for Injuries
- AllSport™ Adhesive Keeps Tape in Place for Up to 7 Days in Nearly Any Condition, Including Water
- Include Easy-to-Follow Instructions

ORDERING DETAILS

- POP-ST-KITS-5 StrengthTape Kit Display, 36 Unit Display Includes: 6 Ankle/Foot, 6 Back/Neck, 6 Knee, 6 Shoulder, 6 Calf/Quad
- POP-ST-5MP StrengthTape Display, 12 Precut Rolls 5 m Includes: 3 Black, 3 Beige, 3 Pink, 3 Light Blue
- POP-ST-5MU StrengthTape Display, 12 Uncut Rolls 5 m Includes: 4 Black, 4 Beige, 2 Red, 2 Royal Blue



STRENGTHBRACE™ HINGED ROM ELBOW BRACE

Equipped with flexible bands that shape around the arms' contours, an adjustable strapping system, and telescopic length adjustments at the bicep and forearm, the StrengthBrace Hinge R.O.M. Elbow Brace offers robust customization to accommodate any patient's needs. The flexion and extension adjustments allow for a lesser or greater range of motion as the patient progresses through recovery.



PRODUCT FEATURES:

- Designed for Those with Soft Tissue Strains/ Repairs, and Post-Op Surgery Recovery
- Plush Padding for Comfortable Fit
- Telescopic Length Adjustments At Forearm and Biceps
- Fully Adjustable Strapping System
- Unique and Flexible Metal Bands Shape Around the Arms' Natural Contours
- Comes with Easily Clipped on Neck Strap
- One-Year Warranty

ORDERING DETAILS

BE3760L	StrengthBrace Hinge R.O.M. Elbow Brace, Left (HCPCS L3760)
BE3760R	StrengthBrace Hinge R.O.M. Elbow Brace, Right (HCPCS L3760)



STRENGTHBRACE™ UNIVERSAL HINGED KNEE BRACE

The StrengthBrace™ Universal Hinged Knee Brace provides stability and support for chronic knee pains and instability. The lightweight, breathable neoprene material with its unique wrap-around design provides optimal compression while the medial and lateral stays provide excellent support and promote proper alignment. The flexion and extension adjustments allow for a lesser or greater range of motion as the patient progresses through recovery.

PRODUCT FEATURES:

- Designed for Mild to Moderate Ligament Joint Pain, Sprain, Strains of the Knee Ligaments, Meniscus Injuries, Chondromalacia of the Patella, and Internal Derangement of the Knee.
- For Treatment of ACL, PCL, MCL, and LCL Instability, and Osteoarthritis
- Universal Wrap-Around Style Fits a Wide Range of Sizes and Requires Minimal Inventory to Satisfy a Diverse Patient Base
- Fits 12"-23": Measure Around Center of the Knee with Leg Extended
- Enclosed Tubular Donut Buttress
- Polycentric Hinge Provides Flexion/Extension and Can Be Easily Set for Desired Range of Motion
- Hinge Breaks At 0°, 10°, 20°, 40°, 60°, and 90°
- Shorter Knee Brace Design, Decreased Material Bulk for Increased Comfort Especially Prolonged Periods of Time
- Can Be Worn on Left or Right Knee
- One-Year Warranty

Knee Compression Sleeve

- Perfect Addition for Patients That Wear a Hinged Knee Brace for Compression Purposes
- Neoprene Lining Provides Extra Comfort and Helps Prevent Brace From Sliding or Moving

ORDERING DETAILS

BK1832-2	StrengthBrace Universal Hinged Knee Brace
BK2397-2	StrengthBrace Universal Knee Compression Sleeve



STRENGTHBRACE™ UNIVERSAL WRIST BRACE

Constructed of lightweight, breathable neoprene material, the StrengthBrace Universal Wrist Brace provides support and stabilization to those with wrist pain.

PRODUCT FEATURES:

- Universal Sizing
- Adjustable Dorsal Stay
- Lightweight, Breathable Material
- Easy Wrap Design for Quick Fitting
- Inset, Stretchable, Soft Sock for Ideal Comfort
- Not Made with Natural Rubber Latex
- One Year Warranty

ORDERING DETAILS

BK3908L	StrengthBrace Universal Wrist Brace, Left
BK3908R	StrengthBrace Universal Wrist Brace, Right



STRENGTHBRACE™ UNIVERSAL ANKLE

Ideal for post cast support or moderate sprains, the StrengthBrace Ankle Brace fits either the left or right foot and is available in multiple sizes. The lace-up enclosure provides a secure medial and lateral fit, while the unique figure-eight strapping provides extra support and stability.

PRODUCT FEATURES:

- Multiple Sizes Available
- Easy Fitting Based on Patient's Shoe Size
- Form Fitting Ankle Stay
- Unique Lacing and Figure-8 Strapping
- Breathable Mesh Under Foot to Prevent Sweating
- Mesh Near Heel to Prevent Chafing/Discomfort
- One-Year Warranty

ORDERING DETAILS

BA1902XS	StrengthBrace Canvas Lace-up Ankle Brace Extra Small (Men's: 6-7, Women's: 7-8)
BA1902S	StrengthBrace Canvas Lace-up Ankle Brace Small (Men's: 7-9, Women's: 8-10)
BA1902M	StrengthBrace Canvas Lace-up Ankle Brace Medium (Men's: 9-11, Women's: 10-12)
BA1902L	StrengthBrace Canvas Lace-up Ankle Brace Large (Men's: 11-13, Women's: 12-14)
BA1902XL	StrengthBrace Canvas Lace-up Ankle Brace Extra Large (Men's: 13-14, Women's: 14-15)



STRENGTHBRACE™ LITE

StrengthBrace LITE is a low-profile lumbar compression brace designed to prevent future injury while providing the support needed for chronic sprains and strains of the lower spine, spinal stenosis, and even post-surgery. This modular brace can conform to a variety of patients, with the ability to adjust to each of their unique support needs throughout their recovery process.

PRODUCT FEATURES:

- Support & Stability for a Variety of Conditions
 - Post-Operative Care
 - Chronic and Mechanical Back Pain
 - Reduction of Pain Associated with - Osteoarthritis and Other Spinal Conditions
 - Compression Fractures
 - Facet Syndrome
 - Degenerative Disc
- Billable Under HCPCS L0627 or L0642
- One Year Warranty

ORDERING DETAILS

BBSBL-0627 StrengthBrace Lite Economy Back Brace
(HCPCS L0627/L0642)

*Prescription Required: This is a Class II Medical Device. Federal law restricts these devices to sale by or on the order of a licensed health care professional.
Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.

OTTOBOCK SPINE BRACE



The SPINE Brace is the first line of defense for the treatment of low back pain without a prescription. Utilizing a patented one-handed rotating type pulley mechanism, it is designed specifically to provide comfortable support and hydrostatic lift to the lower back. A single front fabric fastener makes the brace easy to adjust with one hand in seconds for standing, sitting, walking or driving. The brace is made of a breathable mesh, promoting airflow to keep users cool and clean.



PRODUCT FEATURES:

- One-Step Pulley System Achieves Full Circumferential Compression As Opposed to the Favoring of One Side to Another
- Nylon Mesh Side Panels Provide Vertical Support
- Patented Rotating Type Pulley Mechanism Allows for True 6:1 Mechanical Advantage When Applying the Brace
- The Corset's Rear Power Unit Is Removable for Easy Care and Cleaning
- Retail Packaged
- 5 Year Warranty on the Mechanical Advantage Pulley System; 90 Day Warranty on the Belt/Brace

ORDERING DETAILS

BBSPINE-2025	Ottobock SPINE Brace, X-Small (20" - 25")
BBSPINE-2030	Ottobock SPINE Brace, Small (25" - 30")
BBSPINE-3035	Ottobock SPINE Brace, Medium (30" - 35")
BBSPINE-3040	Ottobock SPINE Brace, Large (35" - 40")
BBSPINE-4045	Ottobock SPINE Brace, X-Large (40" - 45")
BBSPINE-4550	Ottobock SPINE Brace, 2X-Large (45" - 50")
BBSPINE-5550	3 Ottobock SPINE Brace, X-Large (50" - 55")



THE OTTOBOCK® OPTION 1 BRACE

PRODUCT FEATURES:

- A Low-Profile Sagittal (Anterior-Posterior) Control Compression Brace That Is Ideal Low Lumbar Issues
- Covers L1 - S1 of the Spine
- Unique One-Hand Pulley System Reduces Friction and Offers a 5:1 Mechanical Advantage; for Each Pound of Pull, Five Pounds of Intra-Abdominal Compression Are Applied
- Anatomically Contoured Posterior Panel Can Be Heat Molded and Customized to Fit Individual Patients
- 11" Back Panel
- Fits Waist Circumference From 25" - 55"
- 5 Year Limited Warranty

ORDERING DETAILS

BBOPT-L0627	Ottobock Option 1 - AP Control Back Brace - 25" x 55"
BELT-EXT06	6-Inch Extension Belt - 55" x 61"
BELT-EXT12	12-Inch Extension Belt - 55" x 67"

THE OTTOBOCK® OPTION 2 BRACE

PRODUCT FEATURES:

- An Enhanced-Profile Sagittal (Anterior-Posterior) Control Compression Brace
- Covers T9 - S1 of the Spine (Half of the Thoracic Vertebrae, All of Lumbar Vertebrae and Sacrum)
- Unique One-Hand Pulley System Reduces Friction and Offers a 5:1 Mechanical Advantage; for Each Pound of Pull, Five Pounds of Intra-Abdominal Compression Are Applied
- Anatomically Contoured Posterior Panel Can Be Heat Molded and Customized to Fit Individual Patients
- 14" Back Panel
- Fits Waist Circumference From 25" - 55"

ORDERING DETAILS

BBOPT-L0631	Ottobock Option 2 - AP Control + Brace - 25" x 55"
BELT-EXT06	6-Inch Extension Belt - 55" x 61"
BELT-EXT12	12-Inch Extension Belt - 55" x 67"



THE OTTOBOCK® OPTION 3 BRACE



PRODUCT FEATURES:

- An Enhanced-Profile Sagittal (Anterior-Posterior) Control Compression Brace with Two Lateral Support Panels for Patients Requiring Additional Stabilization As Part of Their Treatment
- Covers T9 – S1 of the Spine (Half of the Thoracic Vertebrae, All of Lumbar Vertebrae and Sacrum)
- Unique One-Hand Pulley System Reduces Friction and Offers a 5:1 Mechanical Advantage; for Each Pound of Pull, Five Pounds of Intra-Abdominal Compression Are Applied
- Anatomically Contoured Posterior Panel Can Be Heat Molded and Customized to Fit Individual Patients
- 14" Back Panel
- Fits Waist Circumference From 25" – 55"

ORDERING DETAILS

BBOPT-L0637	OttoBock Option 3 - Control AP/L Brace - 25" x 55"
BELT-EXT06	6-Inch Extension Belt - 55" x 61"
BELT-EXT12	12-Inch Extension Belt - 55" x 67"

PNEU NECK II



PRODUCT FEATURES:

- Helps to Straighten the Neck and Alleviate Pressure on Tight Muscles, Pinched Nerves, Neck Discomfort, Tension Headaches, Osteoarthritis, Swollen Discs and More
- Designed with Three Layers to Prevent, Control and Diminish Painful Neck and Shoulder Discomfort
- Three Air Chambers Gently Inflate Up to 7" with Air Through Attached Inflation Bulb
- Easy-to-Use, Compact and Ideal for Traveling
- Soft Velour Cover Enhances Comfort
- One Size Fits Most
- This Product Contains Natural Rubber Latex Which May Cause Allergic Reactions
- One Year Warranty

ORDERING DETAILS

PNB3550	Pneu Neck II, Portable Cervical Device
---------	--

MEMORY FOAM KNEE SPACER



PRODUCT FEATURES:

- Stabilizes Pelvis and Lower Back, Restoring Alignment of Back, Hips and Knees
- Provides Comfortable Sleep
- Cushions Knee While Relieving Discomfort
- Adjustable Strap Secures Pillow to Leg
- Soft Terry-Velour Cover Removes for Easy Machine Washing
- 90-Day Warranty

ORDERING DETAILS

PC3420 Memory Foam Knee Space, 9" x 7" x 5"

COCCYX SEAT CUSHION



PRODUCT FEATURES:

- Cut-Out Reduces Pressure on Spine and Coccyx (Tail Bone) When Seated
- Provides Comfort and Support When Sitting for an Extended Period of Time
- Helps Relieve Pressure While Distributing Weight Evenly and Providing Proper Spinal Support
- Soft Black Terry-Velour Cover Is Removable for Easy Cleaning
- 90-Day Warranty

ORDERING DETAILS

COC-18 Coccyx Seat Cushion - Memory Foam, 18" x 16" x 3"

JACKSON ROLL-STYLE PILLOW



PRODUCT FEATURES:

- Fiberfill Holds its Shape
- Specially Shaped to Help Provide Comfort and Support
- Versatile Design - Can Be Used for Neck, Back, or Under Your Knees
- Ideal for Traveling
- 90-Day Warranty

ORDERING DETAILS

PP3115

Jackson Roll-Style Pillow, 16.9" x 7.1"



Richmar[®]
a COMPASSHEALTH company

Toll Free +1.888.549.4945
International +1.423.648.7730
4120 South Creek Road
Chattanooga, TN 37406, USA
richmarweb.com