

Vaginal Specula



cusco



graeve

The gynaecologist can select from two vaginal specula models and **different sizes**. Regardless which quality instrument is chosen, all specula excel **through their particular reliability and durability**. The Riester Vaginal Specula is used during gynecological examinations to gently widen the opening of the vagina to allow improved visualization and access to the cervix during a pelvic exam.

graeve vaginal speculum, external polished finish

cusco vaginal speculum, external matt finish

- › No reflections thanks to matt inner surface
- › Made of stainless steel type 18/8
- › Can be sterilized up to 134°C in an autoclave
- › Each speculum is supplied in a cardboard box

No. 5201	graeve, 75 x 20 mm	No. 5230	cusco virgins, 75 x 17 mm
No. 5202	graeve, 95 x 35 mm	No. 5231	cusco, 75 x 32 mm
No. 5203	graeve, 115 x 35 mm	No. 5232	cusco, 78 x 35 mm
		No. 5233	cusco, 100 x 35 mm

Please visit [riester.de](https://www.riester.de) for more information