



## Standard Features

- Painted beam
- Sturdy, enameled steel body
- Easy-to-clean ABS plastic tray
- Dual-read capacity; lb and kg
- Built-in measuring tape 0 to 23 in (0 to 58 cm)

## Specifications

**Capacity:**  
135 lb x 1 oz (62 kg x 20 g)

**Tray Dimensions:**  
(L x W x H)  
24 x 16 x 3 in  
(61 x 41 x 7.5 cm)

**Warranty:**  
Two-year limited

## Part Number/Price

Part #	Model #	Capacity	Price
131470	RL-MIS	135 lb x 1 oz (62 kg x 20 g)	Consult

## Shipping Weight

Part #	Model #	L x W x H	Package Wgt. (lb)	Dim Wgt. (lb)	Shipping Wgt. (lb)
131470	RL-MIS	28 x 20 x 14 in	30	57	57

## Dimensions

	Inches	Centimeters
A	24.0	61.0
B	27.0	68.6
C	21.5	54.5

