

SB6c-26 ROLL-IN BUDDY XXL

For roll in showers



The SB6c-26 mobile shower chair from Showerbuddy is designed for accessible or roll-in showers for larger users needing a wider seat and attendant care.

It can roll straight over the toilet and into the shower, reducing manual handling and transfers for the user and caregivers. The SB6c-26 has a 660mm wide seat cushion and can also be used as a commode chair.

The SB6c-26 Roll-InBuddy XXL is lightweight, durable and strong with a lifetime warranty on the frame.

Product Features:

- Weight capacity - 200kg
- Rolls over toilets Low maintenance
- Arms - removable / lockable / flip-back
- Height adjustable
- Footrests - removable and swing-away, height adjustable
- No rust aluminum frame (Grade 6061 T6)
- Stainless steel parts and hardware (Grade 304)
- Safety belts - lap or chest belts
- Caster Wheels - 100mm all locking
- Commode system
- Tool-free assembly

Product Dimensions:

- Chair Width - over caster wheels 670mm
- Chair Width - over arms 760mm
- Chair Width - Inside Arms 660mm
- Chair length - including footrest 970mm (730mm with footrests removed)
- Backrest Width 650mm
- Backrest Height from Cushion to top 380mm
- Toilet Height Clearance (maximum) 475mm
- Toilet width clearance at narrowest 560mm (at casters)
- Toilet clearance can depend upon the shape of the toilet pan
- Seat Cushion Dimensions 650 W x 430mm D
- Seat Height from Floor 510mm to 585mm (Adjustable)
- Seat Cushion Opening 175mm W x 335mm L
- Box Dimensions L 590mm x W 660mm x H 490mm
- Box Gross Weight 17.5kg

