

# Roscoe Mobility

## Extra Wide Bariatric Wheelchair

The Roscoe Extra Wide Bariatric Wheelchair features a robust design to accommodate patients with a wider seat and higher weight capacity requirement. The heavy-duty, dual cross bar frame has a streamlined look and is the industry standard when it comes to durability and reliability. It is available in multiple widths and weight limits up to 600 pounds. This chair offers a rugged, yet attractive option for the bariatric patient population.

**Available in  
26" and 28"**



### Product Features

- Steel frame with chrome finish
- Heavy-duty black padded reinforced nylon upholstery is durable and easy-to-clean
- Extra large, padded and upholstered armrests are removable
- Chart pocket on backrest provides additional convenience
- Includes tool-free adjustable front riggings with aluminum footplates (swing-away footrests or elevating legrests)
- Elevating padded legrests are equipped with a pneumatic piston for infinite elevating positions
- Chrome handrims with 24" solid rubber tires
- 8" x 2" semi-pneumatic casters are adjustable in three positions
- Precision sealed wheel bearings in front and rear ensure long lasting performance and reliability
- Easy push-to-lock wheel locks
- Limited lifetime warranty

### Product Specifications

- Product Weight <sup>2</sup>: 87 lb. / 89 lb.
- Overall Width with Arm: 34 1/2" / 37"
- Overall Depth: 34"
- Seat Width: 26" / 28"
- Seat Depth: 20"
- Seat-to-Floor Height<sup>1</sup>: 19"
- Back Height: 18"
- Back Style: Fixed
- Weight Capacity: 500 lb./ 600 lb.
- Code: K0007\*



Aluminum Footplates



Adjustable Casters



Chart Pocket



Dual Slide-Tube Frame

Item #	Frame	Seat Width	Armrests	Front Rigging
W7HD26S	Slide Tube	26"	Removable Desk-Length Arms	Swing-Away Footrests
W7HD26E	Slide Tube	26"	Removable Desk-Length Arms	Elevating Legrests
W7HD28S	Slide Tube	28"	Removable Desk-Length Arms	Swing-Away Footrests
W7HD28E	Slide Tube	28"	Removable Desk-Length Arms	Elevating Legrests

\* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.  
<sup>1</sup>NOTE: These heights can vary 1/4-inch due to tire wear and manufacturing tolerances.  
<sup>2</sup>NOTE: Weights based on wheelchair with desk-length arms and no front riggings. Weights will vary depending on how wheelchair is equipped.



**To Order:**  
 (800) 3-ROSCOE  
 fax: (440) 572-4261  
 roscoemedical.com