



WE CARE LIKE FAMILY™



NEW

SecureSafe™ Safety Syringes



Secure • Safe • Easy To Use

- Syringes with Safety Needle
- Safety Hypodermic Needles

New SecureSafe™ Safety Mechanism
 Dynarex New Safety Mechanism protects the end-user from accidental needle stick injuries

It's Easy, It's Safe, and Simple to Use!

SecureSafe™ Syringes with Safety Needle and Safety Hypodermic Needles


The SecureSafe™ safety mechanism has a quick-cover that allows you to easily cover and protect the needle with a press and snap!

1 Open

2 Inject

3 Engage quick cover and snap

See Our Website For Video Demos



See Reverse Side



Find us on:



dynarex® is your one-stop medical supply company since 1967

For more information contact your local Dynarex® distributor sales representative.

Or call **888-396-2739** to find your nearest distributor.

10 Glenshaw Street • Orangeburg, NY 10962 • www.dynarex.com

NEW COMPLETE LINE OF NEW SAFETY SYRINGES

SecureSafe™ Safety Syringes by Dynarex

Hypodermic Needles • Syringes with Needles

8916	SecureSafe Syringe w/ Safety Needle - 10cc - Luer Lock 21G, 1½" needle	6/100/Cs
8919	SecureSafe Syringe w/ Safety Needle - 5cc - Luer Lock 21G, 1½" needle	4/100/Cs
8921	SecureSafe Syringe w/ Safety Needle - 3cc Luer Lock 20G 1½" needle	8/100/Cs
8922	SecureSafe Syringe w/ Safety Needle - 3cc - Luer Lock 20G 1" needle	8/100/Cs
8923	SecureSafe Syringe w/ Safety Needle - 3cc - Luer Lock 21G, 1½" needle	8/100/Cs
8924	SecureSafe Syringe w/ Safety Needle - 3cc - Luer Lock 21G, 1" needle	8/100/Cs
8925	SecureSafe Syringe w/ Safety Needle - 3cc Luer Lock 22G, 1½" needle	8/100/Cs
8926	SecureSafe Syringe w/ Safety Needle - 3cc Luer Lock 22G, 1" needle	8/100/Cs
8927	SecureSafe Syringe w/ Safety Needle - 3cc Luer Lock 23G, 1" needle	8/100/Cs
8928	SecureSafe Syringe w/ Safety Needle - 3cc Luer Lock 25G, 5⁄8" needle	8/100/Cs
8929	SecureSafe Syringe w/ Safety Needle - 1cc Luer Lock 25G, 5⁄8" needle	8/100/Cs
8930	SecureSafe Syringe w/ Safety Needle - 1cc Luer Lock 27G, 1⁄2" needle	8/100/Cs
8931	SecureSafe Syringe w/ Safety Needle - 1cc Luer Lock 25G, 1" needle	8/100/Cs
8932	SecureSafe Syringe w/ Safety Needle - 1cc Luer Lock 25G, 1½" needle	8/100/Cs
8960	SecureSafe Safety Hypodermic Needle - 18G × 1" needle	10/100/Cs
8961	SecureSafe Safety Hypodermic Needle - 18G × 1½" needle	10/100/Cs
8962	SecureSafe Safety Hypodermic Needle - 19G × 1" needle	10/100/Cs
8963	SecureSafe Safety Hypodermic Needle - 19G × 1½" needle	10/100/Cs
8964	SecureSafe Safety Hypodermic Needle - 20G × 1" needle	10/100/Cs
8965	SecureSafe Safety Hypodermic Needle - 20G × 1½" needle	10/100/Cs
8966	SecureSafe Safety Hypodermic Needle - 21G × 1" needle	10/100/Cs
8967	SecureSafe Safety Hypodermic Needle - 21G × 1½" needle	10/100/Cs
8968	SecureSafe Safety Hypodermic Needle - 22G × 1" needle	10/100/Cs
8969	SecureSafe Safety Hypodermic Needle - 22G × 1½" needle	10/100/Cs
8970	SecureSafe Safety Hypodermic Needle - 23G × 1" needle	10/100/Cs
8971	SecureSafe Safety Hypodermic Needle - 23G × 1 ½" needle	10/100/Cs
8972	SecureSafe Safety Hypodermic Needle - 25G × 5⁄8" needle	10/100/Cs
8973	SecureSafe Safety Hypodermic Needle - 25G × 1½" needle	10/100/Cs
8974	SecureSafe Safety Hypodermic Needle - 26G × ½" needle	10/100/Cs
8976	SecureSafe Safety Hypodermic Needle - 27G × ½" needle	10/100/Cs

3/8



See Our Website
For Video Demos

* Our SecureSafe Alligator style safety syringes are comparable to Terumo SurGuard Syringes, Smiths Medical Needle-Pro EDGE, and BD Eclipse. Terumo SurGuard Syringes is a registered trademark of Terumo Corporation. Smiths Medical Needle-Pro EDGE is a registered trademark of Smiths Medical. BD Eclipse is a registered trademarks of Becton, Dickinson and Company

Now With Our New SecureSafe™ Safety Mechanics

See Front Side For Details

See our website for all our medical products! www.dynarex.com