



3pp[®] ThumSling[®] NP



The 3pp[®] ThumSling[®] NP is designed to relieve pain from CMC joint arthritis or instability with adjustable compression and moderate support. The ThumSling NP has a contoured strap that supports the thumb and wrist to control motion and relieve pain at the base of the thumb. Can also be used for ligament sprains.

Features

- The dense foam of the ThumSling NP offers the next level of support and control for those who need more control to reduce pain.
- Fully adjustable, the wearer determines how much compression to apply to stabilize the thumb and relieve pain during activity and rest.
- Resilient, breathable material has light stretch that contours to the hand without being bulky or restrictive.
- The flexible stay (included) can be added wherever extra support is needed.
- Hand or machine washable. Latex-free.

Order Information

Lay hand flat and measure around the MP joints (large knuckles).

	Right	Left	Fits (circumference at MP's)
Sm/Med	P2010-R23	P2010-L23	6" to 7 1/4"
Med/Lrg	P2010-R34	P2010-L34	7 1/4" to 9 1/4"

Watch our **instructional video on the ThumSling** on our website:
www.3pointproducts.com/product-videos-of-splints/3pp-thumsling-video



Call 888-378-7763

Mon-Fri 8:30am - 5:00pm EST for ordering details

www.3pointproducts.com