

Bathing > Bath Boards & Transfer Bench



Homecraft Alton Bath Board. Provides a stable platform for getting into or out of bathtub. Plastic-coated steel board is perforated to allow drainage. Brackets adjust 15"–23¾" apart. 27"W x 12"D. Supports up to 420 lbs.

AA1080



Savanah™ Shower Board. Similar to the Savanah Slatted Bathboard, but with two additional slats for extra comfort and support. May be also be used with the Savanah Bathboard Handle. Seat size 14"D. Supports up to 420 lbs.

557459 27" Internal Adjustment 20"–26"

557798 30" Internal Adjustment 23"–28"

Folding Shower Transfer Bench.

Legs individually adjust for correct seat height and level surface. Phenolic top ensures easy maintenance and a generous surface area makes transfers safe. Available in right or left L-shaped versions. No wall reinforcement needed. Installation by a qualified professional recommended. Measures 32"W x 22½"D x 15"D. Height adjusts 17"–19". Weight capacity: 900 lbs.

920501 Right

920502 Left



920502



Folds up easily



Savanah™ Slatted Bath Board. Smooth reinforced plastic slats provide a strong, comfortable, quick-draining design. Slats are fully sealed against the ingress of water and all fittings are stainless steel to resist corrosion. A serrated face and rubber buffers on brackets help secure the board in place. May also be used with the Savanah Bathboard handle. Depth 9¼". The 28" fits internal widths 21"–25" and the 30" board fits internal widths 23"–27". Supports up to 420 lbs.

557456 28" Fits internal widths 21"–25"

557457 30" Fits internal widths 23"–27"

Savanah Bath Board/Shower Board Handle. For ease of transfer and added user security.

557458



Includes built-in soap dish and drain holes



Portable Plastic Bath Board. Molded-plastic bath board fits securely over bathtub. Portable for storage and travel. Adjusts simply to fit standard bathtubs from 19"–28" wide. Rubber stops help protect tub from scratches. 31"W x 14¼"D. Supports up to 300 lbs.

A83014