

Diaper Scale


Model # NK2000

The NK2000 is the perfect size for diaper weighing applications within the nursery. This fully electronic scale utilizes the latest in load cell technology, providing accurate and repeatable weight data.



Features:

- Three push button design
- Two Mode Operations
- Easy to read $\frac{3}{4}$ " LCD Display
- Auto Shut Off in 2 minutes
- 4 "AA" alkaline battery operated
- Easy to clean surface
- 2 Year Limited Warranty



Three push button design:
Off, On/Tare, Mode

Specifications:

- Maximum Capacity: 2,000g/80oz.
- Dimension: 7" x 9" x 1"
- Weighing Surface: 5 $\frac{3}{4}$ "
- Units: grams/ounces
- Power Source: 4 AA Batteries (included)
- Battery Life: Average - 10,000 readings

