

## Seat Lift Chair Table



HCPCS:



FIGURE A

- Designed for use with a lift chair, standard recliner or couch
- Tabletop swivels to allow for easy access and departure (Figure A)
- Can be used for many activities such as eating, writing, crafts, working on the computer and reading
- Right or left positioning with 180° rotation accommodates seating near any wall
- Tabletop pivots 360° at the touch of a hand
- Teak wood grain laminate finish is attractive
- A built-in magazine rack provides added convenience

### PRODUCTS

### Description

13085LN	Overbed Table,LiftChair,Left,1/ea
13085RN	Overbed Table,LiftChair,Right,1/ea

### SPECIFICATIONS

### Description

Carton Shipping Weight	44 lbs.
Overall Dimensions at Base	N/A
Tabletop Adjustment	23.5" - 33"
Tabletop Surface Dimensions	15" (D) x 30" (W)
Warranty	One Year Limited
Weight	41 lbs
Weight Capacity	25 lbs