



Ritter®

MEDICAL

Seating

Comfort and style for every body

You spend a lot of time on the job—it might be work, but it doesn't have to feel like it. Ritter seating provides you with comfortable seating so you can concentrate on the job at hand rather than the pain in your back. And all of our seating comes with a variety of options, colors and sizes to complement your space (and your patients) just right.






Simple style, classic comfort

Classic Series Stools

Our most popular stools and the foundation of our seating line, the Classic Series is flexible to meet your seating needs and durable enough to withstand the clinical environment with just a touch of style to complement your space. Each stool features a polished, aluminum and black accented five-leg caster base with a hand ring, foot lever or manual height adjustment.



CLASSIC SERIES OPTIONS



- 01 Classic + Value Series Backrest
- 02 Classic Series Glides
- 03 Classic + Value Series Locking Casters
- 04 Classic Series Soft Rubber Casters
- 05 Foot Ring
- 06 6" Lab Extension


Stools you can work with Value Series Stools

Ritter Value Series stools give you the durability and style of our popular Classic Series stools in a more economical, feature-light package. Each stool comes with a black, plastic composite five-leg caster base and a single-lever height control (on pneumatic models).



- 01 Value Series Glides
- 02 Value Series Soft Rubber Casters
- 03 Classic + Value Series Locking Casters
- 04 Classic Series Soft Rubber Casters
- 05 Foot Ring
- 06 6" Lab Extension

VALUE SERIES OPTIONS



Strong,
comfortable
and coordinated



SIDE CHAIRS

We've improved upon our original side chair design with the 280 Side Chair and the 282 Bariatric Side Chair. Both side chairs feature the same contoured seat and back sections, powder-coated black steel tube frame with leveling feet, a wall-saver kickback leg design, an **18-inch seat height**, and high-impact plastic chair arms (optional for the 280). The 282 provides comfortable seating for people of all sizes with a **650-pound weight capacity** and a **30-inch wide seat**.



280
Side Chair




280
Side Chair + Arms



282
Bariatric Side Chair



01



02

01 Adjustable Foot

02 Kickback Leg



400 SERIES AIR-LIFT PHYSICIAN STOOLS

The plush, ergonomic design features a contoured seat with waterfall front edges that molds to the shape of your body for a perfect fit every time. Whether you choose the 425 with a hand ring or the 427 with a foot lever, these stools adjust in height to help you achieve a comfortable working position. (427 shown with optional back and arms)



281 BLOOD DRAW CHAIR

The slightly angled back and adjustable chair arms help improve patient comfort and add a sense of security. Arms are equipped with quick release mechanisms for simple positioning of height and inward/outward angle. The 281 is also available with a drawer, ideal for storing supplies.

SPECIFICATIONS

CLASSIC SERIES STOOLS

Seat Height:
18" - 24.5" (274-278)
19" - 25.5" (279)
Seat Diameter: 16.25"
Seat Cushion Thickness: 3.25"
Back Cushion Dimensions:
14" W x 7.5" H (275, 277)
Caster Base Diameter: 21.5"
Weight Capacity:
350 lb (274, 276)
225 lb (275, 277, 278, 279)

VALUE SERIES STOOLS

Seat Height:
18" - 24.5" (270-273)
Seat Diameter: 16.25"
Seat Cushion Thickness: 3.25"
Back Cushion Dimensions:
14" W x 7.5" H (271, 273)
Caster Base Diameter: 23"
Weight Capacity:
300 lb (270, 272)
225 lb (271, 273)

CLASSIC AND VALUE SERIES OPTIONS

Locking casters
Soft rubber casters
Glides
6" Lab height extension
Foot ring

425/427 AIR-LIFT PHYSICIAN STOOLS

Seat Height:
17.5" - 24" (425)
17" - 23" (427)
Seat Cushion Dimensions:
16.5" W x 17" D
Seat Cushion Thickness: 3"
Back Cushion Dimensions:
13" W x 9.5" H
(with back kit accessory)
Caster Base Diameter: 23"
Weight Capacity: 225 lb

281 BLOOD DRAW CHAIR

Seat Height: 20"
Seat Cushion Dimensions:
20" W x 16.5" D
Seat Cushion Thickness: 2"
Back Cushion Dimensions:
20" W x 16.5" H
Inside Arm Width: 21"
Outside Arm Width: 32"
Storage Drawer Capacity:
8.5" W x 15.25" D x 4.5" H
Weight Capacity: 400 lb

280 SIDE CHAIR

Seat Height: 18"
Seat Cushion Dimensions:
18.5" W x 19.5" D
Seat Cushion Thickness: 2"
Back Cushion Dimensions:
18.5" W x 13.5" H
Inside Arm Width: 17.5"
(on model with arms)
Outside Arm Width: 23.5"
(on model with arms)
Weight Capacity: 400 lb

282 SIDE CHAIR

Seat Height: 18"
Seat Cushion Dimensions:
30" W x 20" D
Seat Cushion Thickness: 2.5"
Back Cushion Dimensions:
30" W x 13.5" H
Inside Arm Width: 30"
Outside Arm Width: 36"
Weight Capacity: 650 lb

STANDARD UPHOLSTERY COLORS

-  Mist 813
-  Stone 814
-  Restful Path 858
-  Iron Ore 851
-  Lunar Gray 845
-  Obsidian 857
-  Oasis 843
-  Soothing Blue 855
-  Dream 847
-  Mineral 844
-  Healing Waters 848
-  Shaded Garden 853
-  Citrus 849
-  Sandy Retreat 850
-  Curative Copper 846
-  Cranberry 859
-  Dark Linen 856
-  Latte 860
-  Robust Brown 852
-  Deep Earth 854

ULTRAFREE[®] UPHOLSTERY COLORS

-  Stone 865
-  Restful Path 868
-  Obsidian 867
-  Harbor 861
-  Wheatgrass 862
-  Spice 863
-  Cranberry 869
-  Dark Linen 866
-  Latte 870
-  Branch 864

*ULTRALEATHER[®] by Ultrafabrics

The color examples shown are the best representation of the original material. Actual color may vary slightly. We strongly recommend that you contact Midmark Customer Experience at 1-800-MIDMARK to request a sample before placing your order.

Midmark is an ISO 13485 and ISO 9001 certified company. For more information or a demonstration, contact your Midmark dealer or call: 1-800-MIDMARK Fax: 1-877-725-6495 Outside the U.S.A. call: 1-937-526-3662 Fax: 1-937-526-8214 or visit our website at midmark.com

© 2018 Midmark Corporation, Dayton, OH. Products subject to improvement changes without notice. Litho in U.S.A.
007-10139-00 Rev. A1 (9/18)

Ritter®